

Rea® Azalea® Minor

Please contact DHCare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@invacare.com
 Or visit our website : www.invacare.eu.com



Rea® Azalea® Minor

Folding Frame Steel Manual Wheelchair - For Smaller Users



Description:

- Seat tilt with maintenance of the center of gravity and the correct position of the legrests: +1°/+24° and of the backrest: 0°/30° by pneumatic cylinders.
- Width-adjustable legrests.
- Heel straps. - Frame colour: Pearl gray (RAL 0515).
- Weight: from 32 kg.
- Maximum user weight: 75 kg.
- 2-year warranty period excluding wear parts



Date :

Purchase Order :

Customer Code :

Distributor Stamp:

DHCare Unique Customization Service

Please specify the quote or special production number :

☐ Please make your selection with a cross in the box

◆ Standard feature

○ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros
BNC0014	☐	Rea® Azalea® Minor for adults		
BNC0014	☐	Rea® Azalea® Minor for young children from 18 months		

PROPULSION TYPE

BNC0025	☐	SELF PROPELLED: 24" / 22" WHEELS WITH USER BRAKES	NC : 1074, 1076	◆
BNC0030	☐	SELF- AND ASSISTANT PROPELLED: 24" / 22" WHEELS WITH HUB BRAKES	NC : 1072, 1074, 1095	

SEAT

SEAT WIDTH

BNC0038	☐	340 mm - Adjustable from 340 to 390 mm		◆
---------	--------------------------	--	--	---

Seat width pre-adjustment

BNC0065	☐	Supplied with adjustment + 50 mm (i.e. 390 mm)		○
---------	--------------------------	--	--	---

Seat depth (SD)

BNC0095	☐	Adjustable from 380 to 450 mm		◆
---------	--------------------------	-------------------------------	--	---

SEAT HEIGHT

BNC0110	☐	400 mm	NC : 0476, 1062	○
BNC0120	☐	450mm	NC : 1055	◆

SEAT ANGLE

BNC0127	☐	NEGATIVE SEAT ANGLE -1 TO +23° (STD +1 TO +24°)		○
---------	--------------------------	---	--	---

SEAT TYPE

BNC0130	☐	DEPTH ADJUSTABLE		◆
---------	--------------------------	------------------	--	---

SEAT CUSHION

BNC0165	☐	CUSHION FLOSHAPE		◆
---------	--------------------------	------------------	--	---

SEAT UPHOLSTERY

BNC0191	☐	BLACK PLUSH TR35		
BNC0201	☐	BLACK PU TR26		◆
BNC0202	☐	BLACK SPACER TR36		○

BACKREST

BACKREST TYPE

BNC0217	☐	FLEX 2 BACKREST (TENSION ADJ. AZALEA MINOR)		○
BNC0218	☐	FLEX 2 BACKREST (FIXED, AZALEA MINOR)		◆



BACKREST COVER/ CUSHION

BNC0250	☐	CUSHION, LAGUNA (LATERAL SUPPORT)		◆
BNC0260	☐	CUSHION, MISTRAL 2 (WAIST SUPPORT)	NC : 0322	

BACKREST UPHOLSTERY

BNC0311	☐	BLACK PLUSH TR35		
BNC0321	☐	BLACK PU TR26		◆
BNC0322	☐	BLACK SPACER TR36		○

SW



0165



0217




0218





0250



0260

PUSHHANDLES					 0340 0345
BNC0340	☐	PUSH HANDLES BRACED		◆	
BNC0345	☐	PUSH BAR ANGLE ADJUSTABLE			
SEAT TILT AND BACKREST REGULATION TYPE					
BNC0350	☐	ASSISTANT OPERATED GAS PISTON SEAT AND BACKREST TILT.		◆	
BNC1263	☐	TILT LOCK (LOCKING OF PRE-SET ANGLE)		◆	
LEGRESTS / FOOTRESTS					
LEGRESTS / HANGERS LH TYPE					
BNC0380	☐	HANGER ANGLE ADJUSTABLE LH. WITH KNEE PAD AND CALF SUPPORT	NC : 0475, 0490, 0942	◆	 0380 / 0430
BNC0390	☐	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, LH		○	
BNC0410	☐	LEGREST AMPUTEE - ALU 2008 LH	NC : 0475, 0476, 0490, 0510, 0520, 0941, 0942		
LEGRESTS / HANGERS RH TYPE					
BNC0430	☐	HANGER ANGLE ADJUSTABLE RH. WITH KNEE PAD AND CALF SUPPORT	NC : 0475, 0550, 0942	◆	 0390 / 0440
BNC0440	☐	HANGER FIXED 80 DEG, SWING IN SWING OUT - ALU 2003, RH		○	
BNC0460	☐	LEGREST AMPUTEE - ALU 2008 RH	NC : 0475, 0476, 0550, 0570, 0580, 0941, 0942		
0410 / 0460					
FOOTRESTS					
FOOTRESTS ONE PIECE					
BNC0475	☐	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE, NO WIDTH ADJUSTEMENT - SHORT TUBE	NC : 0510, 0520, 0570, 0580, 0954, 0956		 0475 / 0476
BNC0476	☐	FOOTREST ONE PIECE DEPTH AND ANGLE ADJUSTABLE, NO WIDTH ADJUSTMENT - LONG TUBE	NC : 0510, 0520, 0570, 0580, 0954, 0956		
FOOTRESTS LH TYPE					
BNC0490	☐	FOOTREST LOWER PART, FIXED WITH SHORT TUBE AND HEEL STRAP, LH.		○	 0510 / 0570
BNC0510	☐	FOOTREST LOWER PART, FIXED WITH LONG TUBE AND HEEL STRAP, LH.	NC : 0954, 0956	○	
BNC0520	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, LH		◆	
FOOTRESTS RH TYPE					
BNC0550	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH SHORT TUBE AND HEEL STRAP, RH		○	 0520 / 0580
BNC0570	☐	FOOTREST LOWER PART, FIXED WITH LONG TUBE AND HEEL STRAP, RH.	NC : 0954, 0956	○	
BNC0580	☐	FOOTREST LOWER PART, DEPTH AND ANGLE ADJUSTABLE WITH LONG TUBE AND HEEL STRAP, RH		◆	
0520 / 0580					
ARMRESTS (Height adjustable, Autolock attachment included)					
ARMRESTS LH TYPE					
BNC0600	☐	STANDARD PAD, LH	NC : 1233	○	 0600 / 0630
BNC0606	☐	CONCAVE GUTTER ADJUSTABLE IN HEIGHT AND DEPTH OVE 70 MM WITHOUT TOOLS (COMPATIBLE WITH SHELF)	NC : 1233		
BNC0610	☐	PADDED, LH.	NC : 0680, 0690, 1230, 1231, 1233	◆	
BNC0615	☐	HEMIPLEGIC ARMRESTS, LH.	NC : 1230, 1231, 1233, 1240, 1245		
BNC0620	☐	WITHOUT ARMREST, LH	NC : 0680, 0690, 1230, 1231, 1240, 1245	○	
0610 / 0640					
ARMRESTS RH TYPE					
BNC0630	☐	STANDARD PAD, RH	NC : 1232	○	 0615 / 0645
BNC0636	☐	CONCAVE GUTTER ADJUSTABLE IN HEIGHT AND DEPTH OVE 70 MM WITHOUT TOOLS (COMPATIBLE WITH SHELF)	NC : 1233		
BNC0640	☐	PADDED, RH.	NC : 0680, 0690, 1230, 1231, 1232, 1240, 1245	◆	
BNC0645	☐	HEMIPLEGIC ARMRESTS, RH.	NC : 1230, 1240, 1245		
BNC0650	☐	WITHOUT ARMREST, RH	NC : 0680, 0690, 1230, 1231, 1240, 1245	○	

SEATING OPTIONS						
ARMREST OPTIONS						
BNC0680	☐	SIDE SUPPORT COVER				
BNC0690	☐	CUSHION TO SIDE REST COVER	ONLY WITH 0680			
SEAT OPTIONS						
BNC0710	☐	INCONTINENCE COVER				
BNC0720	☐	ABDUCTION CUSHION				
BACKREST OPTIONS						
BNC0730	☐	HEADREST, SIDEWISE ADJUSTABLE	NC : 0742, 0744, 0746, 0750, 0770	☐	0720	
BNC0750	☐	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE	NC : 0740, 0742, 0744, 0746, 0770, 0780		0730	
BNC0740	☐	COVER FOR HEADREST			0750	
BNC0770	☐	NECKREST		☒	0770	
BNC0782	☐	UNIVERSAL HEAD/NECKREST HOLDER				
BNC0790	☐	HANDLES CLINIC KIT FOR HEAD-/NECKREST				
TRUNK SUPPORT LH						
BNC0800	☐	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (LH)				
BNC0880	☐	TRUNK SUPPORT SWING-AWAY TR26 LEFT CPLT (V2025)& BRACKET, FIX BACK, PASSIVE		NEW	0782	
TRUNK SUPPORT RH						
BNC0810	☐	TRUNK SUPPORT MULTI ADJUSTABLE WITH STANDARD ARM (RH)			0800 / 0810	
BNC0890	☐	TRUNK SUPPORT SWING-AWAY TR26 RIGHT CPLT (V2025)& BRACKET, FIX BACK, PASSIVE		NEW		
BACKREST OPTIONS (Following)						
BNC0915	☐	TILT SCALE, BACKREST RECLINE				
BNC0920	☐	WEDGES (ONE PAIR)				
BNC0954	☐	Bodypoint Open Buckle Ankle Supports, Side Adjustable, Medium (pair) - Comes loose with bindings				
BNC0956	☐	Bodypoint ankle braces open with buckle, adjustable on the side, Small (pair) - delivered loose with bindings			0880 / 0890	
					0954/0956	
LEGREST OPTIONS						
BNC0941	☐	COVER FOOTPLATE				
BNC0942	☐	CALF STRAP PADDED				
CHASSIS						
CHASSIS COLOUR						
BNC0944	☐	PEARL GREY		☒	0944	
BNC0946	☐	AZURITE		☐	0946	
BNC0947	☐	SAND		☐	0947	
BNC0948	☐	ELECTRIC GREEN		☐	0948	
BNC0949	☐	HAPPY RED		☐	0949	
CASTOR ATTACHMENT TYPE						
BNC0960	☐	Fixed		☒		
FRONT CASTORS						
BNC0982	☐	5,5"x1.3/4" - 140X40 BLACK MEDIUM HARD				
BNC0992	☐	6"x1.1/4" - 150X30 BLACK MEDIUM HARD		☒		
REAR WHEELS						
REAR WHEELS SIZE						
BNC1055	☐	22" / 560 MM - QR REAR WHEEL ATTACHMENT		☐		
BNC1062	☐	24" / 600 MM - QR REAR WHEEL ATTACHMENT		☒		
REAR WHEELS TYPE						
BNC1072	☐	STANDARD (BLACK HUB, COLORLESS RIM & ZINC PLATTED SPOKES)	NC : 1120	☒		
BNC1074	☐	REINFORCED (COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES)		☐		
BNC1076	☐	DRUM BRAKE (22/24" : COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES / 16" : BLACK PLASTIC HUB & RIM)	NC : 1086	☐		
TYRE TYPE						
BNC1082	☐	PNEUMATIC PROFILED BLACK 1" 3/8		☐		
BNC1086	☐	PNEUMATIC LOW PROFILE BLACK 1"		☐		
BNC1092	☐	TYRES SOLID BLACK 1" 3/8		☒		

HANDRIM TYPE					
BNC1100	☐	ALUMINUM		◆	
BNC1120	☐	PLASTIC COATED			
WHEEL OPTIONS					
SPOKEGUARD					
BNC1160	☐	FROST			
SAFETY					
BNC1170	☐	USER MANOUEVERED TYRE BRAKES		◆	 1200
BNC1200	☐	PELVIC BELT - PLASTIC BUCKLE			
POSITIONING OPTIONS					
Think about Bodypoint positioning accessories: belts, harnesses, ankle supports.					
OPTIONS					
TABLE TRAY					
BNC1230	☐	TABLE TRAY INCL MOUNTING KIT	NC : 1237		
BNC1231	☐	TABLE TRAY INCL MOUNTING KIT, WITH ELBOW PADS	NC : 1237		 1230 1231
BNC1232	☐	HALF TABLE TRAY, CLEAR, FLIP UP, LH (ARMREST WITH STRAIGHT NARROW PAD INCL.)	NC : 1235		 1232 / 1233
BNC1233	☐	HALF TABLE TRAY, CLEAR, FLIP UP, RH (ARMREST WITH STRAIGHT NARROW PAD INCL.)	NC : 1235		
BNC1240	☐	SWING AWAY TABLE	NC : 1235, 1237		
TABLE TRAY CUSHION					
BNC1235	☐	TABLE TRAY CUSHION			 1237 1240
BNC1237	☐	PADDING			
OTHER OPTIONS					
BNC1300	☐	HARNESS ATTACHMENT			 1300
BNC1310	☐	ATTACHMENT WHEN USING THE PELVIC BELT FOR			
RAMPS					
1654128	☐	Basic ramp 2 telescopic parts in aluminum and non-slip coating. Folded dimensions: 1210 x 194 mm. 8.8 kg per pair - maximum capacity: 270 kg			 1654128
1654129	☐	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg			 1654129
1654130	☐	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg			 1654130