

# SEAT CUSHIONS ADJUSTABLE SKIN PROTECTION & POSITIONING



ROHO's ISOFLO® Memory Control, located at the front of all Select Series cushions, allows air to flow through the cushion in two ways:

- With the ISOFLO Memory Control open, air can flow freely throughout the entire cushion, similar to a ROHO Single Compartment Cushion.
- With the ISOFLO Memory Control closed, air becomes locked into four separate quadrants offering positioning capabilities.

## ROHO sets the standard for performance in the wheelchair seating industry.

Based on over 40 years of research and individual experience, ROHO® DRY FLOATATION® Technology provides unparalleled performance in facilitating blood flow by fitting and matching the shape of an individual.

Select Series cushions are one of ROHO's most popular product lines because of the revolutionary ISOFLO Memory Control, which offers shape fitting capabilities with the push of a knob. The ISOFLO Memory Control offers a seated individual an adjustable, quick and simple positioning and stability solution.

Additional adjustable skin protection and positioning options, including a series of cushions available with dual valves, provide individuals with increased positioning and stability functions based on their needs.

# ROHO® QUADTRO SELECT® HIGH PROFILE® Cushion

The QUADTRO SELECT HIGH PROFILE Cushion provides maximum skin and soft tissue protection. With the ISOFLO® Memory Control function, the individual has the ability to easily set and save positioning settings by regulating air in four quadrants. The revolutionary ISOFLO Memory Control makes it easy to adjust the individual's posture to the desired sitting position. By pushing the knob to the closed position, the user can maintain the desired sitting posture for maximum function.

The QUADTRO SELECT HIGH PROFILE Cushion is 4.25" tall, providing the maximum immersion for patients.

Based on clinical, scientific and/or engineering evidence, this product may be suitable for individuals who:

- Have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure ulcer (including deep tissue injury), skin grafts, and/or myocutaneous flap surgery
- Currently have skin/soft tissue breakdown and/or any category/stage pressure ulcer (including deep tissue injury)
- Have normal or impaired sensation
- Require accommodation of pelvic asymmetry, up to approximately 3 inches
- Require increased sitting stability
- Require correction of pelvic asymmetry (medial lateral and/or anterior posterior)
- Require different positioning of the left and right femur



## DIMENSIONS

Width: 12.75" to 24.5"  
 Height: 4.25"  
 Weight: Approx. 3.5 lbs.

## MEDICARE CODES

E2624  
 E2625

## WEIGHT LIMIT

No weight limit, yet the cushion must be properly sized to the individual.

## LIMITED WARRANTY

Cushion: 36 Months  
 Cover: 6 Months

■ SKIN PROTECTION ■ POSITIONING ■ ADJUSTABILITY ■ ADAPTABILITY ■ COMFORT

# ROHO® QUADTRO SELECT® LOW PROFILE® Cushion

The QUADTRO SELECT LOW PROFILE Cushion is the lowest profile of all the Select Series cushions with a cell height of 2.25" tall. The LOW PROFILE cushion offers ROHO's DRY FLOATATION® Technology providing maximum skin and soft tissue protection, while offering the benefits of the ISOFLO® Memory Control.

Based on clinical, scientific and/or engineering evidence, this product may be suitable for individuals who:

- Have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure ulcer (including deep tissue injury), skin grafts, and/or myocutaneous flap surgery
- Currently have skin/soft tissue breakdown (including deep tissue injury)
- Have normal or impaired sensation
- Require accommodation of pelvic asymmetry, approximately up to 1 inch
- Require increased sitting stability
- Require correction of pelvic asymmetry (medial lateral and/or anterior posterior)
- Require different positioning of the left and right femur



## DIMENSIONS

Width: 11" to 24.5"  
Height: 2.25"  
Weight: Approx. 2.25 lbs.

## MEDICARE CODES

E2624  
E2625

## WEIGHT LIMIT

No weight limit, yet the cushion must be properly sized to the individual.

## LIMITED WARRANTY

Cushion: 36 Months  
Cover: 6 Months

■ SKIN PROTECTION ■ POSITIONING ■ ADJUSTABILITY ■ ADAPTABILITY ■ COMFORT