



44 **activity chair.**

Invoice address:

Dispatch address: See invoice address

Order
Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The dynamic therapy chair.

- Seat and back unit incl. pad
- Back angle adjustment via gas pressure spring from (-)5° to 20°
- Back height and seat depth adjustment
- Height-adjustable armrests
- Incl. 2-point pelvic harness
- Tilt-in-space adjustment via gas pressure spring
 - for the standard base (-) 15° to 15°
 - for the Hi/Lo base (-) 15° to 25°
- Frame colour: Effect black

HMV-No.: 26.11.05.0xxx



Basic model

activity chair.

Seat unit (standard)

Tilt-in-space adjustment via gas pressure spring (-) 15° to 20°

	Price / piece
Size 0 (72 00 200) (R821)	998 €
Size 1 (72 01 200) (R841)	1.019 €
Size 2 (72 02 200) (R861)	1.165 €

activity chair.




Seat unit (dynamic)

Back angle adjustment via gas pressure spring from (-)5° to 20°; a gas pressure spring enables self-generated movement to be accommodated - the spring option offers the user 10° of movement.


	Price / piece
Size 0 (72 00 300) (R822)	1.139 €
Size 1 (72 01 300) (R842)	1.165 €
Size 2 (72 02 300) (R862)	1.347 €

Dimensions				
		Size 0	Size 1	Size 2
Seat depth		20 - 30 cm	28 - 41 cm	38 - 51 cm
Seat width		30 cm	36 cm	43 cm
Seat width (with hip support pads)		18 - 23 cm	22 - 29 cm	28 - 36 cm
Chest width		17 - 29 cm	17 - 30 cm	24 - 36 cm
Back height		32 - 39 cm	39 - 48 cm	48 - 61 cm
Height headrest - seat		37 - 53 cm	44 - 62 cm	50 - 75 cm
Height armrest - seat		13 - 19 cm	18 - 27 cm	18 - 27 cm
Lower leg length		18 - 30 cm	20 - 38 cm	29 - 54 cm
Seat height above floor	with short / long feet	24 - 32 cm / 47 - 55 cm	32 - 42 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
	with short / long feet and dual castors	34 - 42 cm / 47 - 55 cm	41 - 51 cm / 48 - 58 cm	42 - 52 cm / 51 - 61 cm
	with Hi/Lo base	25,5 - 63,5 cm	30,5 - 58,5 cm	40,5 - 61 cm
Frame width	with short / long feet	53 / 58 cm	59 / 62 cm	68 / 68 cm
	with short / long feet and dual castors	57 / 60 cm	62 / 65 cm	71 / 71 cm
	with Hi/Lo base	66 cm	68,5 cm	73,5 cm
Tilt-in-space	with standard base / with Hi/Lo base			(-)15° - 15° / (-)15° - 25°
Back angle		(-)5° - 20°		
Total length/frame depth		54-73 cm	57-73 cm	77-90 cm
Total width		53-60 cm	59-65 cm	68-70,5 cm
Total height incl. headrest		61-108 cm	76-120 cm	92-136 cm
Total height without headrest		56-94 cm	71-106 cm	90-122 cm
Recommended body height		81-122 cm	102-142 cm	127-188 cm
max. load		34 kg	68 kg	113 kg
Weight of standard base frame		12,2 kg	15,9 kg	24,5 kg
Weight of Hi/Lo frame		25,4 Kg	39 Kg	44,5 kg

The dimensions given apply to both seat units with all bases. If certain dimensions differ between the base frames, these are listed specifically.

Base frame			
<p>Standard base activity chair. incl. height-adjustable, non-slip feet; seat angle adjustment by gas pressure spring from (-) 15° to 15°</p> 	<p>Dynamic base activity chair. incl. height-adjustable, non-slip feet; seat angle adjustment by gas pressure spring from (-) 15° to 15°; a gas pressure spring enables self-generated movement to be accommodated - the spring option offers the user room to move.</p> 	<p>Hi/Lo base activity chair. incl. 4 double castors (100 mm, braked, 1x with directional lock); height adjustment (by gas pressure spring for size 0 and hydraulically for sizes 1+2); seat angle adjustment by gas pressure spring from (-) 25° to 15°</p> 	
Price / piece	Price / piece	Price / piece	Price / piece
Size 0 (72 00 400) (R823) 709 €	Size 0 (72 00 500) (R824) 817 €	Size 0 (72 00 600) (R833) 2.204 €	
Size 1 (72 01 400) (R843) 721 €	Size 1 (72 01 500) (R844) 846 €	Size 1 (72 01 600) (R853) 3.094 €	
Size 2 (72 02 400) (R873) 846 €	Size 2 (72 02 500) (R874) 986 €	Size 2 (72 02 600) (R883) 3.283 €	

Pad colour			
<p>Blue</p> 	<p>Piece</p> <p>for size 0 (72 05 201) (R816) incl.</p> <p>for size 1 (72 00 201) (R836) incl.</p> <p>for size 2 (72 02 201) (R876) incl.</p>	<p>Red</p> 	<p>Piece</p> <p>for size 0 (72 05 301) (R817) incl.</p> <p>for size 1 (72 00 301) (R837) incl.</p> <p>for size 2 (72 02 301) (R817) incl.</p>




Feet			
	Short feet – standard configuration, if nothing else stated	Long feet – for all sizes	
		Price / Set	Price / Set
	for size 0 (72 04 004) (R803) for size 1 (72 00 004) (R893) for size 2 (72 05 004) (R808)	incl. incl. incl.	72 01 004 (R894) incl.

Feet				
	4 Dual castors (75 mm) short version, with brakes, 1 x with direction lock		4 Dual castors (75 mm) long version, with brakes, 1 x with direction lock, for all sizes	
			Upcharge / set	Upcharge / set
	for size 0 (72 06 004) (R805) for size 1+2 (72 02 004) (R895)		382 € 382 €	72 03 004 (R896) 382 €






Accessories								
Headrest (Standard) depth, angle and height adjustable 	Headrest (Shell-shaped) depth, angle and height adjustable 	Headrest (Lateral guide) depth, angle and height adjustable, wing pads can be adjusted from almost flat to 90° 	Price / piece					
			72 00 010 (R898)	342 €	72 00 011 (R897)	367 €	72 00 018 (R881)	471 €

Accessories																							
Lateral pads – height and width adjustable 	<table border="1"> <thead> <tr> <th></th> <th>Depth</th> <th>Height</th> </tr> </thead> <tbody> <tr> <td>for size 0-2</td> <td>13 cm</td> <td>10 cm</td> </tr> <tr> <td>for size 1-2</td> <td>16,5 cm</td> <td>13 cm</td> </tr> </tbody> </table>		Depth	Height	for size 0-2	13 cm	10 cm	for size 1-2	16,5 cm	13 cm	Lateral pads – height and width adjustable, with chest strap 	<table border="1"> <thead> <tr> <th></th> <th>Depth</th> <th>Height</th> </tr> </thead> <tbody> <tr> <td>for size 0+1</td> <td>13 cm</td> <td>10 cm</td> </tr> <tr> <td>for size 2</td> <td>16,5 cm</td> <td>13 cm</td> </tr> </tbody> </table>		Depth	Height	for size 0+1	13 cm	10 cm	for size 2	16,5 cm	13 cm	Push handles – recommended for all activity chair. with dual castors and models with Hi/Lo base 	Price / pair
			Depth	Height																			
		for size 0-2	13 cm	10 cm																			
		for size 1-2	16,5 cm	13 cm																			
	Depth	Height																					
for size 0+1	13 cm	10 cm																					
for size 2	16,5 cm	13 cm																					
for size 0-2 (72 01 005) (R828) for size 1-2 (72 02 005) (R838)	300 € 324 €	for size 0 (72 00 006) (R819) for size 1 (72 01 006) (R829) for size 2 (72 02 006) (R839)	382 € 382 € 409 €	72 00 012 (R878)	243 €																		

Accessories																							
Armrests – height adjustable 	<table border="1"> <thead> <tr> <th></th> <th>Width</th> <th>Length</th> </tr> </thead> <tbody> <tr> <td>for size 0</td> <td>6 cm</td> <td>22 cm</td> </tr> <tr> <td>for size 1+2</td> <td>7,5 cm</td> <td>30,5 cm</td> </tr> </tbody> </table>		Width	Length	for size 0	6 cm	22 cm	for size 1+2	7,5 cm	30,5 cm	Hip support pads – padded and angle adjustable (only for use together with the armrests) 	<table border="1"> <thead> <tr> <th></th> <th>Height</th> <th>Length</th> </tr> </thead> <tbody> <tr> <td>for size 0</td> <td>8 cm</td> <td>14 cm</td> </tr> <tr> <td>for size 1+2</td> <td>10 cm</td> <td>23,5 cm</td> </tr> </tbody> </table>		Height	Length	for size 0	8 cm	14 cm	for size 1+2	10 cm	23,5 cm	Therapy tray – only for use together with the armrests 	Price / pair
			Width	Length																			
		for size 0	6 cm	22 cm																			
		for size 1+2	7,5 cm	30,5 cm																			
	Height	Length																					
for size 0	8 cm	14 cm																					
for size 1+2	10 cm	23,5 cm																					
for size 0 (72 05 007) (R827) for size 1+2 (72 00 007) (R868)	382 € 382 €	for size 0 (72 05 008) (R814) for size 1+2 (72 00 008) (R867)	382 € 382 €	for size 0 (72 00 009) (R667) for size 1 (72 01 009) (R668) for size 2 (72 02 009) (R669)	449 € 449 € 449 €																		

Accessories															
Backrest extension 	Spinal support and tube edge 	Abduction wedge depth adjustable 													
		<table border="1"> <thead> <tr> <th></th> <th>Width</th> <th>Depth</th> <th>Height</th> </tr> </thead> <tbody> <tr> <td>for size 0+1</td> <td>7,5 cm</td> <td>15 cm</td> <td>9 cm</td> </tr> <tr> <td>for size 2</td> <td>10 cm</td> <td>19 cm</td> <td>10 cm</td> </tr> </tbody> </table>		Width	Depth	Height	for size 0+1	7,5 cm	15 cm	9 cm	for size 2	10 cm	19 cm	10 cm	
	Width	Depth	Height												
for size 0+1	7,5 cm	15 cm	9 cm												
for size 2	10 cm	19 cm	10 cm												
	Price / piece	Price / Set	Price / piece												
for size 0+1 (72 01 016) (R849)	193 €	for size 0 (72 00 017) (R886)	167 €												
for size 2 (72 02 016) (R856)	193 €	for size 1 (72 01 017) (R887)	167 €												
		for size 2 (72 02 017) (R888)	167 €												
		for size 0+1 (72 01 013) (R884)	206 €												
		for size 2 (72 02 013) (R891)	206 €												

Accessories			
Knee lateral guide 		Knee support pad with key lock, depth adjustable 	
	Price / pair	Price / piece	Price / piece
72 00 014 (R892)	300 €	for size 0 (72 00 015) (R854)	338 €
		for size 1 (72 01 015) (R858)	338 €
		for size 2 (72 02 015) (R859)	367 €

Harness			
Chest strap 	Chest-shoulder strap  normal  thin	Thigh straps 	Groin harness 
Price / piece	Price / piece	Price / pair	Price / piece
for size 0 (72 00 020) (R802) 149 €	normal for size 0 (72 00 021) (R834) 220 €	for size 0 (72 00 022) (R811) 149 €	for size 0 (72 00 023) (R813) 260 €
for size 1 (72 01 020) (R863) 149 €	normal for size 1 (72 01 021) (R864) 220 €	for size 1 (72 01 022) (R865) 149 €	for size 1 (72 01 023) (R866) 260 €
for size 2 (72 02 020) (R807) 149 €	normal for size 2 (72 02 021) (R889) 220 €	for size 2 (72 02 022) (R813) 149 €	for size 2 (72 02 023) (R846) 260 €
	thin for size 1 (72 04 021) (R882) 220 €		
	thin for size 2 (72 03 021) (R826) 220 €		

Accessories			
Footplate height adjustable, can be folded up and swung away for transfer, recommended for models with Hi/Lo base. 		Ankle positioner – with socket lock, only for use together with the footplate 	
	for size 0	for size 1	for size 2
Lower leg length	23 - 30 cm	30 - 42 cm	38 - 51 cm
	Price / piece		Price / pair
for size 0 (72 00 029) (R847)	682 €		for size 0+1 (72 01 026) (R899)
for size 1 (72 01 029) (R848)	682 €		for size 2 (72 02 026) (R855)
for size 2 (72 02 029) (R857)	682 €		149 €

Accessories			
	Sandals with foot straps only for use together with the footplate		Wedges – only for use together with the footplate and sandals
	Price / pair		Price / pair
	for size 0 (72 00 027) (R944) 149 €		for size 0 (72 00 028) (R697) 118 €
	for size 1 (72 01 027) (R954) 149 €		for size 1 (72 01 028) (R698) 118 €
	for size 2 (72 02 027) (R964) 149 €		for size 2 (72 02 028) (R699) 118 €