

KS DRESSING TROLLEY - small



Product Code	45684
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279456844

Description

- DRESSING TROLLEY - small
Dressing trolley is supplied with upper dispenser (10 tilting bins).
- antimicrobial epoxy painted body with key lock
 - top surface with 4-sided rim
 - top surface with built-in compartment
 - ball-bearing full-extension sliding drawer with self-closing system
 - built-in handle
 - 4 swivelling castors Ø 125 mm - 2 with brake
 - lateral waste container 15 l with lid and catheter holder
 - size: small 760x714xh 1,022 mm

Made in Italy.