

KS DRESSING TROLLEY - standard



Product Code	45685
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279456851

Description

DRESSING TROLLEY - standard
Dressing trolley is supplied with upper dispenser (10 tilting bins).

- antimicrobial epoxy painted body with key lock
- top surface with 4-sided rim
- top surface with built-in compartment
- ball-bearing full-extension sliding drawer with self-closing system
- built-in handle
- 4 swivelling castors Ø 125 mm - 2 with brake
- lateral waste container 15 l with lid and catheter holder
- size: standard 911x704xh 1,075 mm

Made in Italy.