

Start M1

Lightweight wheelchair



Standard equipment

- Aluminium frame
- Tip-assist integrated into the frame
- Individual leg supports, angle-adjustable, plastic foot plates
- Individual calf straps
- Swing-away leg support holders
- Side panel standard with height-adjustable arm support, flip-up
- Seat and back support upholstery, black nylon
- Back support with lumbar curve
- Antibacterial push handles
- 24" drive wheels with PU tyres and quick release axels
- Aluminium handrim
- Knee lever wheel lock
- Frame colour silver metallic

Folding, max. user weight 125 kg/276 lbs

Factory settings:

- Back support height 42 cm (can be changed to 40 cm)
- Anterior seat height 49 cm, posterior seat height 47 cm
- Seat depth 42 cm (adjustable to 40 cm)
- Lower leg length 42 cm

Unit price: EUR X,XXX
(including side panels and leg supports)

All prices in EUR,
valid from 01/01/2022 to 31/12/2022

Provide estimate only **Order**

Billing address/customer number

Delivery address

Sales person

Name

.....
Date

Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Version	Individual order [_HUB]	20' container [_CT] (113 pcs.)	40' container [_CT] (235 pcs.)	40' HC container [_CT] (270 pcs.)
Version 1: 6" caster wheels with PU tyres	Price: xxx,-	Price: xxx,-	Price: xxx,-	Price: xxx,-
Dimensions:				
<input type="checkbox"/> SW 38 cm Art. no. 480F74=20101				
<input type="checkbox"/> SW 40.5 cm Art. no. 480F74=20102				
<input type="checkbox"/> SW 43 cm Art. no. 480F74=20103				
<input type="checkbox"/> SW 45.5 cm Art. no. 480F74=20104				
<input type="checkbox"/> SW 48 cm Art. no. 480F74=20105				
<input type="checkbox"/> SW 50.5 cm Art. no. 480F74=20106				

Version 2: 8" caster wheels (200 mm x 50 mm) with PU tyres	Price: xxx,-	Price: xxx,-	Price: xxx,-	Price: xxx,-
Dimensions:				
<input type="checkbox"/> SW 38 cm Art. no. 480F74=20201				
<input type="checkbox"/> SW 40.5 cm Art. no. 480F74=20202				
<input type="checkbox"/> SW 43 cm Art. no. 480F74=20203				
<input type="checkbox"/> SW 45.5 cm Art. no. 480F74=20204				
<input type="checkbox"/> SW 48 cm Art. no. 480F74=20205				
<input type="checkbox"/> SW 50.5 cm Art. no. 480F74=20206				

Version 3: with drum brake and PU tyres, 8" caster wheels (200 mm x 50 mm) with PU tyres	Price: xxx,-	Price: xxx,-	Price: xxx,-	Price: xxx,-
Dimensions:				
<input type="checkbox"/> SW 38 cm Art. no. 480F74=20301				
<input type="checkbox"/> SW 40.5 cm Art. no. 480F74=20302				
<input type="checkbox"/> SW 43 cm Art. no. 480F74=20303				
<input type="checkbox"/> SW 45.5 cm Art. no. 480F74=20304				
<input type="checkbox"/> SW 48 cm Art. no. 480F74=20305				
<input type="checkbox"/> SW 50.5 cm Art. no. 480F74=20306				

Version	Individual order [_HUB]	20' container [_CT] (113 pcs.)	40' container [_CT] (235 pcs.)	40' HC container [_CT] (270 pcs.)
Version 4 8" caster wheels (200 mm x 50 mm) with PU tyres, height-adjustable side panels, lap belt, standard seat cushion, anti-tipper on the left side, crutch holder on the right side, adaptable back support upholstery "Economy"	Price: xxx,-	Price: xxx,-	Price: xxx,-	Price: xxx,-
Dimensions:				
<input type="checkbox"/> SW 38 cm Art. no. 480F74=20401				
<input type="checkbox"/> SW 40.5 cm Art. no. 480F74=20402				
<input type="checkbox"/> SW 43 cm Art. no. 480F74=20403				
<input type="checkbox"/> SW 45.5 cm Art. no. 480F74=20404				
<input type="checkbox"/> SW 48 cm Art. no. 480F74=20405				
<input type="checkbox"/> SW 50.5 cm Art. no. 480F74=20406				

Version 5: with drum brake and PU tyres, 8" caster wheels (200 mm x 50 mm) with PU tyres, height-adjustable side panels, lap belt, standard seat cushion, anti-tipper on the left side, crutch holder on the right side, adaptable back support upholstery "Economy"	Price: xxx,-	Price: xxx,-	Price: xxx,-	Price: xxx,-
Dimensions:				
<input type="checkbox"/> SW 38 cm Art. no. 480F74=20501				
<input type="checkbox"/> SW 40.5 cm Art. no. 480F74=20502				
<input type="checkbox"/> SW 43 cm Art. no. 480F74=20503				
<input type="checkbox"/> SW 45.5 cm Art. no. 480F74=20504				
<input type="checkbox"/> SW 48 cm Art. no. 480F74=20505				
<input type="checkbox"/> SW 50.5 cm Art. no. 480F74=20506				

Accessories		Quantity per chair	Order quantity	Price
<input type="checkbox"/>	481A53=SK170_CT	Plug-on anti-tipper	1 piece	piece(s)
<input type="checkbox"/>	481A53=SK170_HUB			
<input type="checkbox"/>	481A74=SK065_CT	Crutch holder	1 piece	piece(s)
<input type="checkbox"/>	481A74=SK065_HUB			
<input type="checkbox"/>	481A74=SK025_CT	Transport wheels	2 pieces	piece(s)
<input type="checkbox"/>	481A74=SK025_HUB			
<input type="checkbox"/>	481H74=SK012_CT	Wheel lock lever extension	2 pieces	piece(s)
<input type="checkbox"/>	481H74=SK012_HUB			
<input type="checkbox"/>	481E74=PK101_CT	Side panel with height-adj. arm support, compl. (silver metallic)	1 pair	pair(s)
<input type="checkbox"/>	481E74=PK101_HUB			
<input type="checkbox"/>	481B74=PK162_CT	Elevating leg supports SW 38 – 40.5 cm	1 pair	pair(s)
<input type="checkbox"/>	481B74=PK162_HUB			
<input type="checkbox"/>	481B74=PK163_CT	Elevating leg supports SW 43 – 45.5 cm	1 pair	pair(s)
<input type="checkbox"/>	481B74=PK163_HUB			
<input type="checkbox"/>	481B74=PK164_CT	Elevating leg supports SW 48 cm	1 pair	pair(s)
<input type="checkbox"/>	481B74=PK164_HUB			
<input type="checkbox"/>	481B74=PK165_CT	Elevating leg supports SW 50.5 cm	1 pair	pair(s)
<input type="checkbox"/>	481B74=PK165_HUB			
<input type="checkbox"/>	481B74=PK075_CT	Amputation leg support SW 15" – 19"	1 pair	pair(s)
<input type="checkbox"/>	481B74=PK075_HUB			
<input type="checkbox"/>	481B74=PK076_CT	Amputation leg support SW 20"	1 pair	pair(s)
<input type="checkbox"/>	481B74=PK076_HUB			
<input type="checkbox"/>	476C00=SK1471_CT	Standard seat cushion 15" (cushion thickness 5 cm)	1 piece	piece(s)
<input type="checkbox"/>	476C00=SK1471_HUB			
<input type="checkbox"/>	476C00=SK1472_CT	Standard seat cushion 16" (cushion thickness 5 cm)	1 piece	piece(s)
<input type="checkbox"/>	476C00=SK1472_HUB			
<input type="checkbox"/>	476C00=SK1473_CT	Standard seat cushion 17" (cushion thickness 5 cm)	1 piece	piece(s)
<input type="checkbox"/>	476C00=SK1473_HUB			
<input type="checkbox"/>	476C00=SK1474_CT	Standard seat cushion 18" (cushion thickness 5 cm)	1 piece	piece(s)
<input type="checkbox"/>	476C00=SK1474_HUB			
<input type="checkbox"/>	476C00=SK1475_CT	Standard seat cushion 19" (cushion thickness 5 cm)	1 piece	piece(s)
<input type="checkbox"/>	476C00=SK1475_HUB			
<input type="checkbox"/>	476C00=SK1476_CT	Standard seat cushion 20" (cushion thickness 5 cm)	1 piece	piece(s)
<input type="checkbox"/>	476C00=SK1476_HUB			
<input type="checkbox"/>	481S74=SK010_CT	Lap belt with plastic buckle, 1100 mm, complete (SB 38 – 43 cm)	1 piece	piece(s)
<input type="checkbox"/>	481S74=SK010_HUB			
<input type="checkbox"/>	481S74=SK012_CT	Lap belt with plastic buckle, 1400 mm, complete (SB 45,5 – 50,5 cm)	1 piece	piece(s)
<input type="checkbox"/>	481S74=SK012_HUB			

Note: The accessories are included, not installed. They cannot be returned for credit.