

Start^{M4} XXL

Lightweight wheelchair

Article number: 480F53=40000_K



Standard equipment

- Aluminium tube frame with double crossbrace
- Tip-assist integrated into the frame
- Individual leg supports, angle-adjustable, plastic foot plates
- Swing-away leg support holders
- Side panel with arm support, height-adjustable
- Nylon seat and back support upholstery, black
- Modular design
- Drive wheels with reinforced spokes and quick-release axles
- Knee lever wheel lock
- Back support tube with lumbar curve, reinforced
- Different back support height setting options from 40 – 47.5 cm
- Weight: from 17.5 kg

Folding, max. user weight 160 kg/353 lbs

<p>Base price: EUR X,XXX (including side panels and leg supports)</p> <p>All prices in EUR valid from 01/01/2022 to 31/12/2022</p>	<p>Order no.:</p> <p>-----</p>
---	--------------------------------

<input type="checkbox"/> Provide estimate only <input type="checkbox"/> Order	
Billing address/customer number	
Delivery address	
Sales person	Name

	Date Signature
<p>By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.</p>	

Wheelchair dimensions	
Seat width (SW) _____ cm	Anterior seat height (ASH) _____ cm
Seat depth (SD) _____ cm	Posterior seat height (PSH) _____ cm
Back support height (BH) _____ cm	Lower leg length (LL) (frame measurement) _____ cm
Seat measurements in cm	
Seat width SW	<input type="checkbox"/> 43 <input type="checkbox"/> 45.5 <input type="checkbox"/> 48 <input type="checkbox"/> 50.5 <input type="checkbox"/> 53 <input type="checkbox"/> 55.5 <input type="checkbox"/> 58
Seat depth SD	<input type="checkbox"/> 43 <input type="checkbox"/> 46 <input type="checkbox"/> 49
Seat height measurements in cm (table on the next page)	
Anterior seat height	<input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51
Posterior seat height	<input type="checkbox"/> 38 <input type="checkbox"/> 41 <input type="checkbox"/> 44 <input type="checkbox"/> 47 <input type="checkbox"/> 50

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Information

Seat height table

Anterior seat height in cm


Caster fork, short

Adjustment range


Caster wheel diameter 5"	40 – 44
Caster wheel diameter 5.5"	42 – 44
Caster wheel diameter 6"	43 – 45
Caster wheel diameter 7"	45 – 47

Caster fork, long

Caster wheel diameter 5"	44 – 47
Caster wheel diameter 5.5"	44 – 48
Caster wheel diameter 6"	45 – 49
Caster wheel diameter 7"	47 – 50
Caster wheel diameter 8"	48 – 51




Horizontal drive wheel position



Seat height table

Posterior seat height in cm

Free mounting positions	Drive wheel size	
	22"	24"
1	38	41
2	41	44
3	44	47
4	47	50



The values indicated are rounded values which have been theoretically determined. Max. deviation of 1 cm.

The anterior seat height depends on the selected wheel size, caster fork and mounting position. Please refer to the seat height table!

Specified without seat cushion at 0° seat inclination.

The maximum possible seat inclination is 3.5 cm.

Frame accessories			
MA 50	<input type="checkbox"/>	Transport wheels	Add. charge
MA 57	<input type="checkbox"/>	Crutch holder, right	Add. charge
MA 58	<input type="checkbox"/>	Crutch holder, left	Add. charge
MA 61	<input type="checkbox"/>	Swinging anti-tipper, right	Add. charge
MA 62	<input type="checkbox"/>	Swinging anti-tipper, left	Add. charge

Caster wheels			(dimensions in mm)
MF 03	<input type="checkbox"/>	5" solid rubber	(130 x 33) -
MF 04	<input type="checkbox"/>	5.5" soft	(140 x 40) -
MF 06	<input type="checkbox"/>	6" solid rubber	(150 x 27) -
MF 07	<input type="checkbox"/>	7" pneumatic	(180 x 45) -
MF 08	<input type="checkbox"/>	7" solid rubber	(180 x 45) -
MF 09	<input type="checkbox"/>	8" pneumatic, narrow (long caster fork required)	(200 x 30) -
MF 10	<input type="checkbox"/>	8" solid rubber, narrow (long caster fork required)	(200 x 30) -
MF 12	<input type="checkbox"/>	8" pneumatic, wide (long caster fork required)	(200 x 50) -
MF 13	<input type="checkbox"/>	8" PU tyres, wide (long caster fork required)	(200 x 50) -

Caster forks		
MF 80	<input type="checkbox"/>	Caster fork, long -
MF 81	<input type="checkbox"/>	Caster fork, short -

Drive wheel size		
MG 01	<input type="checkbox"/>	22" -
MG 02	<input type="checkbox"/>	24" -

Drive wheel position		
<input type="checkbox"/> A	<input type="checkbox"/> B	<input type="checkbox"/> C

Drive wheel rim		
MG 05	<input type="checkbox"/>	Standard -
MG 80	<input type="checkbox"/>	Rear wheel with drum brake/attendant Add. charge

Tyres		
MG 30	<input type="checkbox"/>	Pneumatic tyres 1 3/8" -
MG 35	<input type="checkbox"/>	PU tyres, 1 3/8" -

Handrims		
MG 52	<input type="checkbox"/>	Handrim, aluminium, silver anodised, high -
MG 56	<input type="checkbox"/>	Handrim, polished stainless steel, high Add. charge

Drive wheel accessories		
MG 87	<input type="checkbox"/>	High-pressure air pump, 14 bar Add. charge
MG 88	<input type="checkbox"/>	Airman pump, battery-operated Add. charge
MG 84	<input type="checkbox"/>	Silicone handrim cover Add. charge

Wheel lock		
MH 01	<input type="checkbox"/>	Knee lever wheel lock, standard -
MH 06	<input type="checkbox"/>	Knee lever wheel lock for user and attendant Add. charge


Information



MA 50 Transport wheel



MA 62 Swinging anti-tipper



ME 05 Side panel, height-adjustable



MD 58 Stabiliser bar



MB 03 Single-panel leg support



MB 25 Elevating leg support



MB 72 Heel strap




MZ 53 Head support



MZ 54 Clinical installation set



MZ 55 Multi-axis installation set



MG 95 Spoke protector with decor

Wheel lock accessories		
MH 10	<input type="checkbox"/>	Plug-on wheel lock lever extension Add. charge
Back support height in cm		
<input type="checkbox"/> 40	<input type="checkbox"/> 42.5	<input type="checkbox"/> 45 <input type="checkbox"/> 47.5 -
Back support upholstery		
MD 04	<input type="checkbox"/>	Standard -
MD 05	<input type="checkbox"/>	Adaptable Add. charge
MD 06	<input type="checkbox"/>	Without Price reduction
Push handles		
MD 51	<input type="checkbox"/>	Standard -
MD 52	<input type="checkbox"/>	Telescoping Add. charge
Back support accessories		
MD 58	<input type="checkbox"/>	Stabiliser bar for push handles Add. charge
Side panels		
ME 01	<input type="checkbox"/>	Side panel, standard, 300 mm arm support, short mounted -
ME 01	<input type="checkbox"/>	Side panel, standard, 300 mm arm support, long mounted -
ME 06	<input type="checkbox"/>	Side panel, standard, 365 mm arm support, short mounted Add. charge
ME 06	<input type="checkbox"/>	Side panel, standard, 365 mm arm support, long mounted Add. charge
ME 03	<input type="checkbox"/>	Plug-on side panel -
ME 05	<input type="checkbox"/>	Side panel, height-adjustable, 300 mm arm support Add. charge
ME 10	<input type="checkbox"/>	Desk side panel, stepped -
Side panel accessories		
ME 60	<input type="checkbox"/>	Tray Add. charge
Seat upholstery		
MC 13	<input type="checkbox"/>	Upholstery, standard -
MC 06	<input type="checkbox"/>	Without Price reduction
Seat cushion		
MC 01	<input type="checkbox"/>	Seat cushion Standard (thickness 5 cm) Add. charge
MC 02	<input type="checkbox"/>	Cubic Foam seat cushion (thickness 5 cm, up to 150 kg) Add. charge
MC 05	<input type="checkbox"/>	Evolight seat cushion (thickness 6 cm, up to 125 kg) Add. charge
MC 29	<input type="checkbox"/>	Livo (seat cushion thickness 5 cm) Add. charge
MC 30	<input type="checkbox"/>	Evoflair (seat cushion thickness 8 cm) Add. charge
Lower leg length in cm (without seat cushion)		
Lower leg length (LL) in cm		(available adjustment range 27 – 47 cm) -
Leg support bar options		
MB 40	<input type="checkbox"/>	Short (27 – 40 cm) -
MB 41	<input type="checkbox"/>	Medium (32 – 47 cm) -
MB 42	<input type="checkbox"/>	Long (42 – 47 cm) -
Please note that the lower leg lengths vary in case of elevating leg supports.		

Leg supports		
MB 10	<input type="checkbox"/>	Individual leg supports, angle-adjustable, plastic foot plates -
MB 03	<input type="checkbox"/>	Single-panel leg support, angle-adjustable, aluminium Add. charge
MB 25	<input type="checkbox"/>	Elevating leg supports, without automatic length adjustment Add. charge (max. user weight 125 kg), pair
MB 26	<input type="checkbox"/>	Elevating leg support, without automatic length adjustment, right Add. charge (plus standard leg support, left) (max user weight 125 kg)
MB 27	<input type="checkbox"/>	Elevating leg support, without automatic length adjustment, left Add. charge (plus standard leg support, right) (max user weight 125 kg)
MB 21	<input type="checkbox"/>	Transtibial residual limb support, right (plus standard leg support, left) Add. charge (max. user weight 125 kg)
MB 22	<input type="checkbox"/>	Transtibial residual limb support, left (plus standard leg support, right) Add. charge (max. user weight 125 kg)
MB 23	<input type="checkbox"/>	Transtibial residual limb supports (pair) (max. user weight 125 kg) Add. charge

Leg supports accessories		
MB 72	<input type="checkbox"/>	Heel strap for individual leg supports Add. charge

Other accessories		
MZ 51	<input type="checkbox"/>	Lap belt with metal buckle Add. charge
MZ 52	<input type="checkbox"/>	Tool kit Add. charge
MZ 53	<input type="checkbox"/>	Head support (max. user weight 125 kg) Add. charge
MZ 54	<input type="checkbox"/>	Clinical installation kit for head support (max. user weight 125 kg) Add. charge
MZ 55	<input type="checkbox"/>	Head support installation kit, multi-axis, offset (max. user weight 125 kg) Add. charge

Note: Please also use the order form for wheelchair accessories 647F661.

Spoke protector		
MG 90	<input type="checkbox"/>	Spoke protector, transparent Add. charge
MG 95	<input type="checkbox"/>	Spoke protector with decor Add. charge

Frame colours		
002	<input type="checkbox"/>	Anthracite metallic -
003	<input type="checkbox"/>	Jet black RAL 9005 -
004	<input type="checkbox"/>	Ice blue -
016	<input type="checkbox"/>	Silver metallic -
020	<input type="checkbox"/>	Saturn blue metallic -
023	<input type="checkbox"/>	Traffic white RAL 9016 -