

ORDER FORM - mechanical wheelchair**KURY**

Completed and delivered order form is binding order for the production of mechanical KURY wheelchair.
Please, fill out questionnaire legibly and carefully! In case of doubt, please contact your dealer or directly the manufacturer.

RETAIL PRICE	
standard equipped wheelchair (•)	

Product registration	
Serial number	
Other identification	

DEALER or USER	
Company	
User Identification	
Contact	address:
	country
e-mail:	phone:

SHIPPING		DATE of SHIPMENT:	
<input type="checkbox"/>	delivery provided by manufacturer	Product will be delivered by the manufacturer to the address provided by the customer.	shipping costs:
<input type="checkbox"/>	delivery provided by the customer	Shipping is organized by customer after the product is ready to deliver.	-

Form filled by:	date
	signature

Manufacturer only

Výrobek zkompletoval:	datum	Výstupní kontrola:	datum
	podpis		podpis

All measurements are in cm (if not specified by other way).

ZP_IRISX1_01/2024EN

Reset Form

Save Form

Print Form

THE BASIC USER PARAMETERS

MAX. LOAD CAPACITY 100 kg

(including: user's weight, accessories, baggage etc.)

weight

kg

height


cm

- trained active user (interior and exterior use)
- beginner or a passive user with the assistance of accompanying persons
- leg amputation above the knee
- the requirement for a wheelchair's going through the door 60 cm (► for seat width max. 38 cm)
- prognosis of changes in user's body parameters in the future with regard to the construction of the wheelchair (change in weight, height, health condition)

BASIC MEASUREMENT OF USER'S PARAMETERS

Parameters for SEATING

X	cm	parameter for seat depth <i>measurement: user's back > under (back) of knee</i>
Y	cm	<i>measurement: foot > under (back) of knee</i> INCLUDING cushion! (Make sure what shoes used by the user usually.)
P	cm	height of cushion (Consider of cushion squeeze.)
Z1	cm	Height of seat edge in FRONT Without cushion! Cannot be Z1 ≤ Y2!
Z2	cm	Height of seat edge in REAR Without cushion!




V Accompaniment handles height from ground - FIXED ! = Z2 + height of backrest's textile part (see section 5) + cca 5 cm

NASTAVENÍ TĚŽIŠTĚ - umístění sedačky vzhledem k ose zadního kola (nebude-li vyplněno, bude nastavena pozice č. D2)

↳ Změnu těžiště lze servisním zásahem změnit.




D



- position D1 - active** - recommended for active users with good stability and driving wheelchair control
- position D2 - neutral and stable** - suitable for stable drive control and for assistant person drive
- position D3 - neutral and passive** - suitable for stable drive control and for assistant person drive
- position D4 - very passive** - difficult lifting-up of the front part of wheelchair for user and accompanying person, greater strain on the user's arms and shoulders

NOTES FOR MANUFACTURER (Note the limited options for individual adjustment and modification of wheelchair parts. To be consult with manufacturer.)

line	select.	item, part, parameters	pricing
1 frame			
1.2	<input type="checkbox"/>	extended wheelchair frame (individual surcharge for seat depth 47 cm - see 2.2)	▼
1.5	<input type="checkbox"/>	light frame material aluminium 7020 T6	▼
2 seat			
2.1		seat width <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45	•
2.2		seat depth <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45 <input type="checkbox"/> 47 (surcharge - see 1.2)	•
2.3		seat height FRONT (parameter Z1) <input type="checkbox"/> 46 <input type="checkbox"/> 47 ► NOT for front wheels 150 mm (6") <input type="checkbox"/> 48 <input type="checkbox"/> 49	•
2.4		seat height REAR (parameter Z2) <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49	•
2.5		seat upholstery - fabric kortexin, color: <input type="checkbox"/> black <input type="checkbox"/> red <input type="checkbox"/> blue	•
2.7	<input type="checkbox"/>	seat upholstery with straps for multipoint adjustment of seating stiffness	▼
3 antidekubitus cushion			
3.1	<input type="checkbox"/>	antidecubitus cushion GEMINI SPIN - to order GEMINI SPIN cushion, please fill specialized order form	▼
3.2	<input type="checkbox"/>	antidecubitus cushion GEMINI ALPHA - memory foam, 5 cm height, 3D textile fabric, anti-slip bottom, black color, according to seat dimension	▼
4 backrest - TYPE			
4.2	<input type="checkbox"/>	adjustable textile with straps - fixing straps for backrest stiffness adjustment; DOBBY cover (→6.2) or Air-Flow cover (→6.3)	•
4.3	<input type="checkbox"/>	anatomic backrest GEMINI AZO - to order, please fill specialized order form	▼
4.4	<input type="checkbox"/>	anatomic backrest GEMINI TWIN - to order, please fill specialized order form GEMINI TWIN	▼
5 backrest FIXED - height			
		type of construction	upholstery height - take into account height of antidecubitus cushion
5.1	<input type="checkbox"/>	fixed straight <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45	•
5.2	<input type="checkbox"/>	fixed angled (8°) <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45	•
5 backrest ADJUSTABLE HANDLES - height			
		type of construction	upholstery height - take into account height of antidecubitus cushion
5.4	<input type="checkbox"/>	adjustable handles - straight <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45	▼
5.5	<input type="checkbox"/>	adjustable handles - angled (8°) <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45	▼
5.6	<input type="checkbox"/>	prolonged (higher) handle tube for adjust. handles	▼
6 backrest - UPHOLSTERY			
6.2	<input type="checkbox"/>	DOBBY cover - combination DOBBY padding black + kortexin fabric: <input type="checkbox"/> black <input type="checkbox"/> red <input type="checkbox"/> blue	•
6.3	<input type="checkbox"/>	AirFLOW cover breathable, clr. combination: black + → <input type="checkbox"/> green (1) <input type="checkbox"/> turquoise (2) <input type="checkbox"/> red (3) <input type="checkbox"/> yellow (4) <input type="checkbox"/> grey (5) <input type="checkbox"/> pink (6) <input type="checkbox"/> black (7) <input type="checkbox"/> dark blue (8) <input type="checkbox"/> beige (9) <input type="checkbox"/> ruby (10) <input type="checkbox"/> orange (11)	▼
6.4	<input type="checkbox"/>	backside pocket - 30 (w) x 20 (h) x 2 (d) cm	▼

line	select.	item, part, parameters	specification	pricing
8 armrests (tilting and removable) / mudguards (fixed)				
per pair				
8.7	<input type="checkbox"/>	E04 mudguard, fixed	 ABS plastic	•
8.8	<input type="checkbox"/>	E04 M mudguard, fixed	 ABS plastic handle 30 cm adjustable handle - height 18-26 cm	adjust for height: ▼
8.9	<input type="checkbox"/>	E04 KAR mudguard, fixed	 CARBON fibre	▼

9 frame color (sides & cross construction, forks)				
9.2	<input type="checkbox"/>	basic	<input type="checkbox"/> glossy white (1) <input type="checkbox"/> white matt (24) <input type="checkbox"/> yellow (2) <input type="checkbox"/> red (3) <input type="checkbox"/> pink (5) <input type="checkbox"/> violet light (6) <input type="checkbox"/> blue (9) <input type="checkbox"/> glossy black (13) <input type="checkbox"/> black matt (26) <input type="checkbox"/> NAVI blue (29)	•
9.3	<input type="checkbox"/>	metallic (special)	<input type="checkbox"/> violet M (15) <input type="checkbox"/> Nivea blue (16) <input type="checkbox"/> red M (17) <input type="checkbox"/> anthracit M (18) <input type="checkbox"/> silver M (20) <input type="checkbox"/> green M (21) <input type="checkbox"/> KHAKI (28) <input type="checkbox"/> SPRING green (31)	▼
9 color package - set of parts				
9.4	<input type="checkbox"/>	forks	<input type="checkbox"/> white <input type="checkbox"/> black <input type="checkbox"/> red <input type="checkbox"/> blue <input type="checkbox"/> lime	▼

10 footplate				
10.3	<input type="checkbox"/>	one-piece aluminum	angle adjustable, adjust to: <input type="checkbox"/> 90° <input type="checkbox"/> 100° <input type="checkbox"/> other: °	•
10.4	<input type="checkbox"/>	one-piece carbon fibre		▼

11 castor wheels & forks				
11.4	<input type="checkbox"/>	PU 100 (4")	plastic black rim, black full tyre (PU)	•
11.1	<input type="checkbox"/>	PU 125 (5")		•
11.2	<input type="checkbox"/>	PU 150 (6")		•
11.5	<input type="checkbox"/>	ALPU 100 x 35 (4")	aluminum rim, black full tyre (PU)	▼
11.6	<input type="checkbox"/>	ALPU 125 x 35 (5")		<input type="checkbox"/> silver <input type="checkbox"/> black <input type="checkbox"/> red <input type="checkbox"/> blue
11.7	<input type="checkbox"/>	ALPU 150 x 35 (6") Frog Legs	3-spokes aluminum rim, black full tyre (PU)	<input type="checkbox"/> silver <input type="checkbox"/> black
11.10	<input type="checkbox"/>	suspension forks Frog Legs Phase 2 (polymer+carbon fiber)		▼

line	select.	item, part, parameters	specification	pricing
12.1		wheel camber		
12.0	•	<input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 4° ▶ up to 80 kg weight of user and with mudguard E04 usage only!		•

FULLY EQUIPPED REAR WHEELS - cannot change configuration!

12 BASIC				
12.2	<input type="checkbox"/>	24" complete wheel TOP	rim: 24" ultralight black, 24 spokes tyres: pneu Schwalbe MPE black, 24 x 1" hand rim: black anodised	•
12 HIGHER EQUIPMENT / DESIGN WHEELS				
12.3	<input type="checkbox"/>	DYNAMIC COLOR	24 x 1"; Schwalbe MPE, hard anodised hand rim; black hub; colored spokes: spokes: <input type="checkbox"/> black <input type="checkbox"/> white <input type="checkbox"/> red <input type="checkbox"/> blue <input type="checkbox"/> lime	▼
12.4	<input type="checkbox"/>	CAMEL Black	24 x 1"; black rim, Schwalbe MPE tyre, hard anodised hand rim	▼
12.5	<input type="checkbox"/>	CAMEL Silver	24 x 1"; silver rim, Schwalbe MPE tyre, hard anodised hand rim	▼
12.6	<input type="checkbox"/>	EXTREME Light ONE	24 x 1", Schwalbe ONE tyre, black anod. hand rim, black rim, lightweight (1 080 g)	▼
12.7	<input type="checkbox"/>	COUNTRY offroad	24 x 1,95", offroad tyre, black anodised hand rim <input type="checkbox"/> drum brake ready hub	▼

REAR WHEELS TO ORDER - choose: rim + tyres + hand rim

12 rim					
12.8	<input type="checkbox"/>	24" standard	24" simple silver rim, 36 spokes	•	
12.9	<input type="checkbox"/>	24" standard - drum brake ready hub	24" simple silver rim, 36 spokes ▶ select ONLY if (drum) brakes for accompanist (13.4) are selected	▼	
12.10	<input type="checkbox"/>	24" ultralight black	24" black, double rim, black hub, 24 spokes ▶ if FULL TYRE > PU 12 (12.21) only	▼	
12.13	<input type="checkbox"/>	one-hand wheel control: <input type="checkbox"/> right <input type="checkbox"/> left	24" black rim, driving shaft, incl. drive hand rim (black)	▼	
12 tyres					
PNEUMATIC Tyres					
12.14	<input type="checkbox"/>	PNEUMATIC black (Schwalbe RightRun)	24 x 1	Schwalbe RightRun Active black	▼
12.15	<input type="checkbox"/>	PNEUMATIC wide black	24 x 1 3/8	Rubena V69 Hook, standard, wide	▼
12.16	<input type="checkbox"/>	PNEUMATIC Schwalbe MPE	24 x 1	Marathon Plus Evolution, puncture proof	▼
12.17	<input type="checkbox"/>	PNEUMATIC Schwalbe RR black-blue	24 x 1	RightRun	▼
12.18	<input type="checkbox"/>	PNEUMATIC Schwalbe RR black-grey	24 x 1	RightRun	▼
12.19	<input type="checkbox"/>	PNEUMATIC Schwalbe RR black-red	24 x 1	RightRun	▼
12.20	<input type="checkbox"/>	PNEUMATIC Schwalbe ONE ultralight&slim	24 x 1	Schwalbe ONE, extreme light tyre	▼
FULL Tyres					
12.21	<input type="checkbox"/>	FULL TYRE slim (PU 12) black	24 x 1		▼
12.22	<input type="checkbox"/>	FULL TYRE wide (PU 14) black	24 x 1 3/8	▶ NOT with ultralight black rim 12.10	▼
12 hand rim (push rim)					
12.23	<input type="checkbox"/>	aluminium silver (anodised)	<input type="checkbox"/> 10 mm <input type="checkbox"/> 20 mm	24"	▼
12.24	<input type="checkbox"/>	aluminium dark (hard anodised)	<input type="checkbox"/> 10 mm <input type="checkbox"/> 20 mm	24"	▼
12.25	<input type="checkbox"/>	rubberized	<input type="checkbox"/> 10 mm <input type="checkbox"/> 20 mm	24"	▼
12.26	<input type="checkbox"/>	stainless steel	<input type="checkbox"/> 10 mm <input type="checkbox"/> 20 mm	24"	▼
12.27	<input type="checkbox"/>	8 pins		24" ▶ with standard rim ONLY (12.8)	▼
12.44	<input type="checkbox"/>	12 pins		24" ▶ with standard rim ONLY (12.8)	▼
12.28	<input type="checkbox"/>	covered with leather	<input type="checkbox"/> 10 mm <input type="checkbox"/> 20 mm	24", leather type: thread color:	▼

REAR WHEELS TO ORDER - special

12				
12.30	<input type="checkbox"/>	SPINERGY	Specification/type:	▼
13 brakes				
13.1	<input type="checkbox"/>	STANDARD	straight brake lever	•
13.2	<input type="checkbox"/>	ACTIVE	forward curved brake lever ▶ Not for removable footrest! 1.4	▼
13.3	<input type="checkbox"/>	COMPACT	brakes	▼
13.4	<input type="checkbox"/>	brakes FOR ACCOMPANIST (ASSISTANT)	▶ with standard rim only 12.9 or COUNTRY offroad wheel 12.7	▼

line	select.	item	pricing
14 accessories and special equipment			
14.1	<input type="checkbox"/>	anti-tipper (1 pc)	▼
14.2	<input type="checkbox"/>	double leg support belt (front and back) - adjustable, w. 45 mm, textile	▼
14.3	<input type="checkbox"/>	carbon plug of front fork (pair)	▼
14.4	<input type="checkbox"/>	standard pump with clips	▼
14.5	<input type="checkbox"/>	pump with clips 8 ATM (greater pressure fill - Schwalbe MPE, Schwalbe RR etc.)	▼
14.6	<input type="checkbox"/>	hand compressor, wireless, charging battery	▼
14.7	<input type="checkbox"/>	removable pocket (under seat) (black fabric kortexin; buckle for locking of folding wheelchair included)	▼
14.9	<input type="checkbox"/>	silicone handrim covers 24" (pair) <input type="checkbox"/> black <input type="checkbox"/> blue <input type="checkbox"/> red <input type="checkbox"/> yellow <input type="checkbox"/> pink <input type="checkbox"/> green <input type="checkbox"/> violet	▼
14.11	<input type="checkbox"/>	rear wheel spoke guards BASIC (type A) 24" (pair) <input type="checkbox"/> transparent (1) <input type="checkbox"/> black (2)	▼
14.12	<input type="checkbox"/>	rear wheel spoke guards DESIGNED (type B) 24" (pair) <input type="checkbox"/> rings (1) <input type="checkbox"/> night (2) <input type="checkbox"/> eyes (3) <input type="checkbox"/> wheel (4) <input type="checkbox"/> flower I (5) <input type="checkbox"/> flower II (6) <input type="checkbox"/> cascade red-blue (7) <input type="checkbox"/> cascade turquoise (8) <input type="checkbox"/> violet spiral (9) <input type="checkbox"/> blue spiral (10)	▼
14.13	<input type="checkbox"/>	rear wheel spoke guards DESIGNED (type C) 24" (pair) <input type="checkbox"/> floral (1) <input type="checkbox"/> skulls (2) <input type="checkbox"/> strip (3) <input type="checkbox"/> children (4) <input type="checkbox"/> letters (5) <input type="checkbox"/> compass (6) <input type="checkbox"/> spectrum (7) <input type="checkbox"/> flames (8) <input type="checkbox"/> gold wheel (9) <input type="checkbox"/> balloons (10) <input type="checkbox"/> blue sparkle (11)	▼
14.14	<input type="checkbox"/>	rear wheel spoke guards OWN-DESIGNED 24" (pair) - user's own picture or photo will be printed onto rear wheel spoke guards	▼
14.15	<input type="checkbox"/>	original backrest embroidery - price will be calculated according to the range of work	▼
14.26	<input type="checkbox"/>	protective textile covers for front frame part	▼
15 fixing and positioning accessories Addigy line			
15.1	<input type="checkbox"/>	AD 10 Abdominal fixing belt low universal size: 8 x 48 cm	▼
15.2	<input type="checkbox"/>	AD 11 Abdominal fixing belt for chair high <input type="checkbox"/> size 0 (12,5 x 40 cm) <input type="checkbox"/> size 1 (16,5 x 56 cm)	▼
15.3	<input type="checkbox"/>	AD 11 E Abdominal fixing belt for chair, narrow version upholstered part: 12,5 x 56 cm	▼
15.4	<input type="checkbox"/>	AD 12 Perineal fixing belt for chair <input type="checkbox"/> size 0 (12,5 x 40 cm) <input type="checkbox"/> size 1 (16,5 x 56 cm)	▼
15.5	<input type="checkbox"/>	AD 13 Torso fixing belt <input type="checkbox"/> size 0 (12,5 x 40 cm) <input type="checkbox"/> size 1 (16,5 x 56 cm)	▼
15.6	<input type="checkbox"/>	AD 14 Torso and pelvis fixing belt <input type="checkbox"/> size 0 (12,5 x 40 cm) <input type="checkbox"/> size 1 (16,5 x 56 cm)	▼
15.7	<input type="checkbox"/>	AD 18 Standard vest for chair <input type="checkbox"/> size 0 <input type="checkbox"/> size 1	▼
15.8	<input type="checkbox"/>	ADX 1 (pár) Foot support for wheelchair (appl. to footplate) <input type="checkbox"/> size 0 <input type="checkbox"/> size 1	▼
15.9	<input type="checkbox"/>	ADX 2 (pár) Foot belt - ankles strips <input type="checkbox"/> size 0 <input type="checkbox"/> size 1	▼
15.10	<input type="checkbox"/>	ADX 18 Universal elastic band (neoprene) size: 7,5 x 150 cm	▼
15.11	<input type="checkbox"/>	ADX 20 Pelvic/abdominal fixing belt <input type="checkbox"/> size 0 width 18-28 <input type="checkbox"/> size 1 width 23-38	▼
15.12	<input type="checkbox"/>	ADX 22 / <input type="checkbox"/> ADX 22C Trunk Support with zip <input type="checkbox"/> size 0 shoulders 33-43 <input type="checkbox"/> size 1 shoulders 43-53	▼
15.13	<input type="checkbox"/>	ADX 23 / <input type="checkbox"/> ADX 23C Shoulder Harness <input type="checkbox"/> size 0 shoulders 33-43 <input type="checkbox"/> size 1 shoulders 43-53	▼
15.14	<input type="checkbox"/>	ADX 30 Abduction support <input type="checkbox"/> size 0 20-35 cm <input type="checkbox"/> size 1 30-45	▼
15.15	<input type="checkbox"/>	ADX 4PE 4 point pelvic fixing belt <input type="checkbox"/> size 0 w. thighs 18-28 <input type="checkbox"/> size 0 w. thighs 23-38	▼