

<p><b>Base price: EUR X,XXX</b> (including side panels and leg support)</p> <p>All prices in EUR, valid from 01.01.2022 to 31.12.2022</p>	<p>Order no.:</p> <p>_____</p>
---	--------------------------------

## Zenit R

### Active wheelchair

Article number: 480F160=1\_AA04\_C



<input type="checkbox"/> Provide estimate only		<input type="checkbox"/> Order	
<input type="checkbox"/> Attachment for custom fabrication enclosed			
Billing address/customer number			
Delivery address			
Sales person		Name	
		.....	
		Date	Signature
<p>By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products, and you confirm to have read, understood and accepted our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.</p>			


#### Standard equipment

- Open aluminium tube frame
- Choice of front frame offset
- Back support upholstery, adaptable
- Aluminium side panels
- Modular design
- Drive wheels with quick-release axle
- Knee lever wheel lock 4.0
- Calf strap

Max. user weight 120 kg  
Crash-tested according to DIN ISO 7176-19

#### Wheelchair measurements\*

- 1 Anterior seat height (ASH)**  
= ground to top of seat upholstery
- 2 Lower leg length (LL)**  
= top of seat upholstery right/left to top of foot plate centre (without seat cushion)
- 3 Posterior seat height (PSH)**  
= ground to top of seat upholstery
- 4 Back support height (BH)**  
= top of seat tube (crossbrace) to top of back support upholstery
- 5 Seat depth (SD)**  
= front of seat upholstery to back support tube
- 6 Seat width (SW)**  
= outside edge of frame tube to outside edge of frame tube



A manufacturing tolerance of 10 mm applies to all measurements and of 10% to specified weights.

**The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.**

<b>Connection of power add-on drive</b> (only valuet for Germany and Switzerland)			
MZ 20	<input type="checkbox"/>	Mounting set Klaxon 6 km/h	Add. charge
MZ 21	<input type="checkbox"/>	Mounting set Klaxon 10 km/h	Add. charge
MZ 22	<input type="checkbox"/>	Mounting set Klaxon 15 km/h	Add. charge
<b>Road approval for electric drives &gt; 6 km/h</b> (only valuet for Germany)			
MZ 08	<input type="checkbox"/>	Individual approval TÜV	Add. charge
<b>Equipment</b>			
MA 110	<input checked="" type="checkbox"/>	Standard	
<b>Frame tube material</b>			
MA 08	<input type="checkbox"/>	Aluminium	-
MA 12	<input type="checkbox"/>	Carbon (available only with 85° anterior frame angle and 2 cm offset/abduction) (max. load capacity 100 kg)	Add. charge
<b>Anterior frame angle</b>			
MA 20	<input type="checkbox"/>	75°	-
MA 21	<input type="checkbox"/>	85°	-
<b>Front frame offset</b>			
MA 14	<input type="checkbox"/>	0 cm offset	-
MA 15	<input type="checkbox"/>	2 cm offset	-
MA 30	<input type="checkbox"/>	2 cm abduction (possible only up to SW 42 cm)	-
<b>Front frame extension</b>			
MA 23	<input type="checkbox"/>	3 cm	Add. charge
<b>Seat width SW in cm</b>			
<input type="checkbox"/> 32	<input type="checkbox"/> 34	<input type="checkbox"/> 36	<input type="checkbox"/> 38
<input type="checkbox"/> 40	<input type="checkbox"/> 42	<input type="checkbox"/> 44	<input type="checkbox"/> 46
<b>Seat depth SD in cm</b>			
<input type="checkbox"/> 36	<input type="checkbox"/> 38	<input type="checkbox"/> 40	<input type="checkbox"/> 42
<input type="checkbox"/> 44	<input type="checkbox"/> 46	<input type="checkbox"/> 48	<input type="checkbox"/> 50
<b>Anterior seat height (ASH) in cm</b>			
<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41
<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45
<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50
<input type="checkbox"/> 51	<input type="checkbox"/> 52	<input type="checkbox"/> 53	<input type="checkbox"/> 54
<input type="checkbox"/> 46			
<b>Frame accessories</b>			
MA 51	<input type="checkbox"/>	Tip-assist, right	+ 270 g Add. charge
MA 52	<input type="checkbox"/>	Tip-assist, left	+ 270 g Add. charge
MA 55	<input type="checkbox"/>	Anti-tipper, right	+ 517 g Add. charge
MA 56	<input type="checkbox"/>	Anti-tipper, left	+ 42 g Add. charge
MA 60	<input type="checkbox"/>	Frame pad	+ 42 g Add. charge
<b>Caster wheels</b>			
MF 30	<input type="checkbox"/>	4" solid rubber, black	-
MF 35	<input type="checkbox"/>	4" soft caster, aluminium	Add. charge
MF 31	<input type="checkbox"/>	5" solid rubber, black	-
MF 36	<input type="checkbox"/>	5" soft caster, aluminium	Add. charge
MF 04	<input type="checkbox"/>	5.5" soft caster, grey	-
MF 05	<input type="checkbox"/>	6" pneumatic tyre	-
MF 32	<input type="checkbox"/>	6" solid rubber, black	-
MF 37	<input type="checkbox"/>	6" soft caster, aluminium	Add. charge

<b>Caster wheel axle</b>						
MF 90	<input checked="" type="checkbox"/>	Forks with threaded axle				
<b>Posterior seat height (PSH) in cm</b>						
<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43
<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50
<b>Drive wheel size</b>						
MG 02	<input type="checkbox"/>	24"				-
MG 00	<input type="checkbox"/>	25"				-
Drive wheel camber <input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 3° <input type="checkbox"/> 4°						
<b>Drive wheel rim</b>						
			24"	25"		
MG 06	<input type="checkbox"/>	Hollow rim	•	•		-
MG 10	<input type="checkbox"/>	Spinergy SPOX	- 360 g	•		Add. charge
MG 11	<input type="checkbox"/>	Spinergy LX	- 400 g	•		Add. charge
MG 12	<input type="checkbox"/>	Spinergy LXL Blade, black spokes	- 660 g	•		Add. charge
MG 12	<input type="checkbox"/>	Spinergy LXL Blade, white spokes	- 660 g	•		Add. charge
MG 15	<input type="checkbox"/>	Infinity Ultralight (not possible in combination with MG 34)	- 720 g	•	•	Add. charge
MG 150	<input type="checkbox"/>	Infinity Ultralight, coloured (not possible in combination with MG 34)	- 720 g	•		Add. charge
		<b>Infinity Ultralight/caster wheel fork colour package</b>		24"	25"	
	<input type="checkbox"/>	Sparkle light orange	•			Add. charge
	<input type="checkbox"/>	Marine blue	•			Add. charge
	<input type="checkbox"/>	Candy red	•			Add. charge
	<input type="checkbox"/>	Sparkle granny-smith	•			Add. charge
<b>Tyres</b>						
			24"	25"		
MG 31	<input type="checkbox"/>	Schwalbe RightRun 1"	•	•		-
MG 34	<input type="checkbox"/>	PU 1" (not possible in combination with MG 15/MG 150)	+ 344 g	•		-
MG 341	<input type="checkbox"/>	PU 1", black (not possible in combination with MG 15/MG 150)	+ 344 g	•		Add. charge
MG 37	<input type="checkbox"/>	Schwalbe Marathon Plus Evolution 1"	+ 374 g	•	•	Add. charge
MG 371	<input type="checkbox"/>	Schwalbe Marathon Plus Evolution 1.1"	+ 379 g	•		Add. charge
MG 200	<input type="checkbox"/>	Schwalbe ONE 1"	- 226 g	•		Add. charge
<b>Handrims</b>						
			24"	25"		
MG 50	<input type="checkbox"/>	Aluminium, black coated, high	•	•		-
MG 52	<input type="checkbox"/>	Aluminium, silver anodised, high	- 40 g	•	•	Add. charge
MG 53	<input type="checkbox"/>	Aluminium, silver anodised, low	- 120 g	•		Add. charge
MG 54	<input type="checkbox"/>	Rubber-coated, high	+ 500 g	•	•	Add. charge
MG 56	<input type="checkbox"/>	Stainless steel, polished, high	+ 620 g	•	•	Add. charge
MG 57	<input type="checkbox"/>	Stainless steel, polished, low	+ 580 g	•		Add. charge
MG 63	<input type="checkbox"/>	CURVE L (anodised)	+ 200 g	•	•	Add. charge
MG 64	<input type="checkbox"/>	CURVE L Grip for users with limited gripping function	+ 840 g	•	•	Add. charge
MG 65	<input type="checkbox"/>	GEKKO, ergonomic handrim	+ 580 g	•	•	Add. charge

<b>Handrim mounting</b>			
<b>MG 70</b>	<input type="checkbox"/>	<b>Mounting position narrow</b> (15 mm)	-
<b>MG 71</b>	<input type="checkbox"/>	<b>Mounting position wide</b> (25 mm)	-
<b>Tipping point in cm</b> (measured parallel to the ground)			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Alternative measuring method</b> (specification of the horizontal drive wheel position of the Avantgarde)			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
(Tick the drive wheel position of your Avantgarde, we calculate the tipping point)			
<b>Quick-release axle system</b>			
<b>MG 27</b>	<input type="checkbox"/>	<b>Quick-release axle</b>	-
<b>MG 28</b>	<input type="checkbox"/>	<b>Quick-release axle with Tetra release</b>	+ 158 g <b>Add. charge</b>
<b>Drive wheel accessories</b>			
<b>MG 87</b>	<input type="checkbox"/>	<b>High-pressure air pump, 14 bar</b>	<b>Add. charge</b>
<b>MG 88</b>	<input type="checkbox"/>	<b>Airman pump</b>	<b>Add. charge</b>
<b>Wheel locks</b>			
<b>MH 01</b>	<input type="checkbox"/>	<b>Knee lever wheel lock 4.0</b>	-
<b>MH 30</b>	<input type="checkbox"/>	<b>Knee lever wheel lock 4.0, long wheel lock lever</b>	<b>Add. charge</b>
<b>MH 05</b>	<input type="checkbox"/>	<b>Scissor wheel lock</b>	- 142 g <b>Add. charge</b>
<b>Wheel lock accessories</b>			
<b>MH 10</b>	<input type="checkbox"/>	<b>Wheel lock lever extension, fixed</b>	+ 70 g <b>Add. charge</b>
<b>MH 21</b>	<input type="checkbox"/>	<b>Wheel lock lever extension, pluggable</b>	+ 34 g <b>Add. charge</b>
<b>Back support height in cm</b>			
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Back support system</b>			
<b>MD 01</b>	<input type="checkbox"/>	<b>Back support, straight</b> (measured 90° to the seat tube)	-
<b>MD 11</b>	<input type="checkbox"/>	<b>Angle-adjustable, fold-down onto seat tubes</b> (-10°/-5°/0°/+5°/+10°)	+ 240 g <b>Add. charge</b>
<b>MD 13</b>	<input type="checkbox"/>	<b>Back support angle according to customer request</b> (measured 90° to the ground) <input type="checkbox"/> 0° <input type="checkbox"/> -5° <input type="checkbox"/> -10°	<b>Add. charge</b>
<b>Back support upholstery</b>			
<b>MD 03</b>	<input type="checkbox"/>	<b>Adaptable, padded</b>	<b>Add. charge</b>
<b>MD 05</b>	<input type="checkbox"/>	<b>Adaptable, breathable</b>	-
<b>MD 18</b>	<input type="checkbox"/>	<b>Ultra light</b>	<b>Add. charge</b>
<b>MD 06</b>	<input type="checkbox"/>	<b>Without</b>	-
<b>Back support pad widening</b>			
<b>MD 80</b>	<input type="checkbox"/>	<b>+ 2 cm</b>	<b>Add. charge</b>
<b>MD 81</b>	<input type="checkbox"/>	<b>+ 4 cm</b>	<b>Add. charge</b>
<b>MD 82</b>	<input type="checkbox"/>	<b>Without</b>	-
<b>Use of BAXX back support system</b> (include order form, model-specific restrictions are possible)			
<b>MD 44</b>	<input type="checkbox"/>	<b>BAXX aluminium back support</b> (not mounted)	<b>Add. charge</b>


Push handles			
MD 50	<input type="checkbox"/>	Short	-
MD 51	<input type="checkbox"/>	Long	+ 128 g -
MD 55	<input type="checkbox"/>	Height-adjustable/removable	+ 580 g Add. charge
MD 56	<input type="checkbox"/>	Fold-down	+ 108 g Add. charge
MD 59	<input type="checkbox"/>	Without	- 258 g -

Side panels			
ME 02	<input type="checkbox"/>	Aluminium, with splash guard and protection against cold	-
ME 04	<input type="checkbox"/>	Carbon, with splash guard	Add. charge
ME 05	<input type="checkbox"/>	With arm support, short	Add. charge
ME 06	<input type="checkbox"/>	With arm support, long	Add. charge

Seat upholstery			
MC 15	<input type="checkbox"/>	Upholstery with storage compartment	Add. charge
MC 16	<input type="checkbox"/>	Upholstery, readjustable	-
MC 18	<input type="checkbox"/>	Ultra light	Add. charge
MC 06	<input type="checkbox"/>	Without	-

Seat cushions			
MC 01	<input type="checkbox"/>	Standard (thickness of cushion 5 cm)	Add. charge
MC 02	<input type="checkbox"/>	Cubic Foam (thickness of cushion 5 cm)	Add. charge
MC 05	<input type="checkbox"/>	EvoLight (thickness of cushion 6 cm)	Add. charge
MC 20	<input type="checkbox"/>	Terra seat (thickness of cushion 10 cm)	Add. charge
MC 21	<input type="checkbox"/>	Terra Aquos (thickness of cushion 10 cm)	Add. charge
MC 22	<input type="checkbox"/>	Terra Flair (thickness of cushion 10 cm)	Add. charge
MC 23	<input type="checkbox"/>	Cloud (thickness of cushion 10 cm)	Add. charge
MC 29	<input type="checkbox"/>	Livo (thickness of cushion 5 cm)	Add. charge
MC 30	<input type="checkbox"/>	EvoFlair (thickness of cushion 8 cm)	Add. charge

Lower leg length in cm (without seat cushion)		
Lower leg length (LL) in cm ..... (available setting range 28 cm – 51 cm)		
MB 130	<input checked="" type="checkbox"/>	Adjustable



Leg supports			
MB 160	<input type="checkbox"/>	Single-panel, aluminium* (LL 28 – 51 cm)	Add. charge
MB 180	<input type="checkbox"/>	Single-panel, carbon* (LL 28 – 51 cm)	Add. charge
MB 06	<input type="checkbox"/>	Aluminium light, fixed (LL 31 – 51 cm)	-

\* LL can be selected in dependence of wheel size, seat width and anterior seat height.


Other accessories			
MZ 51	<input type="checkbox"/>	Lap belt with metal buckle	Add. charge
MZ 59	<input type="checkbox"/>	Pocket for mobile phone	Add. charge

<b>Spoke protectors</b>		
<b>MG 90</b>	<input type="checkbox"/>	<b>Transparent</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MG 100</b>	<input type="checkbox"/>	<b>Graphics decor</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MG 101</b>	<input type="checkbox"/>	<b>Miracle decor</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MG 102</b>	<input type="checkbox"/>	<b>Network decor</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MG 103</b>	<input type="checkbox"/>	<b>Robots decor</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MG 104</b>	<input type="checkbox"/>	<b>Spectrum decor</b> <span style="float: right;"><b>Add. charge</b></span>

<b>Frame colours</b>		
<b>001</b>	<input type="checkbox"/>	<b>Brilliant silver</b> <span style="float: right;"><b>Add. charge</b></span>
<b>002</b>	<input type="checkbox"/>	<b>Anthracite metallic</b> <span style="float: right;"><b>Add. charge</b></span>
<b>003</b>	<input type="checkbox"/>	<b>Jet black RAL 9005</b> <span style="float: right;"><b>Add. charge</b></span>
<b>004</b>	<input type="checkbox"/>	<b>Ice blue</b> <span style="float: right;"><b>Add. charge</b></span>
<b>005</b>	<input type="checkbox"/>	<b>Sparkle light orange</b> <span style="float: right;"><b>Add. charge</b></span>
<b>009</b>	<input type="checkbox"/>	<b>Vinho sparkle</b> <span style="float: right;"><b>Add. charge</b></span>
<b>010</b>	<input type="checkbox"/>	<b>Marine blue</b> <span style="float: right;"><b>Add. charge</b></span>
<b>012</b>	<input type="checkbox"/>	<b>Sky blue</b> <span style="float: right;"><b>Add. charge</b></span>
<b>048</b>	<input type="checkbox"/>	<b>Jet black matte (texture)</b> <span style="float: right;"><b>Add. charge</b></span>
<b>014</b>	<input type="checkbox"/>	<b>Candy red</b> <span style="float: right;"><b>Add. charge</b></span>
<b>021</b>	<input type="checkbox"/>	<b>Sparkle granny-smith</b> <span style="float: right;"><b>Add. charge</b></span>
<b>023</b>	<input type="checkbox"/>	<b>White RAL 9016</b> <span style="float: right;"><b>Add. charge</b></span>
<b>027</b>	<input type="checkbox"/>	<b>Shimano matte</b> <span style="float: right;"><b>Add. charge</b></span>
<b>039</b>	<input type="checkbox"/>	<b>Black, matte anodised</b> <span style="float: right;">-</span>
<b>040</b>	<input type="checkbox"/>	<b>Anthracite, matte anodised</b> <span style="float: right;">-</span>
<b>041</b>	<input type="checkbox"/>	<b>Gold, matte anodised</b> <span style="float: right;">-</span>
<b>042</b>	<input type="checkbox"/>	<b>Blue, matte anodised</b> <span style="float: right;">-</span>
<b>043</b>	<input type="checkbox"/>	<b>Orange, matte anodised</b> <span style="float: right;">-</span>
<b>044</b>	<input type="checkbox"/>	<b>Green, matte anodised</b> <span style="float: right;">-</span>
<b>045</b>	<input type="checkbox"/>	<b>Violet, matte anodised</b> <span style="float: right;">-</span>
<b>046</b>	<input type="checkbox"/>	<b>Red, matte anodised</b> <span style="float: right;">-</span>

<b>Special accessories</b>		
<b>MZ 52</b>	<input type="checkbox"/>	<b>Tool kit</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MZ 58</b>	<input type="checkbox"/>	<b>Bag for drive wheels</b> (for 2 drive wheels up to 28") <span style="float: right;"><b>Add. charge</b></span>
<b>MZ 70</b>	<input type="checkbox"/>	<b>Kit for transportation in a motor vehicle</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MZ 89</b>	<input type="checkbox"/>	<b>Quick-release axle protection</b> <span style="float: right;"><b>Add. charge</b></span>
<b>MZ 92</b>	<input type="checkbox"/>	<b>Wheelblades</b> (Mini-skis for caster wheels) <span style="float: right;"><b>Add. charge</b></span>
<b>MZ 93</b>	<input type="checkbox"/>	<b>Backpack</b> <span style="float: right;"><b>Add. charge</b></span>

**Information**



**MG 06** Hollow rim



**MG 10** Spinerger SPOX



**MG 11** Spinerger LX



**MG 15** Infinity Ultralight



**MF 35/36/37** 4/5/6"  
soft caster



**MG 100** Graphic



**MG 103** Robots



**MG 104** Spectrum



**MF 90** Caster fork with  
threaded axle



**MH 01** Knee lever wheel  
lock 4.0



**MD 55** Push handles,  
height-adjustable/re-  
movable



**ME 04** Carbon side pan-  
el, with splash guard



**MD 03** Back support up-  
holstery, adaptable,  
padded



**MB 06** Aluminium light leg  
support, fixed



Drive wheel adapter, per-  
manently welded