

BI100 Hydraulic shower table - NEFTI (Made in ITALY)

Purpose

Devices intended to assist patients during washing/showering in a horizontal or inclined position (up to one maximum of 5 °), with the help of an assistant who manages the controls and performs the washing operations.

Main features

Shower stretcher for washing and transferring patients, hydraulically elevatable from 59 to 93 cm, with a tiltable spring up to 5° for proper drainage during and after use. Constructed with epoxi painted steel anti-corrosive treatment of automotive derivation to guarantee reliability and durability.

Padded cover in PVC reinforced with polyester, phthalate-free, sanitizable. High sliding wheels, Ø 125 mm, are equipped with a total block brake. The limited width of the stretcher allows for easy movement without the need for effort even in spaces. Equipped with 180 ° folding handrails to allow an easy positioning of the patient on the stretcher directly from the bed and vice versa.


Image product

REFERENCES

Fabricante en conformidad con 93/42	<i>MORETTI S.P.A.</i>
Παραγωγή	<i>ITALIA</i>
Clase de dispositivo 93/42	<i>Clase I</i>
GMDN Clasification Code	<i>31092</i>
Μέγιστο βάρος χρήστη	<i>180 Kg</i>
Ανύψωση	<i>Υδραυλική</i>

MATERIALES PRINCIPALES

Estructura	<i>Acero pintado</i>
Material	<i>PVC</i>

Dimensions

Dimensions

	<i>A</i>	<i>B</i>	<i>C</i>	<i>D</i>	<i>E</i>	<i>F</i>	<i>H</i>	<i>I</i>	<i>L</i>	<i>M</i>	<i>N</i>
BI100 (cm)	59-93	141	23	20	206	200	67	56	66	73	0-5°

Optional accessories

<i>Model</i>	<i>Dimensions</i>	<i>Materials</i>
BIA110 – electronic	55 x 55 x 2/14 Cm	<i>Polyester - external</i> <i>Polyurethane D.18 - interior</i>
BIA100 – plumber	55 x15 x7 Cm	<i>Polyster – externo</i> <i>Polyurethane D.18 - interior</i>