Etac Swift shower stool/chair

Swift is a prize-winning series of shower stools/chairs designed with functionality and simplicity in mind. The basic shower stool can easily be transformed to a complete shower chair, by adding arm and back supports.

Stable on any floor

Swift has a unique design that adapts to uneven floors thanks to the flexible seat. The soft high quality ferrules adapt to the floor and ensures a good grip. All these features contribute to the users safety.

Changes with the circumstances

Swift is adaptable to different needs. It is possible to start using the shower stool and add arm supports and back support if the circumstances changes. Swift is easy to adjust to the required height with the telescopic legs. It can be set one step lower at the front for a slightly forward tilted seat to facilitate standing up. Soft pads which are specially suited for showering are available should extra comfort be needed.



See the movie about Swift on www.etac.com

Easy to handle

Swift is easy to handle thanks to its light weight. It comes flat-packed and is assembled without the need for any tools. Swift is also corrosion free.



Without arm supports, the attachments holes function as firm grips.

The arm supports have a stable fixation to ensure stability.

The back support has a practical grip for moving the chair



Easy to handle
The flat-packed Swift is easy to handle.



Safe surface The structure of the surface ensures that Swift doesn't become slippery when wet.



Swift is also available in a low version. Seat height 34.5-42 cm (13.6-16.5")



Standard seat height 42-57 cm (16.5-22.4")

Accessories



Back support Green: 81701120 Grev: 84005069 Blue: 84005071



Soap dish Green: 81701140 Grey: 80209266 Blue: 81701340



Back support pad Green: 81701130 Grey: 84005072 The back support pad doesn't become slippery when wet and gives an instantly warm



feeling. Material: Cellular polyethylene. Cleaning at max. 60°C (140°F) Seat pad Green: 81701132

Grey: 84005073 The seat pad doesn't become slippery when wet and gives an instantly warm feeling. Material: Cellular polyethylene. Cleaning at max. 60°C (140°F)



Arm supports pair Green: 81701110 Grey: 84005074 Blue: 84005075

Description

·	Green	Grey	Blue
Swift shower stool	81701010	81701410	81701510
Swift shower chair with arm supports	81701030	81701430	81701530
Swift shower stool low	n/a	81701450	n/a
Swift shower chair low with arm supports	n/a	81701460	n/a

Dimensions

Seat height: 42-57 cm (161/2-221/2")

Seat height Swift low: 34.5-42 cm (13½-16½")

Total height with arm supports: $80-95 \text{ cm} (31\frac{1}{2}-37\frac{1}{2}")$

Total height Swift Low with arm supports: 72.5-80 cm (281/2-311/2")

Seat area: 54 x 41cm (21¼ - 16¼") Total width, stool: 54 cm (211/4")

Total width, with arm supports: 54.5 cm (21")

Depth: 50 cm (19¾") Depth Swift Low: 46.5 cm (18")

Width between arm supports: 45 cm (17¾")

Material

Polypropylene, aluminium

Weight

Shower stool: 3.1 kg (6.8 lbs) Shower stool low: 2.9 kg (6.4 lbs) Shower chair: 4.7 kg (10.4 lbs) Shower chair low: 4.5 kg (8.8 lbs)

Green

Colours





Grey











or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.