



Detailed technical specifications

Panthera Bambino 3



Rev.2021-01

Standard design

Chassis	Chrome molybdenum tubing blue or black
Rear axle	Carbon 6°
Balancing	Ungraded adjustable
Rear wheels	20", 22", 24"
Push rim	Titanium
Brakes	Bambino brake
Castors	120 mm
Footrest	Footplate, adjustable in height, angle, depth
Backrest	Ungraded adjustable angle
Backrest upholstery	Adjustable
Others	Toolkit, Manual

		Bambino 3 Short	Bambino 3 Short	Bambino 3	Bambino 3	Bambino 3	Bambino 3
	Model code	G518	G518	G517	G517	G517	G517
A	Seat width	24 cm	27 cm	24 cm	27 cm	30 cm	33 cm
B	Total width	55 cm	58 cm	55 cm	58 cm	61 cm	64 cm
C	Total length	64 cm	64 cm	69 cm	69 cm	69 cm	69 cm
D	Total height	57.5-65 cm	57.5-65 cm	57.5-65 cm	57.5-65 cm	57.5-65 cm	57.5-65 cm
E	Seat angle	4°	4°	4°	4°	4°	4°
F	Seat height rear	39 cm	39 cm	39 cm	39 cm	39 cm	39 cm
G	Seat height front	41 cm	41 cm	41 cm	41 cm	41 cm	41 cm

H	Seat depth	18-25 cm	18-25 cm	25-30 cm	25-30 cm	25-30 cm	25-30 cm
I	Backrest angle back-forwards	11.5° - (-7.5°)	11.5° - (-7.5°)	11.5° - (-7.5°)	11.5° - (-7.5°)	11.5° - (-7.5°)	11.5° - (-7.5°)
J	Backrest height	20-28 cm alt. 27-35 cm	20-28 cm alt. 27-35 cm	20-28 cm alt. 27-35 cm	20-28 cm alt. 27-35 cm	20-28 cm alt. 27-35 cm	20-28 cm alt. 27-35 cm
K	Footrest height (footrest above ground)	12.5-34.5 cm	12.5-34.5 cm	12.5-34.5 cm	12.5-34.5 cm	12.5-34.5 cm	12.5-34.5 cm
L	Footrest angle	100°	100°	100°	100°	100°	100°
M	Width	26.5 cm	29.5 cm	32.5 cm	35.5 cm	38.5 cm	41.5 cm
N	Length	55 cm	55 cm	60 cm	60 cm	60 cm	60 cm
O	Height	32 cm	32 cm	32 cm	32 cm	32 cm	32 cm
	Total weight	9764 g	9992 g	9764 g	9992 g	10220 g	10448 g
	Transport weight	5574 g	5800 g	5574 g	5800 g	6030 g	6258 g
	Max user weight	60 kg	60 kg	60 kg	60 kg	60 kg	60 kg

Article no. for configuration, not possible to order separately

Seat width	Article no. black	Article no. blue	Article no. no paint
Seat width 24 cm	5170024BLG	5170024G	5170024GP
Seat width 27 cm	5170027BLG	5170027G	5170027GP
Seat width 30 cm	5170030BLG	5170030G	5170030GP
Seat width 33 cm	5170033BLG	5170033G	5170033GP
Seat width short 24 cm	5180024BLG	5180024G	5180024GP
Seat width short 27 cm	5180027BLG	5180027G	5180027GP