


PORTABLE ROLLABLE WHEELCHAIR RAMP
MOBILE ROLLRAMPE ROLLSTUHLRAMPE
RAMPA IN ALLUMINIO AVVOLGIBILE - MOBILE
RAMPA ENROLLABLE PARA SILLAS DE RUEDAS
RAMPE ENROULABLE
РОЛЛ-ПАНДУС ДЛЯ ИНВАЛИДНЫХ КОЛЯСОК
휠체어 경사로
منحدر متنقل
巻取式スロープ

user manual
Bedienungsanleitung
manuale utente
manual de usuario
manuel
руководство пользователя
사용자 설명서
دليل المستخدم
取扱説明書





medlis

innovative ramp systems


Made in Turkey by Medlis Ramps




Ahi Evran Osb Mah. Özbekhanlı Sk. No 33, Sincan,
Ankara, Turkey T: +90 312 315 8795 info@medlis.com.tr




max



350kg



Rec. height : cm
















11kg/m

350kg



SUPPORT FEET CHART

Art No			Includes			Recommended ramp sections where support legs be mounted (Support feet are installed to each of the 6th or 7th ramp sections.)						
	Total Length		Section Qty	Access Profiles Qty	Support Feet Qty							
	cm	inch				1st feet	2nd feet	3rd feet	4th feet	5th feet	6th feet	7th feet
RR..80	80	31,5	3	2	-	-						
RR..100	100	39,4	4	2	-	-						
RR..120	120	47,2	5	2	-	-						
RR..140	140	55,1	6	2	-	-						
RR..160	160	63,0	7	2	-	-						
RR..180	180	70,9	8	2	-	-						
RR..200	200	78,7	9	2	-	-						
RR..220	220	86,6	10	2	-	-						
RR..240	240	94,5	11	2	1	5-6						
RR..260	260	102,4	12	2	1	5-6						
RR..280	280	110,2	13	2	1	6-7						
RR..300	300	118,1	14	2	1	6-7						
RR..320	320	126,0	15	2	1	7-8						
RR..340	340	133,9	16	2	1	7-8						
RR..360	360	141,7	17	2	2	5-6	11-12					
RR..380	380	149,6	18	2	2	6-7	11-12					
RR..400	400	157,5	19	2	2	6-7	12-13					
RR..420	420	165,4	20	2	2	6-7	12-13					
RR..440	440	173,2	21	2	2	6-7	12-13-14					
RR..460	460	181,1	22	2	2	7-8	13-14					
RR..480	480	189,0	23	2	3	6-7	12-13	18-19				
RR..500	500	196,9	24	2	3	6-7	12-13	18-19				
RR..520	520	204,7	25	2	3	6-7	12-13	18-19				
RR..540	540	212,6	26	2	3	6-7	12-13	18-19				
RR..560	560	220,5	27	2	4	5-6	11-12	16-17	22-23			
RR..580	580	228,3	28	2	4	5-6	11-12	16-17	22-23			
RR..600	600	236,2	29	2	4	5-6	11-12	16-17	22-23			
RR..620	620	244,1	30	2	4	6-7	12-13	18-19	23-24			
RR..640	640	252,0	31	2	3	6-7	12-13	18-19	24-25			
RR..660	660	259,8	32	2	4	6-7	12-13	18-19	24-25			
RR..680	680	267,7	33	2	5	5-6	11-12	17-18	22-23	27-28		
RR..700	700	275,6	34	2	5	5-6	11-12	17-18	22-23	27-28		
RR..720	720	283,5	35	2	5	6-7	12-13	18-19	23-24	28-29		
RR..740	740	291,3	36	2	5	6-7	12-13	18-19	24-25	29-30		
RR..760	760	299,2	37	2	5	6-7	12-13	18-19	24-25	30-31		
RR..780	780	307,1	38	2	5	6-7	12-13	18-19	24-25	30-31		
RR..800	800	315,0	39	2	6	5-6	11-12	17-18	22-23	27-28	33-34	
RR..820	820	322,8	40	2	6	5-6	11-12	17-18	22-23	27-28	33-34	
RR..840	840	330,7	41	2	6	6-7	12-13	18-19	24-25	30-31	34-35	
RR..860	860	338,6	42	2	6	6-7	12-13	18-19	24-25	30-31	35-36	
RR..880	880	346,5	43	2	7	5-6	11-12	16-17	21-22	26-27	31-32	37-38
RR..900	900	354,3	44	2	7	5-6	11-12	16-17	21-22	26-27	31-32	37-38
RR..920	920	362,2	45	2	7	6-7	12-13	17-18	23-24	28-29	33-34	38-39
RR..940	940	370,1	46	2	7	6-7	12-13	18-19	24-25	29-30	34-35	39-40
RR..960	960	378,0	47	2	7	6-7	12-13	18-19	24-25	30-31	36-37	41-42
RR..980	980	385,8	48	2	7	6-7	12-13	18-19	24-25	30-31	36-37	42-43
RR..1000	1000	393,7	49	2	7	6-7	12-13	18-19	24-25	30-31	36-37	42-43