

Model

Phoenix i Wheelchair –

Max User Weight – 100kg

Yes I am within max user weight.

CAD

Although, not a requirement it may assist when determining your wheelchair configuration to view a CAD drawing of your order prior to manufacture. You can also order a CAD drawing for quotation purposes.

- Quote only (NO CAD DRAWING) **£0.00**
- Quote with CAD drawing **£300.00 (2 customer amendments ONLY)**
- Wheelchair order with CAD drawing for approval **£250.00 (2 customer amendments ONLY)**


Frame Dimensions

Seat Sizing – Please refer to the diagram below when choosing your seat sizing.

Seat Width

Measurement from left to right seat edge. Max width on diagram (left).


- 14" (35.6cm)
- 15" (38.1cm)
- 16" (40.6cm)
- 17" (43.2cm)



Seat Depth

Measurement from front edge of seat to front of backrest


- 15" (38.1cm)
- 16" (40.6cm)
- 17" (43.2cm)
- 18" (45.7cm)










Control Side

Please choose which side you prefer your control dial on your Phoenix i frame.

- Right
- Left



Front Seat to Floor Height	Rear Seat to Floor Height	Wheel Camber
Measurement from floor to the front of seat.	Measurement from floor to the rear of seat.	Wheel camber is set at 2°
19" (48.3cm) 	17" (43.2cm) 	2° 
Centre of Gravity Measurement from front of backrest to centre of axle. This model has a fixed (non-adjustable) axle.		
<p>Due to the added weight on the front of the wheelchair from the power assist you should select a CG position 1" further forward from what you would do on a conventional wheelchair.</p> <ul style="list-style-type: none"> ▪ i.e. if you normally like an active setup with a +3 CG, selecting +4 on the Phoenix i will feel the same lightness and stability. ▪ If you normally have a stable +2 setting you should select a +3. <p>24" wheels are only viable with +3. If you are seeking an active setup we recommend a 25" or 26" wheel with +4</p> <div style="display: flex; align-items: center; justify-content: center;"> <div style="margin-right: 20px;"> <p>+3" <input type="checkbox"/></p> <p>+4" <input type="checkbox"/></p> </div>  </div>		
Leg Length – Fixed. This is the length from your wheelchair seat (Front edges of seat without cushion) to the back of footrest. The footrest is part of the single piece frame and is factory fixed for height at the time of manufacture. The leg length can be reduced via the adjustable foot plate but can't be extended. This measurement is the max length.	Folding Backrest Height Measurement from rear seat base (excluding seat cushion) to top of backrest. Adjustable angle.	
<p><i>Measure from the back of your knee to floor wearing shoes. Then remove thickness of seat cushion from this measurement to determine wheelchair leg length. Alternatively, you can measure your current wheelchair from the front edge of the seat canvas to the rear edge of the footrest.</i></p> <div style="display: flex; align-items: center; justify-content: center;"> <div style="margin-right: 20px;"> <p>14" (35.6cm) <input type="checkbox"/></p> <p>15" (38.1cm) <input type="checkbox"/></p> <p>16" (40.6cm) <input type="checkbox"/></p> </div>  </div>	<div style="display: flex; align-items: center; justify-content: center;"> <div style="margin-right: 20px;"> <p>14" (35.6cm) <input type="checkbox"/></p> <p>15" (38.1cm) <input type="checkbox"/></p> <p>16" (40.6cm) <input type="checkbox"/></p> <p>17" (43.2cm) <input type="checkbox"/></p> </div>  </div>	

Rear Wheels	
<p>Rear Wheel Size Please select wheel choice.</p> <p><input type="checkbox"/> Spinergy LX <i>All-purpose, daily use, lightweight, travel</i></p> <p><input type="checkbox"/> Spinergy XLX <i>Most rigid. Extra room to hold quick release.</i></p> <p><input type="checkbox"/> Spinergy CLX (24" and 25" only) <i>Lightweight, Aerodynamic, ultra durable carbon wheels</i></p> <p>Please select wheel size: 24" <input type="checkbox"/> (+3" COG only) 25" <input type="checkbox"/> (+3" & +4" COG) 26" <input type="checkbox"/> (+3" & +4" COG)</p> 	<p>Rear Wheel Colour Silver Hubs only available with Black Spokes.</p> <p>X2 Hub Colour</p> <p>Silver Hub <input type="checkbox"/></p> <ul style="list-style-type: none"> • Coloured Spokes with Silver Hub +£55 <p>Black Hub <input type="checkbox"/></p> <p>Spoke Colour</p> <p><input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> White <input type="checkbox"/> Orange <input type="checkbox"/> Green <input type="checkbox"/> Pink <input type="checkbox"/> Purple</p>
<p>Handrims</p> <p><input type="checkbox"/> Standard Anodized – Black</p> <p><input type="checkbox"/> Natural fit LT (<i>Ergonomic, Oval & Light</i>)</p> <p><input type="checkbox"/> Carbolife Quadro (<i>For tetraplegics without triceps function. Larger gripping surface</i>)</p> <p><input type="checkbox"/> Carbolife Gekko (<i>Ergonomic grip, Light, low friction sides reduces heat when braking</i>)</p>	<p>Tires</p> <p><input type="checkbox"/> Schwalbe Right Run (Black/Grey)</p> <p><input type="checkbox"/> Primo Solid HD (Solid tires) – not available with CLX wheels</p> <p><input type="checkbox"/> Schwalbe Marathon Plus Evo (most puncture proof)</p>

Additional Batteries and Chargers	
<p>Batteries All Phoenix i's come with a 4.0AH and 9.0AH battery. Here you can order additional batteries.</p> <p><input type="checkbox"/> 4.0 AH <input type="checkbox"/> Qty +£199.00</p> <p><input type="checkbox"/> 9.0 AH <input type="checkbox"/> Qty +£350.00</p>	<p>Chargers All Phoenix i's come with a standard Travel charger included. Upgrade to a Fast charger or add additional chargers below.</p> <p><input type="checkbox"/> Upgrade to Fast Charger +£100 Recommended. <i>Faster and better charge. Air cooled so battery doesn't heat up.</i></p> <p>Additional Chargers:</p> <p><input type="checkbox"/> Travel Charger +£150.00 (<i>quiet, great for travel, charge in 90 mins</i>)</p> <p><input type="checkbox"/> Fast Charger +£260.00 (<i>fan assisted, fast 40 minute charge</i>)</p>

Options and Accessories

<p>Spare Battery Holders Phoenix i's can have spare battery holders located under the seat. If you would like this feature please select below:</p> <p><input type="checkbox"/> No additional Battery Holders.</p> <p><input type="checkbox"/> 1 (<i>All seat sizes</i>)</p> <p><input type="checkbox"/> 2 (<i>15" seat width and above</i>)</p>	<p>Wheel Locks The Phoenix i comes with electronic downhill built in. For secure parking of your wheelchair we recommend adding manual wheel locks to your order.</p> <p><input type="checkbox"/> No Wheel Locks</p> <p><input type="checkbox"/> Scissor Locks</p>
<p>Lights Super bright, forward illuminating LED lights fitted to the left and right of your wheelchair.</p> <p><input type="checkbox"/> No LED Lights</p> <p><input type="checkbox"/> LED Lights</p>	<p>Anti Tipper An anti-tipper is a safety device that helps prevent a wheelchair from tipping backwards.</p> <p><input type="checkbox"/> No anti-tipper</p> <p><input type="checkbox"/> 1x Anti-Tipper – removable</p> <p><input type="checkbox"/> 2x Anti-Tipper - removable</p>
<p>Push Handles The Phoenix i is designed for an active independent user. The push handles do not extend beyond the height of the backrest. They are intended to be an occasional grab point and for occasional assistance</p> <p>Please choose one:</p> <p><input type="checkbox"/> No Push Handles</p> <p><input type="checkbox"/> Push Handle Ready frame <i>Frame prepared and reinforced to receive push handles at later date. Must be done at manufacture. Adding Push handles will then cost £400.</i></p> <p><input type="checkbox"/> Push Handles +£500 <i>Must be done at manufacture. Frame will be fitted with 2x quick release, detachable handles. Does not fold with backrest</i></p>	<p>Spares The control dial and casters wheels on the Phoenix i are "plug and play" and easy to replace at home should damage occur. We recommend having a spare of each so you can get going again quickly if damage occurs.</p> <p><input type="checkbox"/> Control Unit (<i>control dial</i>)</p> <p><input type="checkbox"/> Control Unit with LED light (<i>control dial – for customers with LED lights fitted to their Phoenix i</i>)</p> <p><input type="checkbox"/> Pair of Motorised Casters</p> <p><input type="checkbox"/> Essential Spares Kit 1 (<i>1x motorised caster wheel and 1x control dial – recommended</i>)</p> <p><input type="checkbox"/> Essential Spares Kit 2 (<i>1x motorised caster wheel and 1x control dial with LED Light Fitting. For customers with LED lights already fitted to their Phoenix i – recommended</i>)</p>

Confirmation

I confirm I understand that the maximum user weight for this chair is 100kg and user weight exceeding this invalidates warranty.

Signature:

Name:

Date:

If you have any questions, please drop us an email at info@phoenixinstinct.com or phone us at 01309 221020