

ORDER FORM / wheelchair

Completed and delivered sighting protocol is **binding order** for the production of mechanical KURY wheelchair. Please, fill out questionnaire legibly and carefully! In case of doubt, please contact your dealer or the manufacturer.

Product registration

Serial number

Date of product shipping

| Dealer or user | | |
|----------------|------------------|------------------|
| Company | | |
| Name | | |
| Contact | address | |
| | | |
| | | |
| | country | |
| | e-mail | phone |
| Other | user height (cm) | user weight (kg) |

| | | | |
|--|-----------|-----------------------------|-----------|
| User agrees with product parameters referred to in this order form questionnaire. According to these parameters, the product will be entered into production. | | Signature of user / dealer: | |
| Questionnaire completed: | date | Output control: | date |
| | signature | | signature |

- ☒ basic equipment included in price
☐ additional item with surcharge (see pricelist)


AS 07/2017 EN

Clear form

Save form

Print form

A PROPORTIONS


X: cmY: cm

- rear wheel position in the holder (change the center of gravity)
- default position 2


K: ☐ ☐ ☐

1 2 3


pict. 1 backrest FIXED




pict. 2 positioning (telescopic) handles



pict. 3 additional adjustable headrest



pict. 4 shortened handles



1 FRAME

- 1.1 ☐ basic wheelchair frame
- 1.2 ☐ extended wheelchair frame + 5 cm (for users of 185 cm high and more)

2 SEAT (Always textile made of cortex fabric, color see 3.3 according to type of backrest)

| | | |
|--|---|---|
| width (S parameter) | <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45 <input type="checkbox"/> 47 | ✓ |
| 2.1 depth (from X parameter) | <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 45 (extended frame) | ✓ |
| seat height from the ground (Z parameter) | <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 | ✓ |
| 2.2A <input type="checkbox"/> anti-decubitus anatomic pad | thickness <input type="checkbox"/> 3 cm <input type="checkbox"/> 5 cm | ○ |
| 2.2B | color: <input type="checkbox"/> black <input type="checkbox"/> grey | ○ |
| | <input type="checkbox"/> colored hem: fill color by 3.3C | ○ |
| 2.3 <input type="checkbox"/> removable pocket (under seat) | color: fill color by 3.3B | ○ |

3 BACKREST

- 3.1A ☐ textilní základní (základní polstrování bez možnosti nastavení tuhosti - barva viz. odd 3.3A)
- 3.1B ☐ fabric textile with fixing straps (for multipoint adjustment of backrest stiffness, upper cushion textile cover select in 3.3)
- 3.1C ☐ anatomic backrest KURY - fill special order form! >

| 3.2 type of backrest | support tube | hight of textile part | |
|---|---|-------------------------------|---|
| 3.2A <input type="checkbox"/> fixed (see pic. 1/ part A) | <input type="checkbox"/> vertical | 26 28 30 32 34 36 38 40 42 45 | ✓ |
| | <input type="checkbox"/> angular | | |
| 3.2B <input type="checkbox"/> telescopic handles (see pic. 2/ part A) | <input type="checkbox"/> vertical | | ○ |
| | <input type="checkbox"/> angular | | ○ |
| | <input type="checkbox"/> for accompanist/assistant of more than 170 cm high | | ○ |

- 3.2C ☐ shortened handles -6 cm - see pict. 4











3.3 upholstery - fabric, color (If the backrest is not ordered, the color applies only to the seat.)

| | | | | |
|---|--|---|---|---|
| 3.3A <input type="checkbox"/> basic for 3.1A | - backrest: kortexin → - seat: kortexin → | kortexin fabric color: | <input type="checkbox"/> black <input type="checkbox"/> dark grey <input type="checkbox"/> yellow <input type="checkbox"/> red <input type="checkbox"/> violet <input type="checkbox"/> blue <input type="checkbox"/> dark green | ✓ |
| 3.3B <input type="checkbox"/> „DOBBY” for 3.1B | - cushion: Dobby black + kortexin → - seat: kortexin → | | | ✓ |
| 3.3C <input type="checkbox"/> „Air-Flow” for 3.1B | - cushion: TEBOX black + supplementary TEBOX → - seat: kortexin black | <input type="checkbox"/> green (1) <input type="checkbox"/> turquoise (2) <input type="checkbox"/> red (3) <input type="checkbox"/> yellow (4) <input type="checkbox"/> grey (5) <input type="checkbox"/> pink (6) <input type="checkbox"/> black (7) <input type="checkbox"/> blue dark (8) <input type="checkbox"/> beige (9) <input type="checkbox"/> ruby (10) <input type="checkbox"/> orange (11) | | ○ |

- 3.4 ☐ textile back pocket 30 (w) x 20 (h) x 2 (d) cm - ONLY with upholstery 3.3B (DOBBY) or 3.3C (Air-Flow)

| | | | | | |
|------|--------------------------|---|--|--|--|
| 3.5 | <input type="checkbox"/> | additional adjustable headrest (see pic. 3/ part A) | headrest (upholstered, TEBOX fabric) | | |
| 3.5A | | | <input type="checkbox"/> small A = 130, B = 210 | <input type="checkbox"/> large A = 170, B = 250 | <input type="checkbox"/> shallow form C = 50mm <input type="checkbox"/> deep form C = 100mm |
| 3.5B | | | hem color: <input type="checkbox"/> black <input type="checkbox"/> grey | | |
| 3.5C | | | <input type="checkbox"/> adjustable and fixing system of headrest (rack) | | |

4 ARMRESTS - tilting and removable

| | | | | |
|------|---|---|---|-----------------|
| 4.1A | <input type="checkbox"/> E 01 long  <input type="checkbox"/> low (15) <input type="checkbox"/> high (20) | <input type="checkbox"/> E 02 low (18) <input type="checkbox"/> E 20 mid (20) <input type="checkbox"/> E 03 high (23)  | <input type="checkbox"/> E 05 long  <input type="checkbox"/> low (18) <input type="checkbox"/> high (23) | ✓ |
| | <input type="checkbox"/> E 06 long divided  <input type="checkbox"/> low (18) <input type="checkbox"/> high (23) | <input type="checkbox"/> E 07 low  | <input type="checkbox"/> E 08 high  | |
| 4.1B | <input type="checkbox"/> E 02M low with extendable handle  | | | ○ |
| 4.2 | Mudguards - non-tilting and non-removable | | | |
| 4.2A | <input type="checkbox"/> E 04 mudguards ABS plastic  | | | ✓ |
| 4.2B | <input type="checkbox"/> E 04M mudguards with extendable handle  | | | ○ |
| 4.2C | <input type="checkbox"/> E 04 CARBON mudguards  | | | <div>CE</div> ○ |

5 COLOR OF FRAME PARTS - SIDE CONSTRUCTIONS

| | | | | | | |
|------|---|--|--|---|---|------|
| 5.1A | <input type="checkbox"/> basic RAL | <input type="checkbox"/> (1) glossy white | <input type="checkbox"/> (24) white matt | <input type="checkbox"/> (2) yellow | <input type="checkbox"/> (3) red | ✓ |
| | | <input type="checkbox"/> (4) ruby | <input type="checkbox"/> (5) light violet | <input type="checkbox"/> (6) violet | <input type="checkbox"/> (7) dark violet | |
| 5.1B | <input type="checkbox"/> metallic | <input type="checkbox"/> (8) light blue | <input type="checkbox"/> (9) dark blue | <input type="checkbox"/> (10) turquoise | <input type="checkbox"/> (11) light green | ○ |
| | | <input type="checkbox"/> (12) green | <input type="checkbox"/> (13) glossy black | <input type="checkbox"/> (26) black matt | | |
| 5.1B | <input type="checkbox"/> metallic | <input type="checkbox"/> (15) violet | <input type="checkbox"/> (16) blue | <input type="checkbox"/> (17) red | <input type="checkbox"/> (18) anthracit | ○ |
| | | <input type="checkbox"/> (19) turquoise | <input type="checkbox"/> (20) silver | <input type="checkbox"/> (21) light green | <input type="checkbox"/> (22) yellow | |
| 5.2 | <input type="checkbox"/> set of colored parts - elox | rear whl. holders + castor whl. forks + brake levers | | | | CE ○ |
| | | <input type="checkbox"/> black | <input type="checkbox"/> silver | <input type="checkbox"/> blue | <input type="checkbox"/> red | |
| 5.3 | <input type="checkbox"/> fastening parts (screws, nuts) | <input type="checkbox"/> white zinc | <input type="checkbox"/> black zinc | | | ✓ |

6 FOOTREST & FOOTPLATE

| | | | | | | |
|------|--|---|----------------------------------|-----------------------|---|---|
| 6.1 | <input type="checkbox"/> one-piece footplate | <input type="checkbox"/> angle 90° <input type="checkbox"/> angle 100° | (Angle can be changed). | side tube of footrest | <input type="checkbox"/> short (according to parameter Y - from 36 to 42 cm) | ✓ |
| 6.2 | <input type="checkbox"/> divided footrest (for width 40 cm only ! see S paramet.) | <input type="checkbox"/> changeless an. 90° <input type="checkbox"/> changeless an. 100° | | | <input type="checkbox"/> long (according to parameter Y - from 44 to 52 cm) | ✓ |
| 6.3A | <input type="checkbox"/> leg support belt | | <input type="checkbox"/> textile | | | ✓ |
| 6.3B | | | <input type="checkbox"/> leather | | | ○ |

7 CASTOR WHEELS

| | wheel rim | sign / size (mm) | rim color | | | tire | |
|------|------------------------------------|--|---|---|-------------|---|---|
| 7.1 | <input type="checkbox"/> plastic | <input type="checkbox"/> PU 100 | black | | | <input type="checkbox"/> Grey Urethane | ✓ |
| | | <input type="checkbox"/> PU 125 | | | | | |
| | | <input type="checkbox"/> PU 150 | | | | <input type="checkbox"/> Black Urethane | |
| | | <input type="checkbox"/> PU 175 | | | | | |
| 7.2A | <input type="checkbox"/> aluminium | <input type="checkbox"/> ALPU 100 x 32 | 5-spoke rim | silver | | Black Urethane | ○ |
| 7.2B | | <input type="checkbox"/> ALPU 100 x 35 | 3-spoke rim | | | | ○ |
| 7.2C | | <input type="checkbox"/> ALPU 125 x 32 | 5-spoke rim | | | | ○ |
| 7.2D | | <input type="checkbox"/> ALPU 125 x 35 | 3-spoke rim | | | | ○ |
| 7.2E | | <input type="checkbox"/> ALPU 100 x 32 C | <input type="checkbox"/> black <input type="checkbox"/> red <input type="checkbox"/> gold | <input type="checkbox"/> blue <input type="checkbox"/> orange <input type="checkbox"/> yellow-green | 5-spoke rim | | ○ |
| 7.2F | | <input type="checkbox"/> ALPU 125 x 32 C | | | | | ○ |

| 8 REAR WHEELS | | | | | | |
|--|---|---|-------------------|--|--|---|
| 8.1 | rear whl. axis angle | <input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 4° (up to 80 kg weight of user and with mudguard E04 usage only) | | | | ✓ |
| 8.2 PREMIUM wheels (complete rear wheels - rim+tire+handrim) | | | | | | |
| | rear wheels marking | specification | wheel tires | | spokes | |
| | | rim, hub, handrim | size | tire (cannot be changed) | | |
| 8.2A | <input type="checkbox"/> Spinergy SPOX Everyday | alu rim, spinergy design handrim: black elox* | 24 x 1" | Schwalbe MPE | vertical, 18x spinergy carbon fiber - black spokes | ○ |
| 8.2B | <input type="checkbox"/> Spider | alu rim + magnesium alloy handrim: silver elox | 24 x 1" | Schwalbe MPE | 6x magnesium spokes, attractive design | ○ |
| 8.2C | <input type="checkbox"/> Dynamic COLOR (KURY design wheels) <div>CE</div> | double rim, black elox design hub | 24 x 1" | Schwalbe MPE | <input type="checkbox"/> white paint | ○ |
| | | handrim: hard elox handrim: <input type="checkbox"/> 15 <input type="checkbox"/> 20mm | | | <input type="checkbox"/> black elox <input type="checkbox"/> red elox <input type="checkbox"/> blue elox <input type="checkbox"/> gold elox <input type="checkbox"/> orange elox <input type="checkbox"/> yellow-green elox | |
| 8.2D | <input type="checkbox"/> Dynamic DYS | double rim, hard elox handrim: <input type="checkbox"/> 15 <input type="checkbox"/> 20mm | 24 x 1" | Schwalbe MPE | desing spokes - 18 stainless steel intersec wires | ○ |
| 8.2E | <input type="checkbox"/> SIX STAR black* | double rim, black elox handrim: black elox | 24 x 1" | Schwalbe MPE | 6x3 vertical stainless spokes (2,6 mm), For rear whls. angle 0° only! | ○ |
| 8.2F | <input type="checkbox"/> Ultra light ONE* | double rim, black elox handrim: black elox | 24 x 0,9" | Schwalbe ONE | 24 spokes, black | ○ |
| 8.2G | <input type="checkbox"/> Black* | double rim, black elox handrim: black elox | 24 x 1" | Schwalbe MPE | 24 stainless spokes black (2 mm) | ○ |
| 8.2H | <input type="checkbox"/> SIX STAR silver* | double rim, silver elox handrim: silver elox | 24 x 1" | Schwalbe RightRun Active | 6x3 vertical stainless spokes (2,6 mm), For rear whls. angle 0° only! | ○ |
| 8.2I | <input type="checkbox"/> Cross* | GrandRock, black elox handrim: <input type="checkbox"/> 15 <input type="checkbox"/> 20mm | 24 x 1,90" | Schwalbe Landcruiser | 36 stainless spokes (2 mm) | ○ |
| 8.2J | <input type="checkbox"/> Country | simple rim., black elox handrim: <input type="checkbox"/> 15 <input type="checkbox"/> 20mm | 24 x 1,95" | PNEU Offroad | 36 stainless spokes | ○ |
| 8.3 STANDARD rim - tires and handrims must be ordered separately (see 8.4 and 8.5) | | | | | | |
| 8.3A | <input type="checkbox"/> ultra light | double rim, white elox handrim see 8.5 | select see 8.4 | • PNEU usage ONLY (!) • NOT for drum brakes! | 24 spokes, stainless | ○ |
| 8.3B | <input type="checkbox"/> standard silver | simple rim handrim see 8.5 | select see 8.4 | • PNEU + PU: select - see 8.4 | 36 spokes, zinc | ✓ |
| 8.3C | <input type="checkbox"/> standard black | simple rim handrim see 8.5 | select see 8.4 | • PNEU + PU: see 8.4 • NOT for drum brakes! | 24 spokes, stainless | ✓ |
| 8.4 PNEU TIRES (marking) | | | size | specification | | |
| 8.4A | <input type="checkbox"/> Schwalbe Marathon Plus Evolution Reflex | | 24 x 1" | black + reflex strip | SmartGuard - protection belt, press. 6-10 Bar; 500g | ○ |
| 8.4B | <input type="checkbox"/> Schwalbe RightRun Active | | 24 x 1" | <input type="checkbox"/> black-grey | puncture protection layer; hand-friendly 2Grip side walls, and attractive two color compound press. 6-10 Bar; 315g | ○ |
| | | | | <input type="checkbox"/> black-blue | | |
| | | | | <input type="checkbox"/> black-red | | |
| | | | | <input type="checkbox"/> black | | |
| 8.4C | <input type="checkbox"/> Schwalbe ONE | | 24 x 1" | black | very light and extremely cut resistant; folding tire; press. 6-9 Bar; 205 g | ○ |
| 8.4D | <input type="checkbox"/> Schwalbe DownTown | | 24 x 1" | grey-black | PNEU with a puncture protection layer and hand-friendly 2Grip side walls, press. 6,0-10,0 Bar; 420g | ○ |
| 8.4E | <input type="checkbox"/> Standard pneu | | 24 x 1" | grey | standard PNEU tire for ultralight and standard rim (see 8.3). | ✓ |
| 8.4G | <input type="checkbox"/> Standard pneu (Rubena V69 Hook) | | 24 x 1 3/8" | <input type="checkbox"/> grey <input type="checkbox"/> black | | |
| URETHAN TIRES (For standard rims ONLY - see 8.3B+C!) | | | | | | |
| 8.4I | <input type="checkbox"/> PU 24 x 1" | | 24 x 1" | <input type="checkbox"/> grey <input type="checkbox"/> black | polyurethane tires for standard rim (see 8.3). | ○ |
| | <input type="checkbox"/> PU 24 x 1 3/8" | | 24 x 1 3/8" | | | |


* Any change of tire or handrim is possible based on the current offer and special price calculations.

✓ basic equipment included in price ○ additional item with surcharge (see pricelist) **CE** - included in COMFORT EDITION packet

| 8.5 | handrim | specification - distance from wheel (It can be changed.) | | |
|------|--|---|--|---|
| 8.5A | <input type="checkbox"/> smooth aluminium | <input type="checkbox"/> 15 <input type="checkbox"/> 20 | | ✓ |
| 8.5B | <input type="checkbox"/> hard anodized | <input type="checkbox"/> 15 <input type="checkbox"/> 20 | | ○ |
| 8.5C | <input type="checkbox"/> rubber coating | <input type="checkbox"/> 15 <input type="checkbox"/> 20 | | ○ |
| 8.5D | <input type="checkbox"/> smooth stainless | <input type="checkbox"/> 15 <input type="checkbox"/> 20 | | ○ |
| 8.5E | <input type="checkbox"/> 8 tabs | <input type="checkbox"/> 15 <input type="checkbox"/> 20 with standard rim only | | ○ |
| 8.5F | <input type="checkbox"/> 12 tabs | <input type="checkbox"/> 15 <input type="checkbox"/> 20 with standard rim only | | ○ |
| 8.5G | <input type="checkbox"/> covered with leather | <input type="checkbox"/> 15 <input type="checkbox"/> 20 leather type: <input type="text"/> Thread color: <input type="text"/> | | ○ |
| 8.6 | brakes | <input type="checkbox"/> standard (straight) <input type="checkbox"/> longer lever <input type="checkbox"/> top lever bent | | ✓ |
| 8.6B | | <input type="checkbox"/> compact brakes | | ○ |
| 8.6C | | <input type="checkbox"/> wheelchair brakes for assistant person | | ○ |
| 8.7A | <input type="checkbox"/> one hand rear wheel control | <input type="checkbox"/> right wheel <input type="checkbox"/> left wheel | | ○ |

| 9 | ACCESSORIES AND SPECIAL EQUIPMENT | | | |
|------|---|--|--|---|
| 9.1A | <input type="checkbox"/> leg support - double belt (front and back) - adjustable, w. 45 mm | <input type="checkbox"/> leather <input type="checkbox"/> textile | | ○ |
| 9.1B | <input type="checkbox"/> backward roll-over prevention support (rear socket tube)* | | | ○ |
| 9.1C | <input type="checkbox"/> quick release axles for quadriplegic | <input type="checkbox"/> 1 pcs <input type="checkbox"/> 2 pcs (price for 1 pcs) | | ○ |
| 9.1D | <input type="checkbox"/> carbon plug of front fork (pair) | (price for 2 pcs) | | ○ |
| 9.1E | <input type="checkbox"/> standard pump with clips | | | ○ |
| 9.1F | <input type="checkbox"/> pump with clips 8 ATM (greater pressure fill - Schwalbe MPE, Schwalbe RR etc.) | | | ○ |
| 9.1G | <input type="checkbox"/> silicone handrim covers | <input type="checkbox"/> red <input type="checkbox"/> violet <input type="checkbox"/> pink <input type="checkbox"/> yellow <input type="checkbox"/> green <input type="checkbox"/> blue <input type="checkbox"/> black (price for 2 pcs) | | ○ |
| 9.1H | <input type="checkbox"/> buckle for locking of folding wheelchair (If selected removable pocket - viz bod 2.3 - this buckle is included.) | | | ○ |
| 9.1I | <input type="checkbox"/> original backrest embroidery - The price will be calculated according to the range of work. | | | |
| 9.1J | <input type="checkbox"/> protective sleeves for rear wheels - waterproof fabric, interior pollution prevention | (price for 2 pcs) | | ○ |
| 9.2 | rear wheel guards | | | |
| 9.2A | <input type="checkbox"/> basic | <input type="checkbox"/> white <input type="checkbox"/> transparent <input type="checkbox"/> black | | ○ |
| 9.2B | <input type="checkbox"/> designed (see swatch) | <input type="checkbox"/> rings <input type="checkbox"/> night <input type="checkbox"/> violet spiral <input type="checkbox"/> flower I <input type="checkbox"/> wheel <input type="checkbox"/> eyes <input type="checkbox"/> blue spiral <input type="checkbox"/> flower II <input type="checkbox"/> red-blue cascade <input type="checkbox"/> turquoise cascade | | ○ |
| 9.2C | <input type="checkbox"/> own design | user's own picture or photo will be print onto rear wheel guard (jpeg, pdf, eps, gif, psd, ai, min. resolution 150 dpi, size 1600x1600 pix - min. size of file (!) 0,7 MB) | | ○ |
| 9.3 | fixing and positioning accessories (see current offer) | approximate dimensions (cm) | | |
| 9.3A | <input type="checkbox"/> AD10 abdominal fixing belt | universal size: (h) 8, (w) 48 | | ○ |
| 9.3B | <input type="checkbox"/> AD11 abdominal fixing belt for chair | <input type="checkbox"/> size 0: (h)12, (w) 40 <input type="checkbox"/> size 1: (h)16, (w) 54 | | ○ |
| 9.3C | <input type="checkbox"/> AD11E abdominal fixing belt thin | universal size: (h)12, (w) 54 | | ○ |
| 9.3D | <input type="checkbox"/> AD12 perineal fixing belt for chair | <input type="checkbox"/> size 0: (h)12, (w) 40 <input type="checkbox"/> size 1: (h)16, (w) 54 | | ○ |
| 9.3E | <input type="checkbox"/> AD13 fixing belt for chair (torso) | <input type="checkbox"/> size 0: (h)12, (w) 40 <input type="checkbox"/> size 1: (h)16, (w) 54 | | ○ |
| 9.3F | <input type="checkbox"/> AD14 fixing belt for chair (torso and pelvis) | <input type="checkbox"/> size 0: (h)12, (w) 40 <input type="checkbox"/> size 1: (h)16, (w) 54 | | ○ |
| 9.3G | <input type="checkbox"/> ADX20 narrow abdominal belt | <input type="checkbox"/> size 0: from w. 37 <input type="checkbox"/> size 1: from w. 48 | | ○ |
| 9.3H | <input type="checkbox"/> ADX22 4 fixpoint belt - vest „X“, zipper 17 cm | <input type="checkbox"/> size 0: (w) 30+(w) 30 <input type="checkbox"/> size 1: (w) 43+(w) 33 | | ○ |
| 9.3I | <input type="checkbox"/> ADX23 shoulder harness | <input type="checkbox"/> size 0: (h)27 <input type="checkbox"/> size 1: (h) 33 | | ○ |
| 9.3J | <input type="checkbox"/> ADX30 abduction support | | | ○ |
| 9.3K | <input type="checkbox"/> ADX2 leg fixation in the ankle area to the footrest | | | ○ |
| 9.4 | wheelchair trays / tables | Tray for wheelchair, 6 mm PC (polycarbonate), with 3 bended sides 30 mm | | |
| 9.4A | <input type="checkbox"/> basic clear (transparent) | | | ○ |
| 9.4B | <input type="checkbox"/> own printed design | | | ○ |
| 9.4C | <input type="checkbox"/> padded elbow parts | | | ○ |
| 9.4D | <input type="checkbox"/> partly padded tray | | | ○ |
| 9.4E | <input type="checkbox"/> full padded tray | | | ○ |

* Back frame part is changed for installation of backward roll-over prevention support (see 9.1B)



✓ basic equipment included in price ○ additional item with surcharge (see pricelist)