

Lisa Rehab folding buggy	Order no.: -----
------------------------------------	-------------------------

Lisa

Rehab folding buggy

Size 1 for children from approx. 3 to 8 years
 Size 2 for children from approx. 6 to 16 years



Standard equipment

Compact folding buggy in a very stable, light-weight aluminium design, back support adjustment for one-hand operation

- Seat upholstery and reinforced back support upholstery
- Adjustable and divisible leg support
- Front and rear wheels with ball bearings
- Wheel lock
- Screw connections guided by sleeves
- Standard lap belt
- Upholstery colour blue
- Front wheels with integrated swivel lock
- Pre-mounted tip-assist, right
- Frame colour: silver

Standard tyres:

- EVA wheels, front, 180 mm
- EVA wheels, rear, 250 mm

<input type="checkbox"/> Provide estimate only <input type="checkbox"/> Order	
Billing address/customer number	
Delivery address	
Sales person	Name Date Signature
By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.	

Standard model			
HR32130001-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 1 with seat upholstery seat depth presetting 30 cm/12"	Add. charge
HR32130002-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 1 with seat upholstery seat depth presetting 35 cm/14"	Add. charge
HR32140001-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 2 with seat upholstery seat depth presetting 35 cm/14"	Add. charge
HR32140002-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 2 with seat upholstery seat depth presetting 40 cm/16"	Add. charge

Accessories			
HR32137100-027	<input type="checkbox"/>	Arm support, piece (please order two pieces)	Add. charge
HR16001102	<input type="checkbox"/>	„Techno“ buggy board	Add. charge
HR32135300-343	<input type="checkbox"/>	Simplified padding, size 1 (coro-blue)	Add. charge
HR32145300-343	<input type="checkbox"/>	Simplified padding, size 2 (coro-blue)	Add. charge
HR32135200-343	<input type="checkbox"/>	Seat pad with multi-functional back, size 1 (coro-blue/dark blue)	Add. charge
HR32145200-343	<input type="checkbox"/>	Seat pad with multi-functional back, size 2 (coro-blue/dark blue)	Add. charge

HR32137700-031	<input type="checkbox"/>	Head support, size 1, matching multifunctional back support, size 1	Add. charge
HR32147700-031	<input type="checkbox"/>	Head support, size 2, matching multifunctional back support, size 2	Add. charge
HR32070060-031	<input type="checkbox"/>	Head support size 1, padded (not available in combin. with seat pad with multifunctional back support)	Add. charge
HR32086500-031	<input type="checkbox"/>	Head support size 2, padded (not available in combin. with seat pad with multifunctional back support)	Add. charge
HR32136100	<input type="checkbox"/>	Five-point belt with neck pad	Add. charge
HR32135500-031	<input type="checkbox"/>	Lateral guides with variable seat width reduction (4-6-8 cm) size 1 (only available in combination with pad for multifunctional back support)	Add. charge
HR32145500-031	<input type="checkbox"/>	Lateral guides with variable seat width reduction (4-6-8 cm) size 2 (only available in combination with pad for multifunctional back support)	Add. charge
HR32137200-031	<input type="checkbox"/>	Flip-down pommel size 1	Add. charge
HR32147200-031	<input type="checkbox"/>	Flip-down pommel size 2	Add. charge
HR32075800-031	<input type="checkbox"/>	Groin strap for multifunctional back support size 1	Add. charge
HR32085800-031	<input type="checkbox"/>	Groin strap for multifunctional back support size 2	Add. charge
HR32137800	<input type="checkbox"/>	Grab rail size 1	Add. charge
HR32147800	<input type="checkbox"/>	Grab rail size 2	Add. charge
HR32135700-031	<input type="checkbox"/>	Padding for grab rail size 1	Add. charge
HR32145700-031	<input type="checkbox"/>	Padding for grab rail size 2	Add. charge
HR32010080	<input type="checkbox"/>	Hook-and-loop foot straps with instep pad size 1	Add. charge
HR32020080	<input type="checkbox"/>	Hook-and-loop foot straps with instep pad size 2	Add. charge
HR32200120-343	<input type="checkbox"/>	Winter slip sack coro blue with detachable front part, with slits for straps and accessories, size 2 (seat depth up to approx. 34 cm)	Add. charge
HR32400120-343	<input type="checkbox"/>	Winter slip sack coro blue with detachable front part, with slits for straps and accessories, size 3 (seat depth up to approx. 39 cm)	Add. charge
HR32100100-074	<input type="checkbox"/>	Rain cape, dark blue size 1	Add. charge
HR32010100-074	<input type="checkbox"/>	Rain cape, dark blue size 2	Add. charge
HR32135800	<input type="checkbox"/>	Folding storage bag size 1	Add. charge
HR32145800	<input type="checkbox"/>	Folding storage bag size 2	Add. charge
HR32136200-094	<input type="checkbox"/>	Canopy in upholstery colour with transparent rain cover	Add. charge

Information

Technical data		
	Size 1	Size 2
Seat width (with/without padding)	32/36 cm	37/41 cm
Seat depth	30 or 35 cm	35 or 40 cm
Back support height	65 cm	70 cm
Back support adjustment	30° (in steps 4 x 7.5°)	30° (in steps 4 x 7.5°)
Lower leg length	27 – 40 cm	27 – 40 cm
User weight (max.)	50 kg	60 kg
Weight	13.6 kg	14.6 kg

Accessories		
HR51040000-027	<input type="checkbox"/>	Robby, size 1 (seat depth 20 – 24 cm, back support height incl. head section 63 – 72.5 cm) Bathing and lying support, includes washable net upholstery, lap belt, leg straps; max. load capacity 30 kg Add. charge
HR51050000-027	<input type="checkbox"/>	Robby, size 2 (seat depth 31 – 42.5 cm, back support height incl. head section 89 – 94 cm) Bathing and lying support, includes washable net upholstery, lap belt, leg straps; max. load capacity 60 kg Add. charge



Robby bathing support