



"WITHIN 1 HOUR OF DELIVERY I WAS IN THE BUSH AND ON THE BEACH! A DREAM COME TRUE AND PROBABLY ONE OF THE BEST FEELINGS I'VE EVER HAD, ABSOLUTE FREEDOM."

Nick
New South Wales

MAGIC MOBILITY EXTREME X8

HIT THE ROAD

Or the sand, or the bush, or the snow.*

*Not to be driven through water

MAKING YOUR CHAIR WORK FOR YOU

Magic Mobility wheelchairs can be found on building sites, farms, stables, national parks - all allowing people to continue to work.



DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help you find them.



**AUSTRALIAN
MADE**
—
BUILT TO LAST
—
**CUSTOMISED
TO YOUR
REQUIREMENTS**

Magic
MOBILITY™

Extreme
X8™



MAGIC MOBILITY EXTREME X8



The Extreme X8 is the flagship wheelchair in the Magic Mobility range. With each release we improve the design with a focus on enhanced performance and reliability. The Extreme X8 is the most compact 4WD wheelchair on the market.

Basic specifications:

- Length: 1150mm
- Width: 710mm
- Seat to floor height: 470mm (490mm with swing away leg rests)
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift or tilt
- Top speed: 10kph (6.2mph)
- Motors: 700W Gear in Line
- Wheels: 14" (145/70-6) black knobby
- Batteries: 73Ahr Gel (40Ahr option)
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-23" wide).*

Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lie down, power elevating legrests
- Seating and backrests: MPS, rehab backrest and adjustable seat pans, foam in place, custom seating (built by our in house team of fabricators), a choice of off the shelf backrests and seat cushions
- Headrests and mounts: MPS, rehab (single and multi axis stem, pivoting, Y-type, L-Type), any custom requirement
- Legrests: centre post, side mount, swingaway, power elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxy and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- WC19 crash test compliant tie downs.



EXTREME X8

"WE'VE BEEN
CREATING
CUSTOM MADE
WHEELCHAIRS
FOR 20 YEARS,
SO THAT YOU CAN
LIVE YOUR LIFE"

Ashley Daff
Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN OWNED, DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

*Various model options excluded

Specifications subject to change. Check online for current version



Visit us online to see how
we change lives every day:
magicmobility.com.au

Postal Address

Magic Mobility Pty Ltd
2/112 Browns Road
Noble Park VIC 3174 Australia

Tel +61 3 8791 5600

Fax +61 3 9701 0144

Email sales@magicmobility.com.au



Join the
Magic Mobility
Facebook and
LinkedIn
communities!

