

LEVO Dealer:

Date:

Dealer reference / Client name:

Anglikerstrasse 20  
CH-5610 Wohlen  
Phone: +41 56 618 44 11  
Fax: +41 56 618 44 10  
e-mail: office@levo.ch


Order number:

Standard configuration 60.001.7010

- Transport dimensions: width 60 cm, length 98 cm, height 45 cm (incl. footplate)
- M/E Wheelchair: Manual driving, seating and powered standing
- Seat plate and backrest upholstery
- Flip-up armrests
- Armrest height and armrest angle adjustable
- Armrests round with built-in control switch
- Fold down backrest
- Standard colors: platinum, red, yellow, blue or black
- Two-piece flip-up footrest
- Pneumatic rear wheels and pump
- Batteries and battery charger
- Chest strap and knee support with wedge



### Anthropometric Data




#### 1 Hip width

Please check the correct box according to the measured hip width in a sitting position:

35 cm  48 cm

#### Seat width:

35 cm	60.004.7010
38 cm	Option
40 cm	Option
42 cm	Option
45 cm	Option
48 cm	Option
50 cm	60.004.7020



#### 2 Upper leg length

There are two seat frame sizes available for the LCEV: S - M and L - XL. The upper leg length determines which frame size and which armrest size to choose for the wheelchair.

Please fill in the measured upper leg length:

Upper leg length: 35 cm  62,2 cm

Example:

Upper leg length 35 cm  
- 6 cm  
**Seat upholstery depth 29.0 cm**

#### Seat upholstery depth:

S-M Frame	29.0 cm	Option
	31.2 cm	Option
	33.4 cm	Option
	35.6 cm	Option
	37.8 cm	Option
40.0 cm	Option	

L-XL Frame	43.0 cm	Option
	45.2 cm	Option
	47.4 cm	Option
	49.6 cm	Option
	51.8 cm	Option
	54.0 cm	Option

#### 3 Lower leg length

The lower leg length determines the footrest position and the seat height of the wheelchair. The seat height determines the wheel size and with that the total height of the wheelchair.

Please fill in the measured lower leg length:

Lower leg length: 36 cm  56 cm

Please check the desired seat height and footrest position, under consideration of the seat cushion thickness.

Example:

Lower leg length 51 cm  
Seat cushion height - 3 cm  
**Footrest position 48 cm**

#### Footrest position:

Seat height 48	36 cm	Option
	38 cm	Option
	40 cm	Option
	42 cm	Option
	44 cm	Option
	46 cm	Option
	48 cm	Option
	50 cm	Option
	52 cm	Option
	54 cm	Option
56 cm	Option	

Now you may choose between seat height 51, 54 or 57 cm

Please measure the anthropometric data as shown. If possible, the user should be in a sitting position and wearing shoes.

**4 Elbow height**

Standard armrest height is 18 cm


14 cm  22 cm

**Wheels and position**

Floor to seat height in the front, measured without seat cushion

Rear wheel 22", castor 5"	48 cm	<input type="checkbox"/>	Option
Rear wheel 24", castor 5"	51 cm	<input type="checkbox"/>	Standard
Rear wheel 24", castor 5"	54 cm	<input type="checkbox"/>	Option
Rear wheel 26", castor 6"	57 cm	<input type="checkbox"/>	60.004.0010

**Armrest**



① Armrest round with hand control		60.005.0040
① Armrest 6 cm/2.4" wide with built-in control switch*		60.005.7010
① Armrest 6 cm/2.4" wide with hand control*		60.005.7020
② Armrest 5 cm/2" longer <i>*Not in combination with swing-away chest support</i>		60.005.7030
③ Dual switch for table tray control, <i>table tray not included</i>		60.007.0010
④ Swing-away chest support <i>Rotatable to the side and flip-up, for increased safety feeling. Distance "D" between backrest tube and the beginning of radius, see photo</i>	<b>D</b> 27 cm 32 cm 37 cm	60.005.7040 60.005.7041 60.005.7042

**Stand-up operation**




Right	<input type="checkbox"/>	Standard
Left	<input type="checkbox"/>	Option

**Frame color**

Platinum		Standard
Red		Option
Yellow		Option
Blue		Option
Black		Option
Blue metallic		60.008.0020
Green metallic		60.008.0020
Red metallic		60.008.0020
Silver metallic		60.008.0020
Orange reflex		60.008.0030
Lemon reflex		60.008.0030
Silver reflex		60.008.0030
Aluminium polished		60.008.0050
Trafic white	RAL: 9016	60.008.0010
Other RAL color	RAL: <input type="text"/>	60.008.0010

**Backrest**



Standard sling backrest height	31 cm 33 cm 35 cm 37 cm 39 cm 41 cm	<input type="checkbox"/>	Standard Option Option Option Option Option
① V-Trak standard backrest, height 40 cm			60.009.0020
② V-Trak medium backrest, height 60 cm			60.009.0030
V-Trak high backrest, height 80 cm			60.009.0040
V-Trak Axxis reduced height backrest, height 31 cm			60.009.0070
V-Trak Axxis low backrest, height 40 cm			60.009.0050
V-Trak Axxis high backrest, height 61 cm			60.009.0060
Without backrest, <i>be aware net price is lower after discount</i>			60.109.6810

**Backrest options**

Backrest extension 15 cm, *for standard backrest*


Seat width ▶	35	38	40	42	45	48	50
Height ▼							
15 cm	32.064.1735	32.064.1738	32.064.1740	32.064.1742	32.064.1745	32.065.1748	32.065.1750

V-Trak backrest extension, 20 cm/8"	<input type="checkbox"/>	32.069.99
-------------------------------------	--------------------------	-----------

Seat cushion with memory foam depending on the seat depth and seat width




Seat width ▶ Seat depth ▼	35	38	40	42	45	48	50
S	30.100.3740	30.100.3770	30.100.3780	30.100.3790	30.100.3800	30.100.3802	30.100.3803
M	30.100.3822	30.100.3820	30.100.3830	30.100.3840	30.100.3850	30.100.3855	30.100.3856
L	30.100.3857	30.100.3860	30.100.3870	30.100.3880	30.100.3890	30.100.3900	30.100.3902
XL	30.100.3907	30.100.3910	30.100.3920	30.100.3930	30.100.3940	30.100.3950	30.100.3956




Additional seat cushion cover

Seat width ▶ Seat depth ▼	35	38	40	42	45	48	50
S	30.100.3741	30.100.3771	30.100.3781	30.100.3791	30.100.3801	30.100.3804	30.100.3805
M	30.100.3823	30.100.3821	30.100.3831	30.100.3841	30.100.3851	30.100.3854	30.100.3953
L	30.100.3858	30.100.3861	30.100.3871	30.100.3881	30.100.3891	30.100.3901	30.100.3903
XL	30.100.3908	30.100.3911	30.100.3921	30.100.3931	30.100.3941	30.100.3951	30.100.3957



Urin proof seat cushion cover

Seat width ▶ Seat depth ▼	35	38	40	42	45	48	50
S	30.100.3743	on request	on request	on request	on request	on request	on request
M	30.100.3824	30.100.3864	on request	on request	on request	on request	on request
L	on request	30.100.3863	30.100.3873	on request	30.100.3892	on request	on request
XL	on request	on request	on request	on request	on request	on request	on request



Headrest



V-Trak headrest standard, *only in combination with V-Trak Axxis backrest*  
 V-Trak headrest winged, *only in combination with V-Trak Axxis backrest*  
 Track extension for V-Trak head support, for V-Trak 60 and 80 cm

	02.067.116
	02.067.117
	02.009.0028

Lateral supports and hardware



V-Trak lateral support, pad and cover

Medium	Right		32.068.0030
Medium	Left		32.068.0040
Large	Right		32.068.0010
Large	Left		32.068.0020

Knee support



Standard knee support

	Standard
--	----------

① Rectangular knee wedge

50x100	Option
60x100	Option
80x100	Option
50x150	Option
80x150	Option

② Conical knee wedge, *gives more support, in particular with spasticity*

50-100x100	60.011.0020
90-140x100	60.011.0020
130-180x100	60.011.0020


③ Sports knee support

	60.011.0050
--	-------------

④ Rigid knee support

Size S (10 - 12 cm)	60.011.0060
Size M (14 - 16 cm)	60.011.0060

Footrest



Two-piece flip-up footrest

	Standard
--	----------

① One piece footrest, fixed

	Option
--	--------

② One piece footrest, flip-up backwards, *not in combination with footrest height adapter*

	60.003.0050
--	-------------

③ Footrest height adapter

For seat width 38, 45, 48	60.003.0070
For seat width 40 and 42 cm*	60.003.0080

*For users with short lower leg length. Footrest will be 3 cm narrower than the seat width.  
 \*Only possible with one piece fixed footrest.*

**Footrest**



- ④ Easy transfer foot plate, *foot plate is closer to the floor which makes transfers easier*

20 cm	<input type="checkbox"/>	60.003.0120
26 cm	<input type="checkbox"/>	60.003.0130

Legrest removable for transport

<input type="checkbox"/>	60.003.0060
--------------------------	-------------

Non-slip footrest surface


<input type="checkbox"/>	60.003.0040
--------------------------	-------------

**Foot plate angle**

*Foot plate is adjustable in all three angles.*

80°	<input type="checkbox"/>	Option
90°	<input type="checkbox"/>	Standard
100°	<input type="checkbox"/>	Option

**Rear wheel position**



*Select the appropriate setting for the rear wheels. The shorter the wheelbase is, the more active the wheelchair becomes and the risk of tipping over gets larger. Anti-tip aids are recommended in that case. Normally the wheels are positioned according to the frame size.*

S	<input type="checkbox"/>	Option
M	<input type="checkbox"/>	Option
L	<input type="checkbox"/>	Option
XL	<input type="checkbox"/>	Option

**Wheel options**




Rear wheels with airless tires  
Rear wheels with black marathon plus tires, *only in combination with air tyres and aluminium hand ring*

<input type="checkbox"/>	60.002.0010
<input type="checkbox"/>	60.002.0090




- ① Spinergy wheels  
*Not available in combination with hub locking brake, drum brakes or one arm drive*
- ② Camber, *not in combination with hub locking brake, drum brakes or one arm drive*

Red	<input type="checkbox"/>	60.002.0015
Black	<input type="checkbox"/>	60.002.0015
White	<input type="checkbox"/>	60.002.0015



- ③ Steel hand rim
- ④ Maxgrepp original hand ring, with rubber cover
- ⑤ Maxgrepp Ergo Para CP hand ring

3°	<input type="checkbox"/>	60.002.0080
6°	<input type="checkbox"/>	60.002.0080
Pair	<input type="checkbox"/>	60.002.0070
	<input type="checkbox"/>	60.002.0060
	<input type="checkbox"/>	60.002.0050



- ⑥ Rear wheels with drum brakes
- ⑦ Rear wheels with hub locking brake
- ⑧ One arm drive, *only in combination with standard wheels*

<input type="checkbox"/>	60.002.0025
<input type="checkbox"/>	60.002.0020
<input type="checkbox"/>	60.002.0030

**Accessories**



- ① Table tray

35	38 - 40	42 - 50
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
32.070.0115	32.070.0116	32.070.0117




- ② Swing away table tray, *only possible with combi armrest and hand control*

35 - 40	42 - 50
<input type="checkbox"/>	<input type="checkbox"/>
32.068.1061	32.068.1062



Combi armrest / handcontrol

Right	<input type="checkbox"/>
Left	<input type="checkbox"/>
	32.070.0110



- ③ Push handles standard
- ④ Push handles 7 cm extended
- ⑤ Push handles quick release height adjustable

Pair	<input type="checkbox"/>	32.069.0080
Pair	<input type="checkbox"/>	32.069.0090
Pair	<input type="checkbox"/>	32.070.0167




- ⑥ Anti-tip rolls

S - M	L	XL 54	XL 57
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
32.069.68	32.069.69	32.069.70	32.069.71

Right	<input type="checkbox"/>
Left	<input type="checkbox"/>

*Please select the size which corresponds with the seat height and frame size chosen on page 1*



- ⑦ Anti-tip / tip aid combination, swing away

Right	<input type="checkbox"/>	32.070.0162
Left	<input type="checkbox"/>	32.070.0162

XL-Adapter, *needed for rear wheel position XL*  
XL-Adapter, *needed for rear wheel position XL*

Right	<input type="checkbox"/>	32.070.0165
Left	<input type="checkbox"/>	32.070.0166



- ⑧ Tip aid

Right	<input type="checkbox"/>	32.069.72
Left	<input type="checkbox"/>	32.069.72

Accessoires



⑨ Spoke guard transparent

22"	24"	26"
02.066.1014	02.066.1015	02.066.1017

Other spoke guards on request

Pair




⑩ Foot guide

small	Right	32.069.95
small	Left	32.069.96
large	Right	32.070.0040
large	Left	32.070.0041



⑪ Brake lever extension (6 cm)


Right	30.100.5030
Left	30.100.5030



⑫ Frame reinforcement kit, *Needed when e-fix or similar systems are mounted. Only in combination with anti-tip swing away, not with anti-tip rolls or tip aid.*

32.065.68
-----------

Vehicle tie-down kit, *fixed attachment points for chair transport. Frame reinforcement kit is included*



35	38	40	42	45	48	50
32.070.0501	32.070.0502	32.070.0503	32.070.0503	32.070.0504	32.070.0505	32.070.0506



⑬ Knee support height adapter, *for users with long lower leg lengths, knee holder gets closer to the joint.*

Pair	32.069.85
------	-----------



⑭ Control switch extension

Vertical for round armrest	32.066.0450
Vertical for wide armrest	32.062.1842
Horizontal	32.070.0156

⑮ Seat depth and footrest quick adjustment kit, *10 clamps with bolt, 4 screws*

60.004.0020
-------------

Remarks