



Joy On the Beach



## TRAVELLING MADE EASY

JoB is a unique beach wheelchair, which makes easy getting around on the beach and get right into the water. Its patented wheels uses the concepts of variable and controlled deformation to minimize the efforts of moving on "difficult" or weak grounds, such as sand, pebbles, gravel, snow or dirt.

JoB is made of high quality, corrosion resistant materials specifically tested for marine environment.

Its characteristics make it a unique traveling companion for tourists with special needs.

## DETAILS

JoB is constructed to maximize ease of use.

The wheelchair is completely disassembled for easy transportation inside the carry bag, and no tools are required to dismount or mount it, since all the joints are fixed by quick locking systems.

A wide range of accessories lets you customize the JoB to suit your needs.




Adjustable Armrest




Adjustable Legrest



Quick lock system



Easy disassembly for transport



Third wheel (optional) for self-propelling



Wide range of accessories