

Address

Fax +49 7261-7351-10

Shipping address

Berollka-aktiv Rollstuhltechnik GmbH
 Jahnstraße 16
 74889 Sinsheim - Germany
Phone: +49 7261-7351-0
 info@berollka.de
 www.berollka.de



Com.:

☐ Order
☐ Estimate

Basic Price 2120,- €

Price valid until 31st of December 2018

F1

Mobility base for seating systems

Capacity up to 120kg

Standard equipment:

user weight + Seating system: _____ kg

- Rigid aluminium frame with abduction
- Castor
- Back with hand bow, angle adjustable, height adjustable and tiltable
- Rear wheels with drum brake (and camber)

max. width
seat shell

basic measure

frame	frame tube	seat frame width B	frame width front A	frame width back C	width between wheels D	total width E
Frame size 0 ☐ Code U700	30mm	26 cm	43 cm	34 cm	37,5 cm 40,5cm*	55,5 cm 52cm*
Frame size 1 ☐ Code U701	30mm	32 cm	49 cm	40 cm	43,5 cm 46,5cm*	61,5 cm 58cm*
Frame size 2 ☐ Code U702	30mm	38 cm	55 cm	46 cm	49,5 cm 52,5cm*	67,5 cm 64cm*
Frame size 3 ☐ Code U703	30mm	44 cm	61 cm	52 cm	55,5 cm 58,5cm*	73,5 cm 70cm*

* Note: for size 2 and 3 we suggest the frame with two gas springs for seat adjustment (see CODE U708) !!!

* without wheel cover and camber

* width with 12" wheel

frame colour

technical data

- ☐ 06 dormantred
☐ 16 RAL 5002 light blue
☐ 17 dormant blue
☐ 19 silver metallic
☐ 20 grey metallic
☐ 22 RAL 9011 graphit black
☐ 23 RAL 9005 black
☐ 24 RAL 1028 melon yellow
☐ 25 RAL 9010 white
☐ 27 mandarin met.
☐ 40 red metallic
☐ 45 RAL 6018 yellow green met.
☐ 46 RAL 9023 pearl dark grey
☐ 53 dormant cassis
☐ 58 satin black
☐ 59 atlantic blue
☐ 60 mojito green

met. = metallized

Because of the material are different colours possible

Edition

"Duo-Color"

Colour 1: _____

- Seat frame
- Footrest holder
- Footrest

Farbe 2: _____

- Frame
- Gas spring bar
- Back part down
- Hand bow

Code U705

☐ OE 60,- €

Seat tilt
-10 to 30°

Incline back
 with gas spring
90 to 120°

Weight:


from 12 kg (without rear wheel)
 from 16 kg (with rear wheel 24")

min. seat height
with rear wheel


12" a. 20" **42 cm**
 22" **43 cm**
 24" **44 cm**

concerning seat frame
(trapezium adapter +2cm)seat can be changed
9cm in height






H total height: 95 - 105cm*
 L total length: 79,5 - 92cm
 T seat depth: size 0 = 32cm
 size 1-3 = 35cm







* according to the angle of the back / hand bow and the size of the rear wheels






frame accessories

Tip shackle with adapter  ☐ right ☐ left Code U711 ☐ Piece 66,- €	Anti-tip wheel with adapter  ☐ right ☐ left Code U712 ☐ Piece 118,- €	Armrest holder with arm pads, fixed at the back tube  (arm pads swing up) Code U746 ☐ Pair 240,- €	Armrest adjustable in height with raster, for mounting on the seating shell  ☐ armrest long ☐ armrest short Code U748 (self installation on the seating shell) ☐ Pair 176,- €	Armrest short adjustable in height and angle adjustable with raster, for mounting on the seating shell  Code U747 (self installation on the seating shell) ☐ Pair 220,- €
--	--	--	--	---




footrest

Single footrest with knee angle adjustable  Code U722 ☐ Piece 244,- €	Single footrest with knee angle and gear segment adjustable  Code U721 ☐ Piece 302,- €	Footrest angle adjustable additional joint for angle adjustable (for single footrest with knee angle CODE U722 or U721)  Code U728 ☐ Pc. OE 85,- €	Footrest folding upwards (for single footrest with knee angle CODE U722 or U721)  Code U725 (not in combination with CODE U728) ☐ Pc. OE 145,- €	Single footrest alu-profile, with holder at the seat tube angle adjustable, folding upwards and LLL adjustable  LLL from upper edge seat frame to footrest ☐ LLL 15 - 32cm (mounted in frame) ☐ LLL 30 - 40cm* (Footrest holder long with footrest bow) Code U726 * Note: only possible at a very high seat height (because of the castor) ☐ Piece 302,- €	Footrest holder height adjustable only in combination with CODE U726 or U727  Code U729 ☐ Pair OE 150,- €
--	---	---	--	---	---

footrest





Spasticity lock for footrest folding upwards for CODE U725 ☐ Code U725S single footrest for CODE U726 ☐ Code U5024 double footrest for CODE U727 ☐ Code U723 OE 60,- € OE 68,- € OE 120,- €	Double footrest alu-profile, with holder on the seat tube angle adjustable, folding upwards and LLL adjustable  LLL from upper edge seat frame to footrest ☐ LLL 15 - 32cm (mounted in frame) ☐ LLL 30 - 40cm* (Footrest holder long with footrest bow) Code U727 * Note: only possible at a very high seat height (because of the castor) ☐ Pair 322,- €	Single footrest alu. flat-profile gear segment knee angle  Code U7025 ☐ Piece 308,- €	Double footrest alu. flat-profile gear segment knee angle  Code U7026 ☐ Pair 390,- €	Heel edge for double footrest for CODE U7026 ☐ Pair Code U7028 for single footrest for CODE U7025 ☐ Piece Code U7027 each 104,- €	Calves pads for double footrest for CODE U7026 ☐ Pair Code U7023 for single footrest for CODE U7025, U721, U722 ☐ Piece Code U7022 each 110,- €
--	--	--	---	---	---

back

Back with hand bow height adjustable, angle adjustable and tiltable ☐ Hand bow U-form (standard) ☐ Hand bow trapezium-form (In comb. with U741 not possible) Code U741 (standard) ☐ with one hand brake Code U744 ☐ with two single brake lever  Back angle: Code U719 ☒ Gas spring Code U741-S	Back with single push handles, height adjustable (2 single brake lever)  ☐ Code U745	Back tubes 10cm higher ☐ Code U742	Back guide in the centre with fork for seating shell  Code U742	Back guide for seating shells with slide bush on back tube  ☐ Code U742R ☐ Code U742RW with CODE U741-S not possible	Back angular adjustment fixed without gas spring ☐ Code U743
Standard	OE 80,- €	☐ Pc. 220,- €	Pc. 224,- €	Pc. 248,- €	without OE

adapter


seat

Trapezium adapter 2-pieces (for seating shell)  Code U731 ☐ Set 264,- €	Trapezium adapter upper part (for seating shell)  Code U732 ☐ Piece 107,- €	Trapezium adapter alu (for mobility base)  Code U733 ☐ Piece 157,- €	HC-Trapezium adapter for seating shell, 2-pieces alu  Code U739 ☐ Set 345,- €	Seat 10cm longer seat frame long (forwards)  Code U734 ☐ Set OE 152,- €	Seat extension 7cm (backwards) Not with Scalomobil  Code U735 ☐ Set OE 144,- €	Seat frame with 2 gas spring for seat angle adjustment ## Special option## for size 2 or 3 Code U708 ☐ Set OE 472,- €
--	--	---	--	--	---	---

Sum

Footrest swing away from the wheelchair SPRINT is in special construction possible II

VD420E-27 Technical changes and print mistakes are reserved.

rear wheel					
Rear wheel with drum brake, hand rims alu., quick release axle and pneumatic tyre (for companion) ☐ 12 inch Code U761 without handrim / with screw axle ☐ 20 inch Code U762 only with 140mm castor possible ☐ 22 inch Code U763 ☐ 24 inch Code U764	Rear wheel with drum brake, handrim alu., quick release axle and solid tyre (for companion) ☐ 20 inch Code U767 ☐ 22 inch Code U768 ☐ 24 inch Code U769	Camber for rear wheel 20", 22" a. 24" ☐ 0° standard ☐ 4° * ☐ 7° * * not possible: - attachment for scalamobil - wheel base long	Wheel cover plastic (with holder)  (for rear wheel 20", 22" or 24") Code U779	Wheel cover aluminium round, coloured coated (with holder)  (for rear wheel 20", 22" or 24") Code U776	Wheel cover PC-Carbon (Clothing guard round with holder) (for rear wheel OFF-ROAD 24")  Code U776PC
Standard			☐ Pair 220,- €	☐ Pair 312,- €	☐ Pair 242,- €

rear wheel				
SAFETY-GRIP wheel rim with integrated hand rim ☐ 22 inch Code U7602 ☐ 24 inch Code U7604 only in combination with rear wheel CODE U763, U764, U768 or U769	"Black Edition" Rear wheel with black hub, black rim and quick release axle with black cap and black tyre Marathon Plus Evolution ☐ 22"x1 pneumatic Code U763S ☐ 24"x1 pneumatic Code U764S	Spoke cover transparent  ☐ 20 inch Code U770 ☐ 22 inch Code U771 ☐ 24 inch Code U772	Handrim cover black  ☐ 20" Code U794 ☐ 22" Code U795 ☐ 24" Code U796	Wheel base long (extension by 9cm) (In combination with Scalamobil not possible) Code U777
☐ Pair OE 60,- €	☐ Pair OE 95,- €	☐ Pair 156,- €	☐ Pair 102,- €	☐ Pair 134,- €

rear wheel	accessories	
Rear wheel 24" with OFF-ROAD tyre (with special size of wheel 507) with aluminium hand rim black instead of normal rear wheels -> only in combination with U764 II!!! ☐ BIG BETTY 24" x 2,40 (62-507) Mountainbike Profil ☐ BIG APPLE 24" x 2,35 (60-507) fine profile (Ballon-tyre) Code U759 Note: only with 0° camber possible II	Rear wheel OFF-ROAD 24" additional (special size of wheel 507) with aluminium hand rim black and quick release axle additional to the normal rear wheel 24" (to change the wheel) ☐ Code U759BB BIG BETTY 24" x 2,40 (62-507) Mountainbike profile ☐ Code U759BA BIG APPLE 24" x 2,35 (60-507) fine profile (Ballon tyre) Note: only with 0° camber possible II	Anchor plate system anchor plates for transport in cars (ISO 7176-19) Assembly kit for mobility base systems including mounting  CRASH TESTED Code U7091 Note: not in combination with E-Fix ☐ Set net price 430,- €
☐ Pair OE 288,- €	☐ Pair 688,- €	

castor	accessories	brake
Castor with screw axle ☐ 140mm black soft (solid) Code U781S ☐ 140mm soft (solid) Code U781 ☐ 6 Zoll solid * Code U784 ☐ 6 Zoll pneumatic * Code U783 ☐ 7 Zoll pneumatic * Code U782 ☐ 7 Zoll solid * Code U788 * 6 inch or 7 inch not by 20" and 22" rear wheel Standard	Attachment for scalamobil only by 24" rear wheel # from size 1 Code U709 (Combination with tip shackle, anti-tip wheel and seat 7cm longer to the back is not possible!!) ☐ Set 260,- €	Fixing brake for rear wheel  ☐ 20 inch Code U791 ☐ 22 inch Code U792 ☐ 24 inch Code U793 ☐ Pair OE 132,- € Brake lever extension straight, foldaway  Code U778 ☐ Pair 50,- €

F1	Mobility base for seating systems	basic price	€ 2120,-
For all orders apply the general terms and conditions.		Sum optinal extras	€
		Shipping and handling	€
		Sum	€
		Tax	€
Purchaser _____ Telephone _____		Total price	€

Attention: Scalamobil-extension:
Wheel base long, anti-tip wheel, tip shackle and seat 7cm longer to the back is not possible !!!