

# NITRO

## Rollator Range

*"This stylish range of lightweight Nitro rollators continues to evolve!"*



Nitro Standard



Nitro Mini



Nitro Tall



Nitro HD

Following on from the success of the Nitro Standard Rollator, we decided to build on this by developing a full range of Nitro Rollators. The Nitro Tall was the first step in extending the range, followed by the introduction of the Nitro HD and finalised with the addition of the Nitro Mini.

[www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)

**drive**  
DeVilbiss<sup>®</sup>  
HEALTHCARE

The Nitro Rollators were developed with a contemporary design in mind, featuring a lightweight aluminium frame and stylish tri-spoke wheels for easy manoeuvrability.



Easily folds with one hand to a compact size for storage and portability

## Superb Comfort, Style & Performance...

- » Lightweight aluminium frame which has been designed for easy manoeuvrability
- » Suitable for use in tight spaces with a caster fork design that enhances turning radius
- » Internally routed brake cables for added safety
- » Push button height adjustment on handles to suit user requirements
- » Cross brace design for added stability
- » Easily folds with one hand to a compact size for storage and portability
- » Integral seat handle serves as a useful handgrip for pulling rollator when used as a shopper or for carrying when folded
- » Large 25cm (10") front casters for easy manoeuvrability
- » Lightweight walking frame is suitable for indoor and outdoor use
- » Removable zippered storage bag for storing shopping items and essentials
- » Plastic clip holds the rollator closed when frame is folded for easy handling
- » Back support height easily adjusts with tool free thumb screw
- » Soft, flexible backrest provides comfort and security while seated
- » Rear stepper tube facilitates safe curb-climbing for both attendant and independent user

### Nitro Mini

- » The Nitro Mini is ideal for users 5'2" and below

### Nitro Tall


- » The Nitro Tall is suitable for users 6'0" and above

### Nitro HD

- » The Nitro HD has a higher weight capacity of 32st (203kg)



Removable zippered storage bag for storing essentials or shopping



Soft, flexible backrest and seat providing great comfort

### Colour Options:



Black (Standard and Heavy Duty)




Red (Standard, Tall and Mini)



Rear stepper tube facilitates safe curb-climbing



Stylish tri-spoke wheels



Integral suspension spring provides weight capacity of 32st (HD only)



Internally routed brake cables for added safety

CODE	DESCRIPTION	WIDTH	DEPTH	HEIGHT MIN/MAX	WEIGHT	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	WEIGHT CAPACITY
RTL10266	Nitro Standard	59cm (23")	48cm (19")	85cm (34") / 97cm (38")	8kg (18lb)	46cm (18")	20cm (8")	52cm (21")	136kg (21st)
RTL10266-H	Nitro Mini	59cm (23")	70cm (28")	71cm (28") / 80cm (31")	8kg (18lb)	46cm (18")	21.5cm (8.5")	46cm (18")	136kg (21st)
RTL10266-T	Nitro Tall	59cm (23")	70cm (28")	91cm (36") / 104cm (41")	8kg (18lb)	46cm (18")	21.5cm (8.5")	60cm (24")	136kg (21st)
RTL10266HD-BLK	Nitro HD	70cm (28")	70cm (28")	91cm (36") / 104cm (41")	13kg (28.5lb)	58cm (23")	26cm (10")	60cm (24")	203kg (32st)

Distributed by: DeVilbiss Healthcare GmbH, Kamenzer Straße 3, 68309 Mannheim, Germany • Tel: +49 (0)621 187 98 230 • [www.drivedevilbiss-int.com](http://www.drivedevilbiss-int.com)

Leaflet Code: LL231 | If you would like to view this product please contact your local stockist.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.\* Depends on user mass, terrain and battery condition.