

Start^{M1} Intro Lightweight Wheelchair

	Order no.:
--	---------------------



☐ Provide estimate only ☐ Order	
Invoice address / Customer number	
Delivery address	
Sales person	Name Date Signature

Standard Equipment

- Aluminum tube frame
- Tip-assist integrated into the frame
- Individual footrests, angle adjustable, plastic footplate
- Swing-away footrest holders
- Desk side panel, different height settings possible; with armrest
- Seat and back upholstery of black nylon
- Modular design
- Rear wheels with quick-release axles
- Wheel lock
- Frame color: silver metallic
- 24" rear wheel with PU tire

Factory presettings:

- Seat depth: 43cm
- Back height: 42.5 cm
- Front seat height: 43-47cm, seat inclination: 3cm

Folding: Max. user weight 125kg

Variant		Price
Start^{M1} Intro Standard Chair		
480F53=10501	[] 40.5 Standard Intro - 24" Rear Wheels & 6" Casters	Add Charge
480F53=10503	[] 45.5 Standard Intro - 24" Rear Wheels & 6" Casters	Add Charge
480F53=10505	[] 50.5 Standard Intro - 24" Rear Wheels & 6" Casters	Add Charge
Start^{M1} Intro with Reclining Back (Special Order)		
480F53=10301	[] 40.5 Intro with Reclining Back - 24" Rear Wheels & 6" Casters	Add Charge
480F53=10303	[] 45.5 Intro with Reclining Back - 24" Rear Wheels & 6" Casters	Add Charge
480F53=10305	[] 50.5 Intro with Reclining Back - 24" Rear Wheels & 6" Casters	Add Charge
* note that the anchor point transportation kit is not transport approved for the reclining back model		
Start^{M1} Intro		
Frame Accessories - To Be Fitted by Dealer		
481A53=SK170	[] Plug-on Anti-Tipper - Right	Add Charge
481A53=SK170	[] Plug-on Anti-Tipper - Left	Add Charge
481A53=SK065	[] Crutch Holder - Right	Add Charge
481A53=SK065	[] Crutch Holder - Left	Add Charge
Non Standard Casters		
481F53=PK018	[] 8" Solid Casters	Add Charge

	Brakes	
481H53=PK056	[] Dual Acting Wheel Lock for User and Attendant	Add Charge
	Wheel Lock Accessories	
481H53=SK012	[] Wheel lock lever extensions [] 1 Pc	Add Charge
	Wheel lock lever extensions [] 2 Pc	Add Charge
	Push Handles	
481D53=SV042	[] Height Adjustable Push Handles	Add Charge
	Back Accessories	
481D53=SK139	[] Adjust Back Upholstery 40.5cm - (Not for Reclining Back Version)	Add Charge
481D53=SK093	[] Adjust Back Upholstery 45.5cm - (Not for Reclining Back Version)	Add Charge
481D53=SK095	[] Adjust Back Upholstery 50.5cm - (Not for Reclining Back Version)	Add Charge
	Side Panels	
481E53=PK003	[] Height Adjustable Side Panel " T Style" (10cm Adjust)	Add Charge
481E53=PK100	[] Desk Side Panel - Stepped	Add Charge
481E53=ST092	[] Longer Armrest Pads 365mm [] 1 Pc	Add Charge
	Longer Armrest Pads 365mm [] 2 Pc	Add Charge
	Side Panel Accessories	
481E53=ST051	[] Tray - SW 40.5cm	Add Charge
481E53=ST053	[] Tray - SW 45.5cm	Add Charge
481E53=ST055	[] Tray - SW 50.5cm	Add Charge
	Footplate Assembly	
481B53=LK100	[] Elevating Footrest Left - Hanger Only	Add Charge
481B53=RK100	[] Elevating Footrest Right - Hanger Only	Add Charge
481B53=PK135	[] Elevating Footrest - SW 40.5 - (Includes Footplates)	Add Charge
481B53=PK140	[] Elevating Footrest - SW 45.5 - (Includes Footplates)	Add Charge
481B53=PK145	[] Elevating Footrest - SW 50.5 - (Includes Footplates)	Add Charge
481B53=LK128	[] Stump Support - Left	Add Charge
481B53=RK128	[] Stump Support - Right	Add Charge
	Spoke Protectors - 24"	
481G00=SK702	[] Spoke Protector Transparent	Add Charge
481G00=SK747	[] Spoke Protector Eye	Add Charge
481G00=SK743	[] Spoke Protector Fire	Add Charge
481G00=SK749	[] Spoke Protector Flower Power	Add Charge
481G00=SK751	[] Spoke Protector Clown	Add Charge
	Other Accessories	
481C08=ST010	[] Tool kit	Add Charge
481S00=SK030	[] Anchor Point Kit for Transportation (Max Load 100 Kg)	Add Charge
481G53=SK071	[] 24" One Arm Drive Set - PU Tyres	Add Charge
Start ^{M1} Intro		
	Seat Cushion	
476C00=SK420	[] Combi 3S Seat Cushion (SW x SD) equivalent to Chair	Add Charge
481C53=SM155	[] Basic Cushion with Incontinence Cover - 40 x 42 x 4.5cm	Add Charge
481C53=SM165	[] Basic Cushion with Incontinence Cover - 42 x 45 x 4.5cm	Add Charge
481C53=SM180	[] Basic Cushion with Incontinence Cover - 49 x 42 x 4.5cm	Add Charge
476C00=SK100	[] CubicFoam Seat Cushion - Required Seat Depth	Add Charge
476C00=SK200	[] CombiFoam Seat Cushion - Required Seat Depth	Add Charge
	New Cushions	
476C00=SK601	[] Terra, 42x42	Add Charge
476C00=SK602	[] Terra, 46x46	Add Charge
476C00=SK603	[] Terra, 46x42	Add Charge
476C00=SK604	[] Terra, 42x46	Add Charge
476C00=SK605	[] Terra, 50x46	Add Charge
476C00=SK606	[] Terra, 38x42	Add Charge
476C00=SK607	[] Terra, 44x46	Add Charge
476C00=SK651	[] Aquos, 42x42	Add Charge
476C00=SK652	[] Aquos, 46x46	Add Charge
476C00=SK653	[] Aquos, 46x42	Add Charge

Anterior seat height-posterior seat height<=3.5 cm

Casters in dependency of the anterior seat height

Caster fork in dependency of the anterior seat height and casters

Rear wheel size in dependency of the posterior seat height

Start M1, M2, M6

Posterior seat height (in cm)

Posterior seat height		Setting range
Rear wheel size	20"	35-45
	22"	37-47
	24"	43-50

Start M4, M5

Posterior seat height (in cm)

Free mounting positions	Rear wheel size	
	22"	24"
1	38	41
2	41	44
3	44	47
4	47	50

With the option "without wider foot space", the anterior seat height values of the M5 apply.

Anterior seat height-posterior seat height<=3.5 cm


Casters in dependency of the anterior seat height

Caster fork in dependency of the anterior seat height and casters

Rear wheel size in dependency of the posterior seat height

Rear wheel adapter

M1, M2, M3, M6



M4, M5

