Invacare® Action®3 NG

Invacare® Action®3 NG
Folding Aluminium Manual Wheelchair



COMPATIBLE alber

Bodypoint

Description:

- Weight: from 14.2 kg.
- Maximum user weight: 125 kg.
- 2-year warranty period excluding wear parts

Date :

Purchase Order :

Customer Code :

Distributor Stamp:

DHCare Unique Customization Service
Please specifiy the quote or special production number :

❑ Please make your selection with a cross in the box ◆ Standard feature ❍ Free of charge option

CODE		DESCRIPTION	NC : Non Compatible NC : 9999 = not compatible with the reference 9999	Euros	
CEB0031	❐	*Invacare® Action®3 NG* - Fixed Backrest			
CEB0036	❐	*Invacare® Action®3 NG* - Folding Backrest			
CEB0041	❐	*Invacare® Action®3 NG* - Reclining Backrest			
FLEXIBLE POSITIONING FILE PACK					
On demand	❐	Asymmetrically adjustable backrest with 8 straps + 2 repositionable trunk supports + **Matrx Libra** seat cushion + a rigid base + a **Bodypoint** 2-point pelvic belt	**Only with 43 cm backrest height**		
ALBER					
CEB0070	❐	Bracket For **Twion / E-Motion M25** - Including Alber Anti-Tipper	NC : 0310, 0315, 0732, 0832, 1064, 1066, 1073, 1074, 1082, 1170, 1185, 1201, 1280, 1320, 1330, 1345, 1365, 1370, 1375, 1430, 1440		
CEB0080	❐	Bracket For **Viamobil V25**	NC : 0300, 0310, 0315, 0455, 0460, 1170, 1201, 1365, 1370, 1375, 1430, 1440		
CEB0090	❐	Bracket For **Viamobil Eco V14**	NC : 0290, 0300, 0310, 0315, 0905, 0906, 0907, 0945, 1185, 1201, 1365, 1370, 1375, 1430, 1440, 1451		
PROPULSION TYPE					
CEB0100	❐	Self Propulsion		◆	
SEAT					
SEAT WIDTH					
CEB0180	❐	380 mm	NC : 0545+0595, 0852, 0855	❍	
CEB0190	❐	405 mm	NC : 0545+0595	❍	
CEB0200	❐	430 mm	NC : 0853	◆	
CEB0210	❐	455 mm	NC : 0853	❍	
CEB0220	❐	480 mm	NC : 0853, 0854	❍	
CEB0230	❐	505 mm	NC : 0853, 0854, 0855, 0856	❍	
SEAT DEPTH					
CEB0240	❐	400 mm	NC : 0879	◆	
CEB0250	❐	450 mm	NC : 0853, 0878	❍	
CEB0252	❐	Adjustable from 450 to 400 mm Delivered at 430 mm (or 450 mm if choice of a seat cushion)	NC : 0330, 0350	31,24	

SEAT HEIGHT				
CEB0260	❒	510 mm	NC : 0738, 0838, 1035, 1050, 1054, 1056, 1066	❍
CEB0270	❒	485 mm	NC : 1054	◆
CEB0280	❒	460 mm		❍
CEB0290	❒	435 mm	NC : 0660/0670/0760/0770 + 24" wheels, 1005, 1010, 1030	❍
CEB0300	❒	410 mm	NC : 24" wheels, 0935, 1005, 1010, 1030, 1054, 1056, 1365, 1370, 1440	❍
CEB0310	❒	385 mm (only with 5" front wheels CEB1054 and 22" rear wheels)	NC : 0660, 0670, 0738, 0760, 0770, 0838, 0935, 1035, 1056, 1068, 1280, 1365, 1370, 1440	❍
CEB0315	❒	360 mm (only with CEB1054 5'' front wheels and 20" rear wheels)	NC : 0600, 0611, 0620, 0660, 0670, 0738, 0760, 0770, 0838, 0935, 1280, 1035, 1056, 1066, 1068, 1365, 1375, 1440	
SEAT UPHOLSTERY TYPE				
CEB0330	❒	NYLON WITHOUT LOOP BANDS SEAT UPHOLSTERY		◆
CEB0340	❒	NYLON WITH LOOP BANDS SEAT UPHOLSTERY		❍
CEB0350	❒	PADDED NYLON SEAT UPHOLSTERY		
BACKREST				
Backrest type - anti-tippers recommended with reclining backrests				
CEB0380	❒	FIXED 0° BACKREST 2 HEIGHTS 400/430 MM (DELIV. 430 MM) - FIXED PUSH HANDLES	NC : 0471, 0945	
CEB0390	❒	FIXED 0° BACKREST 510 MM HEIGHT - FIXED PUSH HANDLES	NC : 0471, 0945	
CEB0395	❒	FIXED 7° BACKREST - HEIGHT ADJUST 380 TO 460 MM (DELIVERED 380 MM)	NC : 0474, 0475, 0478, 0479, 0481, 0905, 0906, 0907, 0940, 1451	
CEB0396	❒	FIXED 7° BACKREST - HEIGHT ADJUST 380 TO 460 MM (DELIVERED 380 MM)	NC : 0474, 0475, 0478, 0479, 0481, 0905, 0906, 0907, 0940, 1451	
CEB0410	❒	FIXED 7° BACKREST - 2 HEIGHTS 400/430 MM (DELIV. 430 MM)	NC : 0471, 0945	◆
CEB0420	❒	FIXED 7° BACKREST - 510 MM HEIGHT	NC : 0471, 0945	
CEB0435	❒	RECLINING (RATCHET) 0/30° BACKREST INCL. BACKRACE & PASSIVE REAR WHEEL POS. ONLY - 3 HEIGHTS 390/410/430 MM (DELIV.430 MM)	NC : 0471, 0738, 0940, 0945, 1185, 1190	◆
CEB0445	❒	RECLINING (GAS SPRING) 0/30° BACKREST INCL. BACKRACE & PASSIVE REAR WHEEL POS. ONLY - 3 HEIGHTS 390/410/430 MM (DELIV.430 MM)	NC : 0471, 0738, 0940, 0945, 1185, 1190	❍
CEB0450	❒	RECLINING (GAS SPRING) 0/30° BACKREST INCL. BACKRACE & PASSIVE REAR WHEEL POS. ONLY - 3 HEIGHTS470/490/510 MM (DELIV.510 MM)	NC : 0471, 0738, 0940, 0945, 1185, 1190	
CEB0455	❒	FOLDING 7° BACKREST - 430 MM HEIGHT	NC : 0471, 0945	◆
CEB0460	❒	FOLDING 7° BACKREST - 510 MM HEIGHT	NC : 0471, 0480, 0481, 0945	
Matrx backrests - See the table "Possible Matrx backrest widths depending on the seat width" at the bottom of the measurement sheet. Find the trunk wedges and all the accessories on the Matrx measurement sheet!				
CEB0493	❒	VICAIR LIBERTY BACKREST - 450 MM BACK HEIGHT INCLUDING BACKBRACE *(Only for authorized countries)*	NC : 0471-0481, 0940, 0945, 0952, 1345	
CEB0451	❒	**Matrx Elite E2** height 160 mm - Curve 80 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0452	❒	**Matrx Elite E2** height 210 mm - Curve 80 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0453	❒	**Matrx Elite E2** height 260 mm - Curve 80 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0454	❒	**Matrx Elite E2** height 310 mm - Curve 80 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0456	❒	**Matrx Elite E2** height 360 mm - Curve 80 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0461	❒	**Matrx Elite E2** height 410 mm - Curve 80 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0462	❒	**Matrx Elite E2** height 460 mm - Curve 80 mm - **Standard** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0457	❒	**Matrx Elite E2** height 510 mm - Curve 80 mm - **Standard** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	



0396




0435




0455 / 0460

CEB0458	❒	**Matrx Elite Deep E2** height 160 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0459	❒	**Matrx Elite Deep E2** height 210 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0466	❒	**Matrx Elite Deep E2** height 260 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0467	❒	**Matrx Elite Deep E2** height 310 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0468	❒	**Matrx Elite Deep E2** height 360 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0463	❒	**Matrx Elite Deep E2** height 410 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0464	❒	**Matrx Elite Deep E2** height 460 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
CEB0469	❒	**Matrx Elite Deep E2** height 510 mm - Curve 160 mm - **Mini** fasteners	NC : 0471, 0474, 0477, 0478, 0481, 0882, 0940, 0945, 1481	
Caution *- loss of approx. 20 mm of seat depth with tiltable backrests - supplied with Mini or Standard fixings except for PMH backrests supplied with Split fixing.*				
Matrx backrest width- see table at the end of the measurement sheet				
CEB0482	❒	360 mm - indicative bust width 265 mm	NC : 1481	❍
CEB0483	❒	380 mm - indicative bust width 290 mm	NC : 1481	❍
CEB0484	❒	410 mm - indicative bust width 315 mm	NC : 1481	❍
CEB0485	❒	430 mm - indicative bust width 340 mm	NC : 1481	❍
CEB0486	❒	460 mm - indicative bust width 370 mm	NC : 1481	❍
CEB0487	❒	480 mm - indicative bust width 390 mm	NC : 1481	❍
CEB0488	❒	510 mm - indicative bust width 420 mm	NC : 1481	❍
BACKREST UPHOLSTERY TYPE				
CEB0471	❒	STANDARD NYLON SLING TYPE (JEMIMA)	NC : 0882, 0905, 0906, 0907, 0882	❍
CEB0474	❒	LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER INCL. BACK POCKET		
CEB0475	❒	LUMBAR SUPPORT (+1 STRAP) PADDED SHAPED COVER		◆
CEB0478	❒	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL. BACK POCKET	NC : 0940, 0945	
CEB0479	❒	CONTOUR ADJUSTABLE (3+1 STRAPS) PADDED SHAPED COVER	NC : 0940, 0945	
CEB0480	❒	CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER	NC : 0940, 0945	
CEB0481	❒	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL. BACK POCKET	NC : 0940, 0945	



0471




0474 / 0475



0480 / 0481

HANGERS/FOOTPLATES

Separate plates delivered with heel straps				
HANGERS RH TYPE				
CEB0510	❒	80° STRAIGHT RH SWING IN/SWING OUT HANGER		◆
CEB0520	❒	ELEVATING LEGREST (0° TO 75°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD)	NC : 0611, 0930, 0935, 0936	
CEB0545	❒	AMPUTEE LEGREST RH	NC : 0600, 0610, 0611, 0930, 0936	
HANGERS LH TYPE				
CEB0560	❒	80° STRAIGHT LH SWING IN/SWING OUT HANGER		◆
CEB0570	❒	ELEVATING LEGREST (0° TO 75°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD)	NC: 0611, 0930, 0935, 0936	
CEB0595	❒	AMPUTEE LEGREST LH	NC : 0611, 0620, 0630, 0930, 0936	
FOOTPLATES RH TYPE				
CEB0600	❒	FLIP UP RH FOOTPLATE - (INCLUDING HEEL STRAP)		◆
CEB0610	❒	FLIP UP DEPTH & ANGLE AJUSTABLE RH FOOTPLATE - (INCLUDING HEEL STRAP)		
CEB0611	❒	FLIP UP ONE PIECE FOOTPLATE - DEPTH & ANGLE ADJUSTABLE	NC : 0620, 0630, 0935, 0936	
FOOTPLATES LH TYPE				
CEB0620	❒	FLIP UP LH FOOTPLATE - (INCLUDING HEEL STRAP)		◆
CEB0630	❒	FLIP UP DEPTH & ANGLE AJUSTABLE LH FOOTPLATE - (INCLUDING HEEL STRAP)		



0510 / 0560



0520 / 0570



0545 / 0595




0610 / 0630



0610 / 0630



0611

ARMRESTS					
ARMRESTS RH FORM IN BLACK COLOR					
CEB0640	❒	REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED - DESK	NC : 0760, 0770, 0772, 0780, 0790, 0795, 0796, 0798, 0838, 1480	◆	 0640 / 0740
CEB0650	❒	REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED - FULL LENGTH	NC : 0760, 0770, 0772, 0780, 0790, 0795, 0796, 0798, 0838		
CEB0656	❒	REMOVABLE AND CONTINUOUSLY HEIGHT-ADJUSTABLE FROM 200 MM TO 280 MM - SHORT NARROW GUTTER (330 MM) WITHOUT ROTATION	NC : 0760, 0770, 0772, 0780, 0790, 0795, 0796, 0798, 0838, 1480		 0656 / 0756
CEB0660	❒	SWING AWAY DETACHABLE RH (N°2) - DESK	NC : 0740, 0750, 0756, 0772, 0780, 0790, 0795, 0796, 0838, 1480	❍	
CEB0670	❒	SWING AWAY DETACHABLE RH (N°2) - FULL LENGTH	NC : 0740, 0750, 0756, 0772, 0780, 0790, 0795, 0796, 0798, 0838		 0660 / 0760
CEB0672	❒	LIFTABLE AND REMOVABLE, HEIGHT-ADJUSTABLE FROM 170 TO 295 MM - SHORT (330 MM) NON-ROTATING NARROW TROUGH WITH CLOTHES GUARD AND ERGONOMIC SUPPORT HANDLE	NC : 0740, 0750, 0760, 0770, 0795, 0796, 0798, 0838, 1480		 0672 / 0772
CEB0680	❒	SWING AWAY DETACHABLE RH (N°3) - DESK	NC : 0740, 0750, 0756, 0760, 0770, 0795, 0796, 0798, 0838, 1480		
CEB0690	❒	SWING AWAY DETACHABLE RH (N°3) - FULL LENGTH	NC : 0740, 0750, 0756, 0760, 0770, 0795, 0796, 0832, 0838		 0680 / 0780
CEB0695	❒	REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED - DESK	NC: 0740, 0750, 0756, 0760, 0770, 0772, 0780, 0790, 0838, 1185, 1480		
CEB0696	❒	REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED - FULL LENGTH	NC : 0740, 0750, 0756, 0760, 0770, 0772, 0780, 0790, 0838, 1185		 0695 / 0795
CEB0698	❒	REMOVABLE AND HEIGHT ADJUSTABLE - SHORT NARROW GUTTER (330 MM) WITHOUT ROTATION	NC : 0740, 0750, 0756, 0760, 0770, 0772, 0780, 0790, 0838, 1185, 1480		
CEB0732	❒	REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED - ARMPAD HEMI BLACK PU TR26	NC : 0838, 1185, 1480		 0698 / 0798
CEB0738	❒	REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 & 435 MM) -	NC : 0740, 0750, 0756, 0760, 0770, 0772, 0780, 0790, 0795, 0796, 0832, 1064, 1066, 1200, 1201, 1480		
ARMRESTS LH FORM IN BLACK COLOR					
CEB0740	❒	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED - DESK	NC : 1480	◆	
CEB0750	❒	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED - FULL LENGTH			
CEB0756	❒	REMOVABLE AND CONTINUOUSLY HEIGHT-ADJUSTABLE FROM 200 MM TO 280 MM - SHORT NARROW GUTTER (330 MM) WITHOUT ROTATION	NC : 1480		
CEB0760	❒	SWING AWAY DETACHABLE LH (N°2) - DESK	NC : 1480	❍	 0732 / 0832
CEB0770	❒	SWING AWAY DETACHABLE LH (N°2) - FULL LENGTH			
CEB0772	❒	LIFTABLE AND REMOVABLE, HEIGHT-ADJUSTABLE FROM 170 TO 295 MM - SHORT (330 MM) NON-ROTATING NARROW TROUGH WITH CLOTHES GUARD AND ERGONOMIC SUPPORT HANDLE	NC : 1480		
CEB0780	❒	SWING AWAY DETACHABLE LH (N°3) - DESK	NC : 1480		 0738 / 0838
CEB0790	❒	SWING AWAY DETACHABLE LH (N°3) - FULL LENGTH			
CEB0795	❒	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED - DESK	NC : 1185, 1480		
CEB0796	❒	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED - FULL LENGTH	NC : 1185		
CEB0798	❒	REMOVABLE AND HEIGHT ADJUSTABLE - SHORT NARROW GUTTER (330 MM) WITHOUT ROTATION	NC : 1480, 1185		
CEB0832	❒	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED - ARMPAD HEMI BLACK PU TR26	NC : 1185, 1480		
CEB0838	❒	REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 & 435 MM) -	NC : 1064, 1066		

SEAT CUSHION TYPE				
Cushion				
CEB0849	❒	BASIC 60 MM FOAM CUSHION - NYLON & ANTI-SLIP COVER	NC : 0853, 0854, 0855, 0856, 0873	
CEB0853	❒	CLASSIC VISCO CUSHION 375X 400 MM -	NC : 0873, 0878, 0879	
CEB0854	❒	CLASSIC VISCO CUSHION 405 X 400 MM -	NC : 0873, 0878, 0879	
CEB0855	❒	CLASSIC VISCO CUSHION 430 X 430 MM -	NC : 0873, 0878, 0879	
CEB0856	❒	CLASSIC VISCO CUSHION 455 X 400 MM -	NC : 0873, 0878, 0879	
CEB0869	❒	FLO-TECH LITE CUSHION THICKNESS 55 MM		
CEB0871	❒	MATRX CONTOUR NG VISCO CUSHION		
CEB0873	❒	FLO-TECH IMAGE CUSHION THICKNESS 55 MM		
CEB0874	❒	FLO-TECH SOLUTION XTRA CUSHION THICKNESS 55 MM		
CEB0876	❒	SEAT CUSHION MATRX LIBRA, BLACK MICRO-VENTILATED MESH COVER, THICKNESS 70 MM		
SEAT CUSHION DEPTH				
CEB0878	❒	FOR 400 MM SEAT DEPTH		❍
CEB0879	❒	FOR 450 MM SEAT DEPTH		❍
OTHER SEATING OPTIONS				
CEB0852	❒	SEAT RIGIDIZER	(Mandatory : choose a seat cushion)	
CEB0882	❒	LATERAL POSITIONING PADS (PAIR)		
CEB0905	❒	NECKREST BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	NC : 0906, 0907, 0940, 0945	
CEB0906	❒	HEADREST BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	NC : 0907, 0940, 0945	
CEB0907	❒	HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA) -	NC : 0940, 0945	
CEB0930	❒	ADJUSTABLE CALF STRAP		
CEB0935	❒	LONGER FOOTREST TUBE (UL= 450-530 MM)		
CEB0936	❒	SHORTER FOOTREST SUPPORT (PAIR) (UL= 230-430 MM)		
CEB0940	❒	BACKBRACE	NC : 0945	
CEB0945	❒	FOLDABLE BACKBRACE	NC : 1451	
COLOURS				
CEB0961	❒	SLATE GREY (RAL7024 FINE TEXTURE)		❍
CEB0964	❒	PEARL WHITE		❍
CEB0966	❒	MATT BLACK		❍
CEB0969	❒	RED MATT		❍
CEB0974	❒	NAVY BLUE		❍
CEB0980	❒	ANTHRACITE BLACK		❍
CEB0991	❒	SPECIAL COLOUR		On Demand
CASTER HOUSING TYPE				
CEB0995	❒	STANDARD		◆
FRONT WHEELS				
CEB1005	❒	200 X 32 MM (8" X 1.1/4") - BLACK PNEUMATIC TYRE		❍
CEB1010	❒	200 X 30 MM (8" X 1.1/4") - BLACK SOLID TYRE		❍
CEB1030	❒	200 X 45 MM (8" X 2") - BLACK SOLID TYRE		❍
CEB1035	❒	150 X 30 MM (6" X 1.1/4") - BLACK SOLID TYRE MEDIUM HARD		❍
CEB1050	❒	150 X 30 MM (6" X 1.1/4") - BLACK SOLID TYRE		◆
CEB1054	❒	125 X 30 MM (5" X 1") - BLACK SOLID TYRE		❍
CEB1056	❒	140 X 40 MM (5" 1/2 X 1" 3/4) - BLACK SOLID TYRE MEDIUM HARD		❍



0849



0853 / 0854




0873



0882




0905



0906 0907



0945



0961 / 0964 / 0966



0974 / 0969 / 0980



1010 1050


REAR WHEELS				
REAR WHEEL SIZE				
CEB1064	❒	20" / 515 MM	NC : 1071, 1073, 1074, 1088, 1170, 1260, 1280	❍
CEB1066	❒	22" / 560 MM	NC : 1074, 1260	❍
CEB1068	❒	24" / 600 MM	NC : 1261	◆
REAR WHEEL TYPE				
CEB1071	❒	STANDARD (BLACK HUB, COLORLESS RIM & ZINC PLATTED SPOKES)	NC : 1170, 1280	◆
CEB1072	❒	REINFORCED (COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES)	NC : 1280	❍
CEB1073	❒	DRUM BRAKE (20/22/24" : COLORLESS HUB, CLEAR ANODIZED RIM & STAINLESS STEEL SPOKES / 12" : BLACK PLASTIC HUB & RIM)	NC : 1270, 1275, 1320, 1330	❍
CEB1074	❒	COMPOSITE (BLACK PLASTIC HUB & RIM)	NC : 1088, 1260, 1261, 1280	❍
TYRE TYPE				
CEB1082	❒	TYRES SOLID BLACK 35 MM (1"3/8)		◆
CEB1084	❒	PNEUMATIC PROFILED BLACK 35 MM (1"3/8)		❍
CEB1088	❒	PNEUMATIC LOW PROFILE BLACK 25 MM (1")		❍
HANDRIM TYPE				
CEB1140	❒	ALUMINIUM PRE-ANODIZED HANDRIM		◆
CEB1150	❒	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC		
REAR WHEELS AXLE TYPE				
CEB1170	❒	FIXED AXLE	NC : 1185, 1365	❍
CEB1180	❒	QUICK RELEASE AXLE		◆
REAR WHEELS POSITION				
CEB1185	❒	Active position (37.5 mm axles in front of the backrest tubes) - Specific rear plates	NC : 1320, 1330, 1345, 1365	
CEB1190	❒	Standard position (25 mm axles in front of the backrest tubes)		◆
CEB1200	❒	Passive position (25 mm axles behind the backrest tubes - standard with reclining backrest)		❍
CEB1201	❒	Very rearward position (axles 60 mm behind the backrest tubes)	NC: 1275, 1280, 1320, 1330, 1365	
Pump				
CEB1230	❒	PUMP AND SUPPORT		
CHASSIS OPTIONS TYPE				
CEB1241	❒	FRONT WHEELS EXTENSION (+ 50 MM OUTWARD)		❍
CEB1260	❒	SPOKE GUARDS PAIR 600 MM (24") BLACK	NC : 1261	
CEB1261	❒	SPOKE GUARDS PAIR CLEAR		
BRAKES TYPE				
CEB1270	❒	PUSH TO LOCK BRAKES		◆
CEB1275	❒	PULL TO LOCK BRAKES		❍
CEB1280	❒	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)		
CEB1320	❒	OAD WIRE CONTROLED BRAKE RH	NC : 1390	
CEB1330	❒	OAD WIRE CONTROLED BRAKE LH	NC : 1390	
SAFETY OPTIONS TYPE				
CEB1345	❒	TRANSPORTATION BRACKETS		
CEB1360	❒	POSTURE BELT - PLASTIC BUCKLE	NC : 1361, 1362, 1363, 1364	
CEB1361	❒	HIP BELT, 2-PT, CTR-PULL, PB - M X 460MM PADS, FLAT-MOUNT	NC : 1362, 1363, 1364	
CEB1362	❒	HIP BELT, 2-PT, CTR-PULL, PB - L X 620MM PADS, FLAT-MOUNT	NC : 1363, 1364	
CEB1363	❒	HIP BELT, 4-PT, CTR-PULL, PB - M X 460MM PADS, FLAT-MOUNT	NC : 1364	
CEB1364	❒	HIP BELT, 4-PT, CTR-PULL, PB - L X 620MM PADS, FLAT-MOUNT		
CEB1365	❒	TRANSFER KIT (100 MM (4") WHEELS) - ANTI-BARRIER		
CEB1370	❒	REAR ANTI-TIPPERS PAIR	NC : 1375, 1440	
CEB1375	❒	REAR ANTI-TIPPERS PAIR WITH TIPPING AID - (REA TYPE)	NC : 1440	
CEB1380	❒	REFLECTORS KIT		
CEB1390	❒	BRAKE EXTENSION PAIR		



Standard



Composite



1" 1.3/8" 1.3/8"

1088 1084 1082



1241



1270



1345



1365



1370




1375



1361 / 1362



1363 / 1364

OPTIONS					
CEB1430	❒	CANE HOLDER			 1440  1480
CEB1440	❒	TIPPING AID RH			
CEB1451	❒	DRIP STAND SUPPORT WITH DRIPSTAND (WITH BACKBRACE 0940 ONLY)			
OTHERS OPTIONS					
CEB1480	❒	ACRYLIC TRAY			 1483/ 1484
CEB1480	❒	ACRYLIC TRAY			
CEB1480	❒	ACRYLIC TRAY			
ACCESSORIES (delivered separtly)					
AP1493219	❒	GUTTER WITH INDEPENDENT SIDE ROTATION, ADJUSTABLE IN HEIGHT AND DEPTH, RIGHT	Only compatible with armrest 0680, 0690, 0695, 0696, 0732		
AP1493220	❒	GUTTER WITH INDEPENDENT SIDE ROTATION, ADJUSTABLE IN HEIGHT AND DEPTH, LEFT	Only compatible with armrest 0780, 0790, 0795, 0796, 0832		
RAMPS					
1654128	❒	**Basic ramp 2 telescopic parts** in aluminum and non-slip coating. Folded dimensions: 1210 x 194 mm 8.8 kg per pair - maximum capacity: 270 kg			 1654128
1654129	❒	**Quick2go Compact 3-part aluminum telescopic ramp** with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg			 1654129
1654130	❒	**Aluminium folding mono ramp** Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg			 1654130