LIW CARE TECHNOLOGY
ORDER FORM

ORTHOPEDIC SEATS
LIWCARE.COM

BAFFIN NEOSIT HL

Baffin neoSIT HL is a multi-adjustable orthopedic seat. Available in 3 sizes: S, M, L and with with the chassis frame HL21 with the wheels of 100 mm diameter or chassis frame HL21 with the wheels of 150 mm front and 310 mm rear diameters.



○ **BAFFIN neoSIT HL INDOOR Seat with the chassis frame HL21 with the wheels of 100 mm diameter**



○ **BAFFIN neoSIT HL OUTDOOR Seat with the chassis frame HL21 with the wheels of 150 mm front and 310 mm rear diameters**


LIW's order number:	Date:
Distributor name:	Customer's order number:
Invoice information: Street: Postal Code:	
Shipping information: Street: Postal Code: Tel.	
Ordered For:	
Quantity: 1	
PRICE FOR STANDARD CONFIGURATION BAFFIN neoSIT HL INDOOR (STD): 1690 EUR	
PRICE FOR STANDARD CONFIGURATION BAFFIN neoSIT HL OUTDOOR (STD): 1790 EUR	

COLOUR

- ○ Blue
- ○ Red
- ○ Grey
- ○ Purple
- ○ Black

BASIC PARAMETERS (cm/inch)	SIZE	○ Baffin neoSIT HL SIZE S	○ Baffin neoSIT HL SIZE M	○ Baffin neoSIT HL SIZE L
Seat depth	A	16-31 / 6.3" - 12.2"	20-31 / 7.9" – 12.2"	27-40 /10.6"-15.7"
Back length	B	42 / 16.5"	50 / 19.7"	58 / 22.8"
Seat width	C	15-30 / 5.9" - 11.8"	15-30 / 5.9" - 11.8"	20-32 / 7.9"-12.6"
User's height	D	80-100 / 31.5"- 39.3"	100-120 / 39.3"- 47.2"	115-140 / 45.2"- 55.1"
Maximum User's weight		60 kg /132 lb	60 kg /132 lb	60 kg /132 lb

DIMENSIONS			**ADDITIONAL NOTES**
A		cm / inches	
B		cm / inches	
C		cm / inches	
D		cm / inches	
E		cm / inches	
F		cm / inches	
G		cm / inches	
Weight		kg / lb	



CLEAR

Code no.	Pict. no.	Name	Comments	Qty.	Price (EUR)	
BAFHL-I-S BAFHL-I-M BAFHL-I-L	-	Orthopedic indoor seat including upholstery [Siedzisko z tapicerką]		1	STD	
BAFHL-O-S BAFHL-O-M BAFHL-O-L	-	Orthopedic outdoor seat including upholstery [Siedzisko z tapicerką]		1	STD	
HEADRESTS – optional [Zagłówki]						
BAFHL-H-01	1	Butterfly headrest with multi adjustable headrest hardware [Zagłówek z mocowaniem]		1		☐
BAFHL-H-02	2	Headrest enlarged on both sides [Zagłówek powiększony dwustronnie]		1		☐
BAFHL-H-03	3	Headrest right side enlarged [Zagłówek powiększony prawostronnie]		1		☐
BAFHL-H-04	4	Headrest left side enlarged [Zagłówek powiększony lewostronnie]		1		☐
BAFHL-H-05	5	Anatomical headrest [Zagłówek anatomiczny]		1		☐
LATERAL SUPPORTS – choose one option [Peloty piersiowe]						
BAFNS-L-01-S BAFNS-L-01-M BAFNS-L-01-L	6	Lateral supports [Peloty piersiowe] size S ○ M ○ L ○	2 pcs.	1	STD	
BAFNS-L-91	-	Additional foam for lateral supports [Dodatkowe wypełnienie pelot piersiowych]	Upgrade option	2		☐
TRUNK AND HIP SUPPORTS [Wsparcie tułowia I bioder]						
BAFHL-V-01-S BAFHL-V-01-M BAFHL-V-01-L	7	Vest [Kamizelka] Size S ○ M ○ L ○		1	STD	
BAFHL-B-01-S BAFHL-B-01-M BAFHL-B-01-L	8	Pelvic belts [Pasy stabilizujące miednicę] Size S ○ M ○ L ○		1	STD	
BAFHL-L-02	9	Hip supports with armrests [Peloty biodrowe z podłokietnikami]	2 pcs.	1	STD	
BAFHL-L-92	-	Additional foam for hip supports [Dodatkowe wypełnienie pelot biodrowych]	Upgrade option	2		☐
BAFHL-L-03	10	Knee supports with legs abduction belts [Peloty z pasami odwodzącymi uda]	2 pcs.	1	STD	
BAFHL-A-01	11	Spine lock [Blokada kręgosłupa]		1		☐
FOOT SUPPORT ACCESSORIES [Akcesoria stóp]						
BAFHL-F-05-S/M BAFHL-F-05-L	16	Footplate Size S/M ○ L ○	1 pc	1	STD	
BAFHL-F-01-S BAFHL-F-01-M BAFHL-F-01-L	12	3D rotational foot platforms [Platformy stóp 3d] Size S ○ M ○ L ○	2 pcs.			
BAFHL-F-04-1 BAFHL-F-04-2 BAFHL-F-04-3	14	Foot platforms height adjustment [Regulacja wysokości platform/podnóżka] Range 1 ○ 2 ○ 3 ○	2 pcs.	1	STD	
BAFHL-F-02	13	Foot straps [Paski do stóp]	4 pcs.	1	STD	
BAFHL-F-03-S BAFHL-F-03-M BAFHL-F-03-L	15	3D rotational foot platforms with heel supports [Platformy stóp 3d z zapiętkami] Size S ○ M ○ L ○	Upgrade option	1		☐
ADDITIONAL ACCESSORIES [Dodatkowe akcesoria]						
BAFHL-T-01-S BAFHL-T-01-M/L	17	Tray (acrylic) [Stolik] Size S ○ M/L ○		1		☐
BAFHL-A-02	18	Polyamide adapter		1		☐
BAFHL-A-03	-	Additional seat filling [Nakładka skracająca siedzisko]		1		☐
BAFHL-A-04-S BAFHL-A-04-M BAFHL-A-04-L	-	Extra set of standard upholstery [Tapicerka kompletna] Size S ○ M ○ L ○		1		☐

CLEAR


PICTURES

1




Butterfly headrest
22 x 18 cm/8.6'' x 7''

2




Headrest enlarged on both sides
Up horn + 5 cm / + 1.9''
Bottom horn + 7 cm / +2.7''

3




Headrest right side enlarged
Up horn + 5 cm / + 1.9''
Bottom horn + 7 cm / +2.7''

4




Headrest left side enlarged
Up horn + 5 cm / + 1.9''
Bottom horn + 7 cm / +2.7''

5




Anatomical headrest

6



Lateral supports
Spacing 13 - 30 cm / 5.1"– 11.8"
Size S
Pad: 9 x 8 cm / 3.5"x 3.1"
Size M
Pad: 13 x 11 cm / 5.1"x 5.1"
Size L
Pad: 15 x 14 cm / 5.9"x 5.5"


7



Vest

Size S	Size M	Size L
A - 26 cm / 10.2''	A - 31 cm / 12.2''	A - 39 cm / 15.3''
B - 17 cm / 6.6''	B - 20 cm / 7.8''	B - 25 cm / 9.8''
C - 14 cm / 5.5''	C - 16 cm / 6.2 ''	C - 20 cm / 7.8''
D - 8,5 cm / 3.3''	D - 8,5 cm / 3.3''	D - 8,5 cm/ 3.3''
E - 15 cm / 5.9''	E - 19 cm / 7.4''	E - 29 cm / 11.4''


8



Pelvic belts

Size S:	Size M:	Size L:
A - 31 cm / 12.2''	A - 35 cm / 13.77''	A - 45 cm / 17.7''
B - 8 cm / 3.1 ''	B - 8 cm / 3.1''	B - 10 cm / 3.9''

9



Hip supports with armrests


Size S/M	Size L
Armrest height adjustment: 13-19 cm / 5.1"-7.5"	Armrest height adjustment: 17-23 cm / 6.6"- 9"

10



Knee supports with legs abduction belts 9x18 cm / 3.5"x7"
Spacing:
Size S/M: 20-32 cm / 7.8"-12.6"
Size L: 25-37 /9.8" – 14.6"

11



Spine Lock


12



3D rotational foot platforms

Size S	Size M	Size L
Length: 18 cm / 7''	Length: 21,5 cm / 8.4''	Length: 24 cm/ 9.4''
Width: 9 cm / 3.5''	Width: 12 cm/ 4.7''	Width: 12 cm / 4.7''

13




Foot straps
(set of 4 pcs.)

14




Dimension G	Foot platforms (12)	Footplate (16)
Range 1 (cm/inch)	10-21 / 3.9"- 8.3"	17-27 / 6.7"- 10.6"
Range 2 (cm/inch)	13-24 / 5.1" – 9.4"	28-38 / 11"- 15"
Range 3 (cm/inch)	22-35 / 8.7" – 13,8"	

15



3D rotational foot platforms with heel supports
Heel supports height: 12 cm / 4.7"

16	17	18
		
Footplate **Size S/M** 26 x 18 cm/10.2”x7” **Size L** 27 x 21 cm/10.6”x8.2”	**Tray** **Size S** 43 x 42 cm/16.9”x16.5” **Size M/L** 56 x 53 cm/22”x20.8”	**Polyamide adapter**